

Product datasheet for **RG216772**

FBXO48 (NM_001024680) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: FBXO48 (NM_001024680) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: FBXO48
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG216772 representing NM_001024680
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCATAAAAACCTCCAAGAGAAACAATAATTTAAGAGTTTCTCACACAGAAGCGAACTCTGTGGATGCTG
AGAAGGAAAAAATGAGAGTCAAACAACCTTTTTGAAGTCTGCCTGCAGAAATCACTTTTAAATTTT
CAGTCAGCTGGACATTCGGAGTCTGTGCAGGGCTTCATTGACATGCAGGAGCTGGAATGACACAATAAGA
AACAGTGACTCTTTATGAAAACCTCACTGCATGACTGTAAGAGCTGTGTGCCGAAGAGAAATAGATGATG
ATCTAGAAAAGTGTTATTCCTGGAGGGTGATACTGCTGAGGAATTACCAGAAGAGTAAAGTGAAACACGA
ATGGCTAAGTGGCAGATACAGCAACATTTGTTCTCCATTAGCCTACCAGAAAAAATCATGTACCCAATG
GATGCAGATACATGGGGGAAATTTCTAGAAGCAGAAGTGGAAAGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216772 representing NM_001024680
Red=Cloning site Green=Tags(s)

MHKNSKRNNLRVSHTEANSVDAEKEKNESQNNFFELLPAEITFKIFSQDIRSLCRASLTCRSWNDTIR
NSDSLWKPCHMTVRAVCRREIDDDLESGYSWRVILLRNYQKSKVKHEWLSGRYSNICSPIISLPEKIMYPM
DADTWGEILEAELER

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001024680

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001024680.3](#)

RefSeq Size: 2732 bp


RefSeq ORF: 468 bp

Locus ID: 554251

UniProt ID: Q5FWF7

Cytogenetics: 2p13.3

Product images:



Circular map for RG216772