

Product datasheet for **RG216796**

Aspartate beta hydroxylase (ASPH) (NM_004318) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aspartate beta hydroxylase (ASPH) (NM_004318) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Aspartate beta hydroxylase
Synonyms:	AAH; BAH; CASQ2BP1; FDLAB; HAAH; JCTN; junctin
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG216796 representing NM_004318
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCCAGCGTAAGAATGCCAAGAGCAGCGGCAACAGCAGCAGCAGCGGCTCCGGCAGCGGTAGCACGA
 GTGCGGGCAGCAGCAGCCCCGGGCCCCGAGAGAGACAAAGCATGGAGGACACAAGAATGGGAGGAAAGG
 CGGACTCTCAGGAACCTCATTCTTCACGTGTTTTATGGTATTGCATTGCTGGGCGTCTGGACATCTGTA
 GCTGTCGTTTTGGTTTTGATCTTGTGACTATGAGGAAGTTCTAGGAAAAGTCTAGGAATCTATGATGCTGATG
 GTGATGGAGATTTTGTGATGGATGATGCCAAAGTTTTATTAGGACTTAAAGAGAGATCTACTTCAGAGCC
 AGCAGTCCCGCCAGAAGAGGCTGAGCCACACACTGAGCCCCGAGGAGCAGGTTCTGTGGAGGCAGAACCC
 CAGAATATCGAAGATGAAGCAAAAGAACAATTCAGTCCCTTCTCCATGAAATGGTACACGCAGAACATG
 TTGAGGGAGAAGACTTGAACAAGAAGATGGACCCACAGGAGAACCACAACAAGAGGATGATGAGTTTCT
 TATGGCGACTGATGTAGATGATAGATTTGAGACCTGGAACCTGAAGTATCTCATGAAGAAACCGAGCAT
 AGTTACCACGTGGAAGAGACAGTTTACAAAGACTGTAATCAGGATATGGAAGAGATGATGTCTGAGCAGG
 AAAATCCAGATCCAGTGAACCAAGTAGTAGAAGATGAAAGATTGCACCATGATACAGATGATGTAACATA
 CCAAGTCTATGAGGAACAAGCAGTATATGAACCTCTAGAAAATGAAGGGATAGAAATCACAGAAGTAACT
 GCTCCCCCTGAGGATAATCCTGTAGAAGATTCACAGGTAATTGTAGAAGAAGTAAGCATTTCCTGTGG
 AAGAACAGCAGGAAGTACCACCAGAAAACAATAGAAAAACAGATGATCCAGAACAAAAGCAAAGTTAA
 GAAAAAGAAGCCTAACTTTTAAATAAATTTGATAAGACTATTAAGCTGAACCTGATGCTGCAGAAAAA
 CTCCGTAAGGGGAAAAATTGAGGAAGCAGTGAATGCATTTAAGAAGTACGACAAAATCCCTCAGA
 GTCCACGAGCAAGATATGGGAAGGCGCAGTGTGAGGATGATTTGGCTGAGAAGAGGAGAAGTAATGAGGT
 GCTACGTGGAGCCATCGAGACCTACCAAGAGGTGGCCAGCCTACCTGATGTCCCTGCAGACCTGCTGAAG
 CTGAGTTTGAAGCGTCGCTCAGACAGGCAACAATTTCTAGGTCATATGAGAGGTTCCCTGCTTACCTGC
 AGAGATTAGTTCAACTATTTCCAATGATACTTCTTAAAAATGACCTTGGCGTGGGATACCTCTTGAT
 AGGAGATAATGACAATGCAAAGAAAGTTTATGAAGAGGTGCTGAGTGTGACACCTAATGATGGCTTTGCT
 AAAGTCCATTATGGCTTCATCCTGAAGGCACAGAACAAAATTGCTGAGAGCATCCCATATTTAAAGGAAG
 GAATAGAATCCGGAGATCCTGGCACTGATGATGGGAGATTTTATTTCCACCTGGGGGATGCCATGCAGAG
 GGTGGGAACAAGAGGCATATAAGTGGTATGAGCTTGGGCACAAGAGAGGACACTTTCATCTGTCTGG
 CAACGCTCACTCTACAATGTGAATGGACTGAAAGCACAGCCTTGGTGGACCCAAAAGAAACGGGCTACA
 CAGAGTTAGTAAAGTCTTTAGAAAGAAACTGGAAGTTAATCCGAGATGAAGGCCTTGCAGTATGGATAA
 AGCCAAAGGTCTTCTTCTGCCTGAGGATGAAAACCTGAGGGAAAAAGGGGACTGGAGCCAGTTACAGCTG
 TGGCAGCAAGGAAGAAGAAATGAAAATGCCTGCAAAGGAGCTCCTAAAACCTGTACCTTACTAGAAAAGT
 TCCCCGAGACAACAGGATGCAGAAGAGGACAGATCAAAATATCCATCATGCACCCCGGGACTCACGTGTG
 GCCGCACACAGGGCCCAAACTGCAGGTTCCGAATGCACCTGGGCTTGGTGATTCCAAGGAAGGCTGC
 AAGATTGATGTGCCAACGAGACCAAGACCTGGGAGGAAGGCAAGGTGCTCATCTTTGATGACTCCTTTG
 AGCAGGAGGTATGGCAGGATGCCTCATCTTTCCGGCTGATTTTCATCGTGGATGTGTGCCATCCGGAAC
 GACACCACAGCAGAGACGCAGCCTTCCAGCAATT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG216796 representing NM_004318
 Red=Cloning site Green=Tags(s)

```

MAQRKNAKSSGNSSSSSGSGSGSTSAGSSSPGARRETKHGGHKNGRKGGLSGTSFFTFWVMIALLGVWTSV
AVVWFDLVDYEEVLGKGLGIYDADGDGDFDVKVLLGLKERSTSEPAVPPEEAEPHTEPEEQVPVEAEP
QNIIEDEAKEQIQSLLHEMVHAEHVEGEDLQQEDGPTGEPQQEDEFMATDVDDRFETLEPEVSHEETEH
SYHVEETVSQDCNQDMEEMMSEQENPDSSSEPVEDERLHHDTDDVYQVYEEQAVYEPLNEGIEITEVT
APPEDNPVEDSQVIVEEVSIFPVEEQEVPPEPNRKTDDPEQKAKVKKKKPKLLNKFDKTIKAELDAAEK
LRKRKIEEAVNAFKELVRKYPQSPRARYGKAQCEDDLAEKRRSNEVLRGAIETYQEASLPDVPADLLK
LSLKRRSDRQQFLGHMRGSLTLQRLVQLFPNDTSLKNDLGVGYLLIGDNDNAKKVYEEVLSVTPNDGFA
KVHYGFIKKAQNKIAESIPYLKEGIESGDPGTDGGRFYFHLGDAMQRVGNKEAYKWYELGHKRGHFASVW
QRSLYNVNLKAQPWWTPKETGYTELKSLERNWKLIRDEGLAVMDKAKGLFLPEDENLREKGDWSQFTL
WQQGRRNENACKGAPKTCTLLEKFPETTGCRGQIKYSIMHPGTHVWPHTGPTNCRFRMHLGLVIPKEGC
KIRCANETKTWEEGKVLIFDDSFHEVWQDASSFRLIFIVDVWHPPELTPQRRSLPAI
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_004318

ORF Size: 2274 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004318.2](#), [NP_004309.2](#)

RefSeq Size: 2452 bp

RefSeq ORF: 2277 bp

Locus ID: 444

UniProt ID: [Q12797](#)


Cytogenetics: 8q12.3

Domains: TPR, Asp_Arg_Hydrox, Asp-B-Hydro_N

Protein Families: Druggable Genome, Transmembrane

Gene Summary: This gene is thought to play an important role in calcium homeostasis. The gene is expressed from two promoters and undergoes extensive alternative splicing. The encoded set of proteins share varying amounts of overlap near their N-termini but have substantial variations in their C-terminal domains resulting in distinct functional properties. The longest isoforms (a and f) include a C-terminal Aspartyl/Asparaginyl beta-hydroxylase domain that hydroxylates aspartic acid or asparagine residues in the epidermal growth factor (EGF)-like domains of some proteins, including protein C, coagulation factors VII, IX, and X, and the complement factors C1R and C1S. Other isoforms differ primarily in the C-terminal sequence and lack the hydroxylase domain, and some have been localized to the endoplasmic and sarcoplasmic reticulum. Some of these isoforms are found in complexes with calsequestrin, triadin, and the ryanodine receptor, and have been shown to regulate calcium release from the sarcoplasmic reticulum. Some isoforms have been implicated in metastasis. [provided by RefSeq, Sep 2009]

Product images:



Circular map for RG216796