

Product datasheet for **RG216887**

CLDN19 (NM_148960) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLDN19 (NM_148960) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CLDN19
Synonyms:	HOMG5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG216887 representing NM_148960 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCA**ACTCAGGCCTCCAGCTCCTGGGCTACTTCTTGGCCCTGGGTGGCTGGGTGGGCATCATTGCTA**
GCACAGCCCTGCCACAGTGAAGCAGTCTTCTACGCAGGCGACGCCATCACTACTGCCGTGGGCCCTA
TGAAGGGCTCTGGATGTCCTGCGCCTCCAGAGCACTGGCAAGTGCAGTGAAGCTCTACGACTCGCTG
CTCGCCCTGGACGGTCACATCCAATCAGCGCGGGCCCTGATGGTGGTGGCCGTGCTCTGGGCTTCGTGG
CCATGGTCTCAGCGTAGTTGGCATGAAGTGTACGCGGGTGGGAGACAGCAACCCATTGCCAAGGGCCG
TGTGGCATCGCCGGGGAGCCCTTTCATCCTGGCAGGCTCTGCACTTGGACTGCTGTCTCGTGGTAT
GCCACCCTGGTGACCCAGGAGTTCTTCAACCCAAGCACACCTGTCAATGCCAGGTATGAATTTGGCCAG
CCCTGTTCTGGTGGGCTGGGCTCAGCTGGCCTGGCCGTGCTGGGCGGCTCCTTCTCTGCTGCACATGCC
GGAGCCAGAGAGACCCAACAGCAGCCACAGCCCTATCGGCCTGGACCCTCTGCTGCTGCCCGAGAACCA
GTTGTTAAATTGCCCGCCTCCGCCAAGGGCCCCCTGGGTGTG

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG216887 representing NM_148960
 Red=Cloning site Green=Tags(s)

MANSGLLQLLGYFLALGGWVGIIASTALPQWKQSSYAGDAIITAVGPYEGLWMSQSTGQVQCKLYDSL
 LALDGHIQSARALMVAVLLGFVAMVLSVVGKCTRVGDSNPIAKGRVAIAGGALFILAGLCTLTAVSWY
 ATLVTQEFFNPSTPVNARYEFGPALFVGWASAGLAVLGGSFLLCCTCEPERPNSSQPYPYRGPSSAAAREP
 VVKLPASAKGPLGV

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_148960

ORF Size: 672 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_148960.1](#), [NP_683763.1](#)

RefSeq Size: 2859 bp

RefSeq ORF: 675 bp

Locus ID: 149461

UniProt ID: [Q8N6F1](#)


Cytogenetics: 1p34.2

Protein Families: Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction

Gene Summary: The product of this gene belongs to the claudin family. It plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity. Defects in this gene are the cause of hypomagnesemia renal with ocular involvement (HOMGO). HOMGO is a progressive renal disease characterized by primary renal magnesium wasting with hypomagnesemia, hypercalciuria and nephrocalcinosis associated with severe ocular abnormalities such as bilateral chorioretinal scars, macular colobomata, significant myopia and nystagmus. Alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jun 2010]

Product images:



Circular map for RG216887