

Product datasheet for **RG216959**

CDV3 (NM_017548) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDV3 (NM_017548) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CDV3
Synonyms:	H41
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG216959 representing NM_017548 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGAGACGGAGGAGCGGAGCCTGGACAACCTCTTTGCCAAGAGGGACAAGAAGAAGAAGGAGC
GGAGCAACCGGGCGCGAGTGCCGCGGGCGCAGCGGGCAGCGCCGGCGGAAGCAGTGGAGCCGGGTGC
GGCGGGCGGGCGGGCGGGGACCCGGCCGGGTGACGGCGGGACCGCCAGCGCGGGGCTGCGGGC
CCAGGGCGGCCACCAAGGCTGTGACGAAGGACGAAGATGAATGGAAAGAATTGGAGCAAAAAGAGTTG
ATTACAGCGGCTCAGGGTTCAGGCAATGCAAAT AAGCAGTGAAAAGGAAGAAGACGATAATGAAAAGAG
ACAAGATCCAGGTGATAACTGGGAAGAAGTGGAGGTGGTGGAGGTATGGAAAAATCTTCAGGTCCC
TGGAAATAAACAGCTCCAGTACAAGCACCTCCTGCTCCAGTAATTGTTACAGAAACCCAGAACCCAGCGA
TGACTAGTGGTGTGTATAGGCCTCCTGGGGCCAGGTAAACCACAACAAGGAAAACACCACAAGGACCACC
AGAAATCTACAGTGATACACAGTTCCTATCCCTGCAGTCAACTGCCAAGCATGTAGAAAGCCGGAAGGAT
AAAGAAATGGAGAAGAGCTTTGAAGTAGTAAGACACAAAATAGAGGTAGGGATGAGGTTTCAAAAACC
AGGCCCTAAACTTCAGCTAGACAACCAATATGCTGTGCTTAAAATCAGAAAAGCAGCCACTCACATA
CAAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG216959 representing NM_017548
Red=Cloning site Green=Tags(s)

MAETEERSLDNFFAKRDKKKKKERSNRAASAAGAAGSAGGSSGAAGAAGGGAGAGTRPGDGGTASAGAAG
 PGAATKAVTKDEDEWKELEQKEVDYSGLRVQAMQISSEKEEDDNEKRQDPGDNWEEGGGGGGMEKSSGP
 WNK TAPVQAPPAPVIVTETPEPAMTSGVYRPPGARLTTTRKTPQGPPEIYSDTQFPSLQSTAKHVESRKD
 KEMEKSFVVRHKNRGRDEVSKNQALKQLDNQYAVLENQKSSHSQYN

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_017548

ORF Size: 3473 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017548.5
RefSeq Size:	3338 bp
RefSeq ORF:	777 bp
Locus ID:	55573
UniProt ID:	Q9UKY7
Cytogenetics:	3q22.1