

Product datasheet for **RG216980**

GGPS1 (NM_001037278) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: GGPS1 (NM_001037278) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: GGPS1
Synonyms: GGPPS; GGPPS1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG216980 representing NM_001037278
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGCATAATGCCAGTTTACTCATCGATGATATTGAAGACAACCTCAAACTCCGACGTGGCTTCCAG
TGGCCACAGCATCTATGGAATCCCATCTGTCAATCTGCCAATTACGTGTATTTCTTGGCTTGA
GAAAGCTTAACCTTGATCACCCAGATGCAGTGAAGCTTTTACCCGCCAGCTTTTGGAACTCCATCAG
GGACAAGGCTAGATATTTACTGGAGGGATAATTACACTTGTCCCACTGAAGAAGATATAAAGCTATGG
TGCTGCAGAAAACAGGTGGACTGTTGGATTAGCAGTAGGTCTCATGCAGTTGTTCTCTGATTACAAAGA
AGATTTAAAACCGCTACTTAATACACTTGGGCTCTTTTTCCAAATTAGGGATGATTATGCTAATCTACAC
TCCAAAGAATATAGTGAAAACAAAAGTTTTTGTGAAGATCTGACAGAGGGAAAGTTCTCATTTTCTACTA
TTCATGCTATTTGGTCAAGGCCTGAAAGCACCCAGGTGCAGAATATCTTGCGCCAGAGAACAGAAAACAT
AGATATAAAAAAATACTGTGTACATTATCTTGAGGATGTAGGTTCTTTGAATACACTCGTAATACCCCTT
AAAGAGCTTGAAGCTAAAGCCTATAAACAGATTGATGCACGTGGTGGAAACCCTGAGCTAGTAGCCTTAG
TAAACACTTAAGTAAGATGTTCAAAGAAGAAAATGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG216980 representing NM_001037278
 Red=Cloning site Green=Tags(s)

MLHNASLLIDDIEDNSKLRGFPVAHSIYGIPSVINSANYVYFLGLEKVLTDHPDAVKLFTRQLLELHQ
 GQGLDIYWRDNYTCPTEEEYKAMVLQKTGGLFGLAVGLMQLFSDYKEDLKPLLNTLGLFFQIRDDYANLH
 SKEYSENKSFCEDLTEGKFSFPTIHAIWSRPESTQVQNILRQRTEIDIKKYCVHYLEDVGSFEYTRNTL
 KELEAKAYKQIDARGGPELVALVKHLSKMFKEENE

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001037278

ORF Size: 738 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001037278.2
RefSeq Size:	2828 bp
RefSeq ORF:	741 bp
Locus ID:	9453
Cytogenetics:	1q42.3
Protein Pathways:	Metabolic pathways, Terpenoid backbone biosynthesis
Gene Summary:	This gene is a member of the prenyltransferase family and encodes a protein with geranylgeranyl diphosphate (GGPP) synthase activity. The enzyme catalyzes the synthesis of GGPP from farnesyl diphosphate and isopentenyl diphosphate. GGPP is an important molecule responsible for the C20-prenylation of proteins and for the regulation of a nuclear hormone receptor. Alternate transcriptional splice variants, both protein-coding and non-protein-coding, have been found for this gene. [provided by RefSeq, Sep 2010]