

## Product datasheet for **RG216996**

### **ANKRD20A3 (NM\_001012419) Human Tagged ORF Clone**

#### **Product data:**

|                                  |                         |
|----------------------------------|-------------------------|
| <b>Product Type:</b>             | Expression Plasmids     |
| <b>Tag:</b>                      | TurboGFP                |
| <b>Symbol:</b>                   | ANKRD20A3               |
| <b>Synonyms:</b>                 | MGC176486; MGC198508    |
| <b>Mammalian Cell Selection:</b> | Neomycin                |
| <b>Vector:</b>                   | pCMV6-AC-GFP (PS100010) |
| <b>E. coli Selection:</b>        | Ampicillin (100 ug/mL)  |



[View online »](#)

ORF Nucleotide Sequence: >RG216996 representing NM\_001012419  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGAAGTTATTCGGCTTCGGGAGCCGAGGGGCCAGACGGCCAGGGCTCCATAGACCAGTCTACACGG  
 GTTCCGGATACCGAATCCGGGACTCCGAATGCAGAAGATCCACAGGGCAGCTGTCAAAGGCGACGCCGC  
 GGAGGTGGAGCGCTGCCTGGCGCGCAGGAGCGGAGACCTGGACGCCCTGGACAAGCAGCACAGAAGTCT  
 CTACACTTGGCCTGTGCCAGTGGCCATGTGCAAGTGGTCACTCTCCTGGTTAACAGAAAAATGCCAGATTG  
 ATGCTGTGACAAAGAAAACAGAACGCCTTTGATACAGGCTGTCCATTGCCAGGAAGAGGCTTGTGCCGT  
 TATTCTGCTGGAACATGGCGCAATCCAAACCTTAAGGATATCTACGGCAACACTGCTCTCCATTATGCC  
 GTGTATAGTGAGAGCACCTCACTGGCAGAAAACTGCTTCCCATGGTGACATATTGAAGCACTGGACA  
 AGGACAATAATACCCACTTTTATTGCTATAATTTGCAAGAAAGAGAAAAATGGTGAATTTTTATTGAA  
 AAGGAAAGCAAGTTCACATGCCGTTGATAGGCTGAGACGGTCAGCTCTCATGCTTGTGTATACTATGAC  
 TCACCAAGGTATTGTCAATATCCTTCTTAAGCAAAATATTGATGTCTTCGCTCAAGACATGTGTGGACGAG  
 ATGCAGAAGATTATGCTATTTCTCATCTTTGACAAAAATTCACAAACAAATTTTGGACATAAAAAAGAA  
 GATACTTAAAAAGGAGAAAAATCAGATGTTGGAAGTTCGATGAATCTGCAGTCAGCATTTTCCATGAAGT  
 CGTGTGGATTCACTTGCCTGCATCGGATGACAAAGACTTGAATGTTGCTACTAAGCAGTGTGTCCCGGAGA  
 AAGTGTGAGAGCCTTTACCTGGATCTTCGATGAAAAAGGAAACAGAATAGTCAATGGACAAGGAGAAGG  
 GCCTCTGCAAAACATCCTTCTTGAAGCCTAGCACTGAAGTAGAAGATCCTGCTGTGAAAGGAGCAGTA  
 CAAAGAAAGAAATGTACAGACATTGAGAGCAGAACAAGCCTTACCAGTGGCTTCAGAGGAAGAGCAAGAAA  
 GGCATGAAAGAAAGTAAAGAAAGCAACACAGGTCAAAGAAAGAAATAATACAAACAAAAGTGAAAAAT  
 ACAACTTTCAGAAAATATATGTGATAGTACATCTTCTGCTGCTGCTGGCAGATTAACCCAACAAAGAAAG  
 ATTGGGAAAACGTATCCTCAGCAATTTCCCAAGAAGCTGAAGGAAGAGCATGATAGATGCACCTTAAAC  
 AAGAAAATGAAGAAAAACAATGTTAATATGCTGTACAAAAAATAGAGAAGAAATTAGAAAGGAAGA  
 GAAACAATATAAGAAAGAAAGTTGAAGCAAAACAACCTGAACCAACTGTTCACTAGAGATGAATCA  
 AAGACTGCAAGAAATACTCCAAATCGGGATTTTCATAATCATGAAGAAATGAAAGGTCTGATGGATGAAA  
 ATTGCATTTTGAAGGCAGATATTGCTATACTCAGACAGGAAATATGTACAATGAAAAATGACAACCTTGA  
 AAAAGAAAAATAATCTTAAGGACATTAATTTGTTAAAGAAACAATGCTGCCCTTGAAGATATATA  
 AAACCTCAATGAGGAAATGATAACAGAAACAGCATTCCGGTATCAACAAGAGCTTAATTATCTCAAGGCTG  
 AGAATACAAGGCTCAATGCCGAACCTGTTGAAGGAAAAAGAAAGCAAGAAAGACTGGAAGCTGACATTGA  
 ATCTTATCAGTCTAGACTGGCTGCTGCTATAAGCAAAACACAGTGAAAGTGTGAAAACAGAAAGAAACCTA  
 AAACCTGCTTTAGAGAGAACACGAGATGTTTCTGTACAAGTAGAATGAGTTCGCTATTTCCAAAGTAA  
 AAGATGAGAATGAGTTTCTTACTGAACAACTTTCTGAAACACAAATTAATTCATGCCTTAAAGATAA  
 GTTCCGTAAAGACAAGAGATAGTCTCAGAAAAAGTCATTGGCTTTAGAACTGTACAAAACGACCTAAGC  
 CAAACACAGCAGCAAAACAGGAAATGAAAGAGATGTATCAAAATGCAGAAGCTAAAGTGAATAATTCCA  
 CTGGAAGTGAAGTGTGTAGAAGAGAGGATATGTCACCTCCAACGTGAAAATGCGTGGCTTGTACAGCA  
 ACTAGATGACGTTATCAGAAAGAGGATCATAAGAGATAGTAACATAATCCAAAGAGGCTTTATTGAG  
 AGTGGAAAGAAAGACCTCGTGCTAGAAGAGAAAAAGTAAGAAGCTAATGAATGAATGTGATCATTTAAAG  
 AAAGTCTCTTTAGTATGAGAGAGAGAAAAACAGAGGAGTAGTAAGTATCAAGGAAGATAAATTTTCA  
 AACTCTAGAAAGACAATT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG216996 representing NM\_001012419

Red=Cloning site Green=Tags(s)


MKLFGFGSRRGQTAQGSIDHVTGSGYRIRDSELQKIHRAAVKGDAAEVERCLARRSGDLDALDKQHR  
LHLACASGHVQVVTLLVNRKCQIDVCDKENRTPLIQAVHCQEEACAVILLEHGANPNLKDIYGNALHYA  
VYSESTSLAEKLLSHGAHIEALDKDNNTPLLFAIICKKEKMEVFLKRRASSHAVDRLRRSALMLAVYYD  
SPGIVNILLKQIDVFAQDMCGRDAEDYAISSHLLTKIQQILEHKKKILKKEKSDVGSSDESAVSIFHEL  
RVDSLPAASDDKDLNVATKQCVPEKVSEPLPGSSHEKGNRIVNGQGEGPPAKHPSLKPSTEVEDPAVKGAV  
QRKNVQTLRAEQALPVASEEEQERHERSEKKQPQVKEGNNNTNKSEKIQLSNICDSTSSAAAGRLTQQRK  
IGKTPQQFPKLLKEEHRCTLKQENEEKTNVNMLYKKNREELERKEKQYKKEVEAKQLEPTVQSLEMKS  
KTARNTPNRDFHNHEEMKGLMDENCILKADIAILRQEICTMKNDNLEKENKYLKDIKIVKETNAALEKYI  
KLNEEMITETAFRYQQELNYLKAENTRLNAELLKEKESKKRLEADIESYQSRLAAAIKHSSESVKTERNL  
KLALERTRDVSQVEMSSAISKVKDENEFLEQLSETQIKFNALKDKFRKTRDSLRRKSLALETQVNDLS  
QTQQQTQEMKEMYQNAEAKVNNSTGKWCVEERICHLORENALVQQLDDVHQKEDHKEIVTNIQRGFIE  
SGKKDLVLEEKSKKLMNECDHLKESLFQYEREKTEGVVSIKEDKYFQTSRKTI

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM\_001012419

ORF Size:

2469 bp

OTI Disclaimer:


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <u>NM_001012419.2</u> , <u>NP_001012419.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq Size:</b>           | 3964 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 2472 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 441425                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>            | <u>Q5VUR7</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 9q21.11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

**Product images:**



Circular map for RG216996