

Product datasheet for **RG217025**

PDE9A (NM_001001575) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PDE9A (NM_001001575) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PDE9A
Synonyms: HSPDE9A2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG217025 representing NM_001001575
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACGCATT CAGAAGAGCATT CAAAA TCAATGA AACTGAAAGCTGAAGTTGCAAATCACTTGGCTGTCC
TAGAGAAACGCGTGAATTGGAAGGACTAAAAGTGGTGGAGATTGAGAAATGCAAGAGTGACATTAAGAA
GATGAGGGAGGAGCTGGCGGCCAGAAGCAGCAGGACCAACTGCCCTGTAAGTACAGTTTTTTGGATAAC
CACAAGAAGTTGACTCCTCGACGCGATGTTCCCACTTACCCCAAGTACCTGCTCTCCAGAGACCATCG
AGGCCCTGCGGAAGCCGACCTTTGACGTCTGGCTTTGGGAGCCCAATGAGATGCTGAGCTGCCTGGAGCA
CATGTACCACGACCTCGGGCTGGTCAGGGACTTCAGCATCAACCCGTGCACCCCTCAGGAGGTGGCTGTT
TGCCTCCACGACAACACAGAAACAACCCCTTCCACAACCTCCGGCACTGCTTCTGCGTGGCCCAGATGA
TGTACAGCATGGTCTGGCTCTGCAGTCTCCAGGAGAAGTTCTCACAAACGGATATCCTGATCCTAATGAC
AGCGGCCATCTGCCACGATCTGGACCATCCCGGCTACAACAACACGTACCAGATCAATGCCCGCACAGAG
CTGGCGGTCCGCTACAATGACATCTCACCGCTGGAGAACCACCACTGCGCCGTGGCTTCCAGATCCTCG
CCGAGCCTGAGTGCAACATCTTCTCCAACATCCCACTGATGGTTCAAGCAGATCCGACAGGGAATGAT
CACATTAATCTTGCCACTGACATGGCAAGACATGCAGAAATATGGATTCTTTCAAAGAGAAAAATGGAG
AATTTTGACTACAGCAACGAGGAGCACATGACCCTGCTGAAGATGATTTTATAAAATGCTGTGATATCT
CTAACGAGGTCCGTCCAATGGAAGTCGACAGACCTTGGGTGGACTGTTTATTAGAGGAATATTTTATGCA
GAGCGACCGTGAGAAGTCAGAAGGCCTTCTGTGGCACCGTTCATGGACCAGACAAAGTACCAAGGCC
ACAGCCCAGATTGGGTTTCATCAAGTTTGTCTGATCCCAATGTTTAAAACAGTGACCAAGCTCTTCCCCA
TGGTTGAGGAGATCATGCTGCAGCCACTTTGGGAATCCCGAGATCGCTACGAGGAGCTGAAGCGGATAGA
TGACGCCATGAAAGAGTTACAGAAGAAGACTGACAGCTTGACGTCTGGGGCCACCGAGAAGTCCAGAGAG
AGAAGCAGAGATGTGAAAAACAGTGAAGGAGACTGTGCC

ACGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)


Protein Sequence: >RG217025 representing NM_001001575
 Red=Cloning site Green=Tags(s)

MDAFRRAFKINELKAEVANHLAVLEKRVELEGLKVVEIEKCKSDIKKMREELAARSSRTNCPCKYSFLDN
 HKKLTTPRRDVPTYPKYLLSPETIEALRKPTFDVWLWEPNEMLSCLEHMYHDLGLVRDFSINPVTLRRWLF
 CVHDNYRNNPFHNRHCFCAVQMMYSMVWLCSLQEKFSQTDILILMTAAICHDLDHPGYNNTYQINARTE
 LAVRYNDISPLENHHCVAFAFQILAEPECNIFSNIPDPGFKQIRQGMITLILATDMARHAEIMDSFKEKME
 NFDYSNEEHMTLLKMILIKCCDISNEVRPMEVAEPWVDCLLEEYFMQSDREKSEGLPVAPFMDRDKVTKA
 TAQIGFKFVLIIPMFETVTKLFPMVEEIMLQPLWESRDREELKRIDDAMKELQKKTDSLTSGATEKSRE
 RSRDVKNSEGDCA


TRPLE - GFP Tag - V

Restriction Sites: SgfI-NotI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001001575

ORF Size:	1299 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001575.1 , NP_001001575.1
RefSeq Size:	1675 bp
RefSeq ORF:	1302 bp
Locus ID:	5152
UniProt ID:	O76083
Cytogenetics:	21q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism
Gene Summary:	The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]