

## Product datasheet for **RG217220**

### **KHSRP (NM\_003685) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	KHSRP (NM_003685) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	KHSRP
Synonyms:	FBP2; FUBP2; KSRP; p75
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RG217220 representing NM\_003685  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGTCGGACTACAGCACGGGAGGACCCCGCCCGGGCCGCGCCCGCCCGCGGGGGCGGGGAGCCG  
GAGGCGCCGGGGAGGCCCTCCGCGGGCCCGCAGGCGGGGGACCGGGCGGCGCGGTCCCTGCGG  
CGGCGGCCGGGGCGGGGTGCGCCGGGGCCCTCTCAGCCACCCGGGAGGCGGGCCCGGAATCCGC  
AAGGACGCTTTCGCCGACGCGTGCAGCGGGCCCGCCAGATTGCAGCCAAAATTGGAGGCGATGCTGCCA  
CGACAGTGAATAACAGCACTCCTGATTTTGGTTTTGGGGCCAAAAGAGACAGTTGGAAGATGGAGATCA  
ACCGGAGAGCAAGAAGCTGGCTTCCAGGGAGACTCAATCAGTTCTCACTTGGACCCATCCATCCTCCC  
CCAAGGACTTCAATGACAGAAGAGTACAGGGTCCAGACGGCATGGTGGGCTGATCATTGGCAGAGGAG  
GTGAACAAATTAACAAAATCCAACAGGATTCAGGCTGCAAAGTACAGATTTCTCCAGACAGCGGTGGCT  
ACCGGAGCGCAGTGTGCTTTCAGAGGACCCAGAAATCTGTCCAGAAAGCAAGATGATGCTGGATGAC  
ATTGTGTCTCGGGTCTGTGGGGCCCCCAGGACAGTTCCACGACAACGCCAACGGGGCCAGAACGGCA  
CCGTGCAGGAGATCATGATCCCCGCGGCAAGGCCGCGCTGGTCAATTGGCAAGGGCGGGGAGACCATTAA  
GCAGCTGCAGGAACGCGCTGGAGTGAAGATGATCTTAATTGAGGACGGATCTCAGAATACGAATGTGGAC  
AAACCTCTCCGCATATTGGGGATCCTTACAAAGTGCAGCAAGCCTGTGAGATGGTATGGACATCCTCC  
GGGAACGTGACCAAGGCGGCTTTGGGGACCGGAATGAGTACGGATCTCGGATTGGCGGAGGCATCGATGT  
GCCAGTGGCCAGGCATTTCTGTTGGCGTGGTCAATGGCCGAGTGGAGAGATGATCAAGAAGATCCAGAAT  
GATGCTGGCGTGGGATACAGTTCAAGCAAGATGACGGGACAGGGCCCGAGAAGATTGCTCATATAATGG  
GGCCCCAGACAGGTGCGAGCACGACGCCGGATCATCAACGACCTCCTCCAGAGCCTCAGGAGTGGTCC  
CCAGGTCTCCAGGGGTCCAGGCATGCCCCGGGGGGCCGAGGCCGAGGAAGAGGCCAAGGCAATTGG  
GGTCCCCCTGGCGGGGAGATGACCTTCTCCATCCCCACTCACAAGTGTGGGCTGGTATCGGCCGAGGTG  
GCGAGAATGTGAAAGCCATAAACAGCAGACGGGAGCCTTCGTAGAGATCTCCGGCAGCTGCCACCCAA  
CGGGGACCCCACTTCAAGTTGTTTCATCATCCGGGTTTACCCAGCAGATTGACCACGCCAAGCAGCTT  
ATCGAGGAAAAGATCGAGGGTCTCTCTGCCAGTTGGACCAGGCCAGGTGGCCAGGCCCTGCTGGCC  
CAATGGGGCCCTTCAATCCTGGGCCCTTCAACCAGGGGCCACCCGGGCTCCCCACATGCCGGGGGCC  
CCCTCCTCACCAGTACCCACCCAGGGTGGGGCAATACCTACCCAGTGGCAGCCGCTGCTCCTCAT  
GACCCAAGCAAAGCAGCTGCAGCGGCCGCGGACCCCAACGCCGCTGGGCCGCTACTACTCACACTACT  
ACCAGCAGCCCCCGGCCCGTCCCCGGCCCGCACCGGCCCTGCGGCCACCCGGCTCAGGGTGAAGCC  
CCCTCAGCCCCACCCACCGGCCAGTCGGACTACACTAAGGCCTGGGAAGAGATTACAAAAAGATCGGC  
CAGCAGCCCCAGCAGCCCGGAGCGCCCCACAGCAGGACTACACGAAGGCTTGGGAGGAGTACTACAAGA  
AGCAAGCGCAAGTGGCCACCGGAGGGGTCCAGGAGCTCCCCAGGCTCCAGCCAGACTACAGTCCCGC  
CTGGGCGGAATATTACAGACAGCAGGCCGCTTACTACGGACAGACCCAGGTCTGGCGGCCCCAGCCG  
CCGCCACGCAGCAGGGACAGCAGCAGGCTCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG217220 representing NM\_003685  
 Red=Cloning site Green=Tags(s)


```
MSDYSTGGPPPPPPAGGGGGAGGAGGGPPPPGAGDRGGGGPCGGGPGGGSAGGPSQPPGGGGPGIR
KDAFADAVQRRARQIAAKIGGDAATTVNNSTPDFGFGGQKRQLEDGDQPE SKKLASQGDSSQLGPIHPP
PRTSMTEEYRVPDGMVGLIIGRGGEQINKIQDMSGCKVQISPDSSGGLPERSVSLTGAPESVQKAKMMLDD
IVSRGRGGPPGQFHDNANGGQNGTVQEIMIPAGKAGLVIGKGGETIKQLQERAGVKMILIQDGSQNTNVD
KPLRIIGDPYKVVQACEMVMDILRERDQGGF GDRNEYGSRI GGGIDVPVPRHSVGVVIGRSGEMIKKIQN
DAGVRIQFKQDDGTGPEKIAHIMGPPDRCEHAARIINDLLQSLRSGPPGPPGMPGPPGGRGRGQGNW
GPPGGEMTFSIPTHKCLVIGRGGENVKAINQQTGAFVEISRQLPPNGDPNFKLFIIRGSPQQIDHAKQL
IEEKIEGPLCPVGPGGGPGPAGPMGPFNP GPFNQGGPAGPHAGGPPPHQYPPQGWNTYPQWQPPAPH
DPSKAAAAADPNAAWAAYYSHYYQQPPGPVGPAPAPAAPPAQGEPPQPPTGQSDYTKAWEEYKKIG
QQPQQGAPPQQDYTKAWEEYKKQAQVATGGGPGAPPGSQPDYSAAWAEYYRQQAAYYQTPGPGGPQP
PPTQQGQQAQ
```


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



**Plasmid Map:**


**ACCN:** NM\_003685

**ORF Size:** 2133 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_003685.1</a> , <a href="#">NP_003676.1</a>
<b>RefSeq Size:</b>	3009 bp
<b>RefSeq ORF:</b>	2136 bp
<b>Locus ID:</b>	8570
<b>UniProt ID:</b>	<a href="#">Q92945</a>
<b>Cytogenetics:</b>	19p13.3
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>Gene Summary:</b>	The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization (Min et al., 1997 [PubMed 9136930]; Gherzi et al., 2004 [PubMed 15175153]). [supplied by OMIM, Apr 2010]