

Product datasheet for **RG217539**

RPL31 (NM_001099693) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPL31 (NM_001099693) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: RPL31
Synonyms: L31
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG217539 representing NM_001099693
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCCCGCAAAGAAGGGTGGCGAGAAGAAAAAGGCCGTTCTGCCATCAACGAAGTGGTAACCCGAG
AATACACCATCAACATTCACAAGCGCATCCATGGAGTGGGCTTCAAGAAGCGTGCACCTCGGGCACTCAA
AGAGATTCGAAATTTGCCATGAAGGAGATGGAACTCCAGATGTGCGCATTGACACCAGGCTCAACAAA
GCTGTCTGGCCAAAGGAATAAGGAATGTGCCATACCGAATCCGTGTGCGGCTGTCCAGAAAACGTAATG
AGGATGAAGATTCACAAATAAGCTATATACTTTGGTTACCTATGTACCTGTTACCACTTCAAAGTAA
GTTCTCCATCCCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG217539 representing NM_001099693
Red=Cloning site Green=Tags(s)

MAPAKKGGEEKKGRSAINVVVREYTIINIKRIHGVGFKKRAPRALKEIRKFAMKEMGTPDVRIDTRLNK
AVWAKGIRNVPYRIRVRLSRKRNEDEDSNKLYTLVTVYVPVTTFKSKFSIP

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001099693

ORF Size: 363 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001099693.1 , NP_001093163.1
RefSeq Size:	801 bp
RefSeq ORF:	366 bp
Locus ID:	6160
UniProt ID:	P62899
Cytogenetics:	2q11.2
Protein Pathways:	Ribosome
Gene Summary:	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L31E family of ribosomal proteins. It is located in the cytoplasm. Higher levels of expression of this gene in familial adenomatous polyps compared to matched normal tissues have been observed. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]