

Product datasheet for **RG217571**

ZDHHC16 (NM_198043) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZDHHC16 (NM_198043) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZDHHC16
Synonyms:	APH2; DHHC-16
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG217571 representing NM_198043 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGAGGCCAGCGGAGCCTGCTGCTGGGCCCGCCCGCCTCTGCCTCCGCCTCCTTCTGCTGCTGGGT
ACAGGGCGCGCTGTCCACCTCTACTCCGGGTCTAGTACAGCGCTGGCGCTACGGCAAGGTCTGCCTGCC
CTCCCTGCTCTACAACCTCTTTGGGGCAGTGACACCGCTGTTGATGCTGCCTTTGAGCCTGTCTACTGG
CTGGTAGACAACGTGATCCGCTGGTTGGAGTGGTTCGTGGTCTGGTATCGTCTGACAGGCTCCA
TTGTAGCTATCGCTACCTGTGTCTGCCTCTCATCCTCCGAACCTACTCAGTGCCACGACTCTGCTG
GCATTTCTTCTATAGCCACTGGAATCTGATCCTGATTGTCTTCCACTACTACCAGCCATCACCCTCCG
CCTGGGTACCCACCCAGGGCAGGAATGATATCGCCACCGTCTCCATCTGTAAGAAGTGCATTTACCCCA
AGCCAGCCCAACACACCACTGCAGCATCTGCAACAGGTGTGTGCTGAAGATGGATCACCCTGCCCTG
GCTAAACAATTGTGTGGGCCACTATAACCATCGTACTTCTTCTTTCTGCTTTTTTCTGACTCTGGGC
TGTGTCTACTGCAGTATGGAAGTTGGGACTTTTCCGGGAGGCTTATGCTGCCATTGAGACTTATACC
AGACCCACCACCCACCTTCTCCTTTGAGAAAGGATGACTCACAAGAGTCTTGTCTACCTCTGGTTCCT
GTGCAGTCTGTGGCACTTGCCCTGGGTGCCCTAACTGTATGGCATGCTGTTCTCATCAGTCGAGGTGAG
ACTAGCATCGAAAGGCACATCAACAAGAAGGAGAGAGCTCGGCTACAGGCCAAGGGCAGAGTATTTAGGA
ATCCTTACAACCTACGGCTGCTTGGACAACCTGGAAGGTATTCTGGGTGTGGATACAGGAAGGCACTGGCT
TACTCGGGTGTCTTACCTTCTAGTCACTTGCCCATGGGAATGGAATGAGCTGGGAGCCCCCTCCCTGG
GTGACTGCTCACTCAGCCTCTGTGATGGCAGTG

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG217571 representing NM_198043
 Red=Cloning site Green=Tags(s)

MRGQRLLLLGPARLCLRLLLLLGYRRRCPLLRGLVQRWRYGKVCRLSLLYNSFGGSDTAVDAAFEFVYV
 LVDNVIWFGVVFVVLVIVLTGSIVAIAYLCVLPILRTYSVPRLCWHFFYSHWNLILIVFHYYQAITTP
 PGYPPQGRNDIATVSICKKCIYPKPARTHHCSICNRCVLKMDHHCPLNNCVGHYNHRYFFSFCFFMTLG
 CVYCSYGSWDLFREAYAAIETYHQTPPTFSFRERMTHKSLVYLWFLCSSVALALGALTVWHAVLISRGE
 TSIERHINKERRRLQAKGRVFRNPYNYGCLDNWKVFLGVDTRHWLTRVLLPSSHLPHGNGMSWEPPPV
 VTAHSASVMAV


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_198043

ORF Size:	1083 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198043.2 , NP_932160.1
RefSeq Size:	1768 bp
RefSeq ORF:	1086 bp
Locus ID:	84287
UniProt ID:	Q969W1
Cytogenetics:	10q24.1
Protein Families:	Transmembrane
Gene Summary:	Palmitoyl acyltransferase that mediates palmitoylation of proteins such as PLN and ZDHHC6 (PubMed:28826475). Required during embryonic heart development and cardiac function, possibly by mediating palmitoylation of PLN, thereby affecting PLN phosphorylation and homooligomerization (By similarity). Also required for eye development (By similarity). Palmitoylates ZDHHC6, affecting the quaternary assembly of ZDHHC6, its localization, stability and function (PubMed:28826475). May play a role in DNA damage response (By similarity). May be involved in apoptosis regulation (By similarity). Involved in the proliferation of neural stem cells by regulating the FGF/ERK pathway (By similarity).[UniProtKB/Swiss-Prot Function]