

Product datasheet for **RG217639**

SEPP1 (SELENOP) (NM_001093726) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SEPP1 (SELENOP) (NM_001093726) Human Tagged ORF Clone
Symbol:	SEPP1
Synonyms:	SELP; SeP; SEPP; SEPP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG217639 representing NM_001093726 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGCTCTGGGCTGGACGCGGTGGCTCACGCTTGTAAATCCCAGCACTTTGGGAGGCCAAGGCAGGCGGA
TCACGAGGACAACCCAGCAATGTGGAGAAGCCTGGGGCTTGCCCTGGCTCTCTGTCTCCATCGGG
AGGAACAGAGAGCCAGGACCAAGCTCCTTATGTAAGCAACCCAGCCTGGAGCATAAGAGATCAAGAT
CCAATGCTAACTCCAATGGTTCACTGACTGTGGTTGCTCTTCTCAAGCCAGCTGATACCTGTGCATAC
TGCAGGCATCTAAATTAGAAGACCTGCGAGTAAACTGAAGAAAGAGGATATTCTAATTTCTTATAT
TGTTGTTAATCATCAAGGAATCTTCTCGATTAAAATACACACATCTTAAGAATAAGGTTTCAGAGCAT
ATTCCTGTTTATCAACAAGAAGAAAACCAACAGATGTCTGGACTCTTTAAATGGAAGCAAAGATGACT
TCCTCATATATGATAGATGTGGCCGTCTGTATATCATCTTGGTTTGCCTTTTTCTTCCCTAACTTTCCC
ATATGTAGAAGAAGCCATTAAGATTGCTTACTGTGAAAAGAAATGTGGAACTGCTCTCACGACTCTC
AAAGATGAAGACTTTTGTAAACGTGTATCTTTGGCTACTGTGGATAAAACAGTTGAACTCCATCGCCTC
ATTACCATCATGAGCATCATACAATCATGGACATCAGCACCTTGGCAGCAGTGAGCTTTCAGAGAAATCA
GCAACCAGGAGCACCAATGCTCCTACTCATCCTGCTCCTCCAGGCCTTCATCACCACATAAGCACAAG
GGTCAGCATAGGCAGGGTCAACCAGAGAACCGAGATATGCCAGCAAGTGAAGATTTACAAGATTTACAAA
AGAAGCTCTGTCGAAAGAGATGTATAAATCAATTAATCTGTAAATGCCCCACAGATTAGAGTTGGCTCC
TAGGAGCTGATGCTGCCATTGTGCGACATCTGATATTTGAAAAACAGGGTCTGCAATCACCTGACAGTGT
AAAGAAAACCTCCATCTTTATGTAGCTGACAGGGACTTCGGGCAGAGGAGAACATAACTGAATCTTGTG
AGTGACGTTTGCCTCCAGCTGCCTGACAAATAAGTCAGCAGCTTATACCCACAGAAGCCAGTGCCAGTTG
ACGCTGAAAGAATCAGGCAAAAAGTGAGAATGACCTTCAAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG217639 representing NM_001093726
 Red=Cloning site Green=Tags(s)

MCSGLDAVAHACNPSTLGGQGRITRTPAMWRSLLGLALALCLLPSGGTESQDQSSLCKQPPAWSIRDQD
 PMLNSNGSVTVVALLQAS*YLCILQASKLEDLRVKKKEGYSNISIYIVVNHQGISSRLKYTHLKNKVSEH
 IPVYQEEENQTDVWTLNGSKDDFLIYDRCGRLVYHLGLPFSFLTFPYVEEAIKIAYCEKKCGNCSLTTL
 KDEFDFCKRVSLATVDKTVETPSPHYHHEHHNHGHQHLGSSELSENQQPGAPNAPTHPAPPGLHHHKKHK
 GQHRQGHPENRDMPAEQLQDLQKQLCRKRCINQLLCKLPTDSELAPRS*CCHCRHLIFEKTGSAIT*QC
 KENLPSLCS*QGLRAEENITESCQ*RLPPAA*QISQQLIPTEASAS*R*KNQAKK*E*PSN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001093726

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone encodes a selenoprotein containing the rare amino acid selenocysteine (Sec). Sec is encoded by UGA codon, which normally signals translational termination. Expression of this clone is not guaranteed due to the nature of selenoproteins.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001093726.2](#)

RefSeq Size: 2252 bp

RefSeq ORF: 1236 bp

Locus ID: 6414


UniProt ID: [P49908](#)

Cytogenetics: 5p12

Protein Families: Secreted Protein

Gene Summary: This gene encodes a selenoprotein that is predominantly expressed in the liver and secreted into the plasma. This selenoprotein is unique in that it contains multiple selenocysteine (Sec) residues per polypeptide (10 in human), and accounts for most of the selenium in plasma. It has been implicated as an extracellular antioxidant, and in the transport of selenium to extra-hepatic tissues via apolipoprotein E receptor-2 (apoER2). Mice lacking this gene exhibit neurological dysfunction, suggesting its importance in normal brain function. Sec is encoded by the UGA codon, which normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, designated the Sec insertion sequence (SECIS) element, that is necessary for the recognition of UGA as a Sec codon, rather than as a stop signal. The mRNA for this selenoprotein contains two SECIS elements. The use of alternative polyadenylation sites, one located in between the two SECIS elements, results in two populations of mRNAs containing either both (predominant) or just the upstream SECIS element (PMID:27881738). Alternatively spliced transcript variants have also been found for this gene. [provided by RefSeq, Oct 2018]

Product images:



Circular map for RG217639