

Product datasheet for **RG217667**

LDLRAD2 (NM_001013693) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LDLRAD2 (NM_001013693) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LDLRAD2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG217667 representing NM_001013693 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGCTTGTGTCTTCTGCAGTTGCCCAAAGGTTGCTCTTGCTGGGGGCAGCCGCCCTGACTGCAA
CTGCTTTGGAGACAGCCGACCTGGCGGAAGTGTGCGGGCAGAGCTGGCAGGGGGACGGGCTGCTGCTGCG
CTCGCACGCCGCATCGCGCAGGTTCTACTTCGTGGCTCCGGACACCGACTGCGGGCTCTGGGTGCAGGCG
GCAGCCCCCGCGACCGGATCCGCTTCCAGTTCGCTTCTTCTGGTCTACAGCCTGACCCCCGCGCCCC
CGGCGCTCAACACCTCTCCCCGGCCCCGGCCGACCCGTGCGCCCCGGCTCCTACCTGCGGTTCTACGA
GGGCCCGCGGGGCGCCCCGGCCCTGGGGTCCCCACTGTGCGGCCTGAACATCCCGGTGCTGTGGCA
TCCTCCGACCCCTTTCTAGGCCTGCGCCTGGTCACGAGAGGCCGCCAGCCCCGCTGGACTTCGTGGGCG
AAGTCACCTCTTCCGTCTGGGACCTTGTGGTGCCTACTTCCGCTGCCAGAATGGCAGGTGCATCCCTC
AAGCCTCGTGTGTGACCCTGGGCGATGGACAAGTGTGGCGATGGCAGTGACCAGGGCTCCTGGTACCA
GCTGACTGCAGAGGTCCCTCTCCGGTGGCCAGCCAGACAGGAAGTACAGATGCCCATACCTCCAGATCCC
TGACTCCCTCCCAGCTCTCGGTCTGCAGGATCCCTCTGGATTGCAGCTGAGAGGAGTTCCCCAGCAGG
CAGGGACCCACGAGACAAGACGCAGCTTGAAGGCTCCACTGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG217667 representing NM_001013693
 Red=Cloning site Green=Tags(s)

MEACCLLQLPQRLLLLGAAALTATALETADLAELCGQTWQGDGLLLRSHAASRRFYFVAPDTCGLWVQA
 AAPGDRIRFQFRFFLVYSLTPAPPALNTSSPAPADPCAPGSYLRFYEGPPGAPRPLGSPLCGLNIPVPVA
 SSGPFLGLRLVTRGRQPRVDFVGEVTSFRLGPCGAYFRCQNGRCIPSSLVCDPWGMDNCGDGSQGSWSP
 ADCRGPSVPVPSQTGSTDAHTSRSLTPSPALGSAGSLWIAAERSSPAGRDPTRQDAALEGSTE

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001013693

ORF Size: 816 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001013693.1, NP_001013715.1

RefSeq Size: 4012 bp

RefSeq ORF: 819 bp


Locus ID: 401944

UniProt ID: Q5SZI1

Cytogenetics: 1p36.12

Protein Families: Druggable Genome

Product images:



Circular map for RG217667