

Product datasheet for **RG217756**

BRAP (NM_006768) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRAP (NM_006768) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	BRAP
Synonyms:	BRAP2; IMP; RNF52
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG217756 representing NM_006768
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGAGTGTGCTACTGTTTGTATCCGATTGGAGCTCGCGAACACTCGCCTGTCCCCGCCGGCTTCGGCT
 TCAGCGCCGCGCGCGGGAAATGTCTGATGAGGAGATAAAAAAGACGACACTAGCCTCAGCTGTAGCCTG
 TTAGAAGGCAAGTCACCAGGAGAGAAAGTAGCGATTATCCATCAGCATCTCGGCCGTCGAGAAATGACA
 GATGTGATCATTGAGACCATGAAGTCCAACCCAGATGAACAAAACTACAGTGAAGAAAGGAAGTCTT
 CAGAAGCCTCCCCACTGCGCAAAGAAGTAAAGATCACAGTAAGGAATGCATAAACGCTGCCCCAGATTC
 TCCGTCCAACAGCTTCCAGACCAGATTTATTCTTCAGTGGAATCCATCAGTTGAAATAGTTCATGGT
 ATTATGCACCTATATAAGACAAATAAGATGACCTCCTTAAAGAAGATGTGCGGCGCAGTGCCATGCTGT
 GTATTCTCACAGTCCCTGCTCAATGACCAGTCATGACCTTATGAAGTTTGTGCCCATTTAACGAAGT
 AATTGAACAAATGAAAATTATCAGAGACTCTACTCCCAACCAATATATGGTGCTGATAAAGTTTCGTGCA
 CAGGCTGATGCGGATAGTTTTATATGACATGCAATGGCCGCCAGTTCAACTCAATAGAAGATGACGTTT
 GCCAGCTAGTGTATGTGGAAGAGCTGAAGTGCTCAAATCTGAAGATGGCGCCAGCTCCAGTGATGGA
 CCTGACTGAACCTCCCAAGTGACCGGTGTGTCTGGAGCGCATGGACGAGTCTGTGAATGGCATCCTCACA
 ACGTTATGTAAACACAGCTTCCACAGCCAGTGTCTACAGCGCTGGGACGATACCAGTGTCTGTTTGCC
 GGTACTGTCAAACGCCCGAGCCAGTAGAAGAAAATAAGTGTTTTGTAGTGTGGTGTTCAGGAAAATCTTTG
 GATTTGTTTAATATGCGGCCACATAGGATGTGGACGGTATGTGAGTCGACATGCTTATAAGCACTTTGAG
 GAAACGCAGCACACGTATGCCATGCAGCTTACCAACCATCGAGTCTGGGACTATGCTGGAGATAACTATG
 TTCATCGACTGGTTGCAAGTAAAACAGATGGAAAAATAGTACAGTATGAATGTGAGGGGGATACTTGCCA
 GGAAGAGAAAAATAGATGCCTTACAGTTAGAGTATTATATTTTACTAACCAAGCCAGCTGGAATCTCAGCGA
 ATCTACTGGGAAAAACAAGATAGTTCGGATAGAGAAGGACACAGCAGAGGAAATTAAACAATGAAGACCA
 AGTTTAAAGAAACAATTGAGAAGTGTGATAATCTAGAGCACAACTAAATGATCTCTCTAAAGAAAAAGCA
 GTCTGTGGAAAGAAAGTCACTCAGCTAAACACAAAAGTGGCCAACTCACCACGAGCTCAAAGAGGAG
 CAGGAAATGAACAAGTGTTCGAGGCCAACCAAGTCTCTGCAGAACAGCTAAAGAGGAGGAGAGGG
 TGCTGAAGGAGACCTGTGACCAAAAAGATCTGCAGATCACCGAGATCCAGGAGCAGCTGCGTGACGTCAT
 GTTCTACCTGGAGACACAGCAGAAGATCAACCATCTGCCTGCCGAGACCCGGCAGGAAATCCAGGAGGGA
 CAGATCAACATCGCCATGGCCTCGGCCTCGAGCCTGCCTCTTCGGGGGCGAGTGGGAAGTTGCCCTCCA
 GGAAGGGCCGCGCAAGAGGGGCAAG

ACGCGTACGCGGCGCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:


>RG217756 representing NM_006768
 Red=Cloning site Green=Tags(s)

MSVSLVVIRLELAHSPVPAGFGFSAAAGEMSDEEIKKTTLASAVACLEGKSPGEKVAIHQHLGRREMT
 DVIIETMKNPDELKTTVEERKSSEASPTAQRSDHSKECINAAPDSPSKQLPDQISFFSGNPSVEIVHG
 IMHLYKTNKMTSLKEDVRRSAMLCLTVPAAMTSHDLMKFVAPFNEVIEQMKIIRDSTPNQYMLIKFRA
 QADADSFYMTGNRQFNSEDVQCLVYVERAEVLKSEDGASLPVMDLTELPKCTVCLERMDESUNGILT
 TLCNHSFHSQCLQRWDDTTPVCRYCQTPPEVEENKCFECGVQENLWICLICGHIGCGRYVSRHAYKHFE
 ETQHTYAMQLTNHRVWDYAGDNYVHRLVASKTDGKIVQYCEGDTQEEKIDALQLEYSYLLTSQLESQR
 IYWENKIVRIEKDTAEEINNMKTKFKETIEKCDNLEHKLNDLLKEQSVKRCQTLNKTAKLNLNELKEE
 QEMNKCLRANQVLLQNKLEEEVRLKETCDQKDLQITEIQEQLRDVFMFLETQQKINHLPAETRQEIQEG
 QINIAMASASSPASSGGSKLP SRKGRSKRGK

TRTRPLE - GFP Tag - V

Restriction Sites:


SgfI-MluI

Cloning Scheme:


ACCN:	NM_006768
ORF Size:	1776 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_006768.5
RefSeq Size:	2063 bp
RefSeq ORF:	1779 bp
Locus ID:	8315

UniProt ID:	<u>Q7Z569</u>
Cytogenetics:	12q24.12
Protein Families:	Druggable Genome
Gene Summary:	The protein encoded by this gene was identified by its ability to bind to the nuclear localization signal of BRCA1 and other proteins. It is a cytoplasmic protein which may regulate nuclear targeting by retaining proteins with a nuclear localization signal in the cytoplasm. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG217756