

## Product datasheet for **RG217893**

### **C17orf77 (NM\_152460) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C17orf77 (NM\_152460) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** C17orf77  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG217893 representing NM\_152460  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGATGAATTGGCGCTGTCATTTCCCTGACATGTCTCCTTCTGAGAACAGAGCGTCTGTCTCCAA  
GCCAGCCGCTGTCATTCCAGTGCCTAAAGCTCCAGCCACACTCACCTGGGAGGATGAGAAGCAGCAGAG  
ATGGGGACAGCCACATGGTCCTGTCTCCTCTCCTTGGTGATCACAGGTGTCTGGTGCCCTTCAGG  
GACTTGAATCCAAGCTCGGAAGTCAACACCGCAAACCTACTCGAATCTCCTTCAAGTCTCTTGCTCACTT  
CCTGTTACATCTGTAGTACTTCTCTTTTACATATTAGGAGAAAAGAGGTGTCATTCTTTGAAAAGACT  
AAGGTACAGTGTATGTTGCAAGGTGTGCCAAACTTTTGTGCATGTGAAAAGAAAATGTGTCCGGGACT  
GGTCAGGTTTGCACGGGAGTGACGTAGGAGCCAAGGAGCAGGAGGAGCCAGGAGGAACCAAGCATTGA  
GGAGCTGCGGGATCTACTGCCTGGAGGAGAGGACGGACAAGGCCAGCCATGAGGAGTGACGGGAGAGGAG  
CACTTTGGGAAGACCACAGTGCACAGGACTCACACCCAGCCTTGCTGGAGAAAGCCCATGCCCCAGGCTT  
CTTCCAGGCAGCCCACTGTAAGACATCTCATTGCCTCTTCTGCCCTGGGCTGTCTGACCCACTGCCCC  
TTCCCCAGGAACACTTCCACTGGGATCA

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >RG217893 representing NM\_152460  
**Red**=Cloning site **Green**=Tags(s)

MDELALSFSLTCLLPENRASLSPSQPLSFQCLKAPATLTWEDEKQQRWGQPHGPVSSPLLGDHRCCLVPFR  
DLNPSSEVNTANLLESPSLLLTSYICSYFSFYILGEKRCHSLKRLRYSVCKVCPNFCACGKENVSGT  
GQVCTGVHVGAKQEPEGGTQALRSCGIYCLEERTDKASHEECRERSTLGRPQCTGLTPSLAGESPCPRL  
LPGSPTVRHLIASSCPGLSDPLPLPPGTLPGLS


**TRTRPLE** - GFP Tag - V



[View online »](#)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_152460

**ORF Size:** 729 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_152460.2](#), [NP\\_689673.2](#)

**RefSeq Size:** 2525 bp

**RefSeq ORF:** 732 bp


**Locus ID:** 146723

**UniProt ID:** [Q96MU5](#)

Cytogenetics: 17q25.1

Protein Families: Secreted Protein

### Product images:



Circular map for RG217893