

## Product datasheet for **RG217911**

### SEC23B (NM\_006363) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SEC23B (NM_006363) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SEC23B
Synonyms:	CDA-II; CDAIL; CDAN2; CWS7; HEMPAS; hSec23B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide  
Sequence:**

>RG217911 representing NM\_006363  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCGACATACCTGGAGTTCATCCAGCAGAATGAAGAACGGGATGGTGTGCGTTTTAGTTGGAACGTGT  
 GGCTTCCAGCCGCTGGAGGCTACAAGAATGGTTGTACCCTGGCTTGCTCCTTACTCCTTTGAAAGA  
 ACGTCCAGACCTACCTCCTGTACAATATGAACCTGTGCTTTCAGCAGGCCAACTTGTAAAGCTGTTCTC  
 AACCCACTTTGTGAGTTGATTACGAGCAAACTTTGGGCTGTAAATTTCTGTTTTCAAAGAAATCAGT  
 TTCCTCCAGCTTATGGAGGCATATCTGAGGTGAATCAACCTGCCGAATTGATGCCCCAGTTTTCTACAAT  
 TGAGTACGTGATACAGCGAGGTGCTCAGTCCCCTCTGATCTTTCTCTATGTGGTTGACACATGCCTGGAG  
 GAAGATGACCTTCAAGCACTCAAAGAGTCCCTGCAGATGTCCCTGAGTCTTCTCCTCCAGATGCTCTGG  
 TGGGTCTGATCACATTTGGAAGGATGGTGCAGGTTTCATGAGCTAAGCTGTGAAGGAATCTCCAAAAGTTA  
 TGTCTTCCGAGGGACCAAGGATTTAACTGCAAAGCAAATACAGGATATGTTGGGCTGACCAAGCCAGCC  
 ATGCCCATGCAGCAAGCAGCACCTGCACAACCACAGGAGCACCTTTTGCTTCAAGCAGATTTCTGCAGC  
 CTGTTCAAGATTGATATGAACCTCACTGATCTTCTTGGGGAGCTACAGAGGGACCCATGGCCAGTAAC  
 TCAGGGGAAGAGACCTTTGCGATCCACTGGTGTGGCTTTGTCCATTGCTGTTGGCTTGTGGAGGGCACT  
 TTTCAAACACAGGAGCCAGGATCATGTGTTTACTGGAGGTCCCCTACCCAAGGGCCTGGCATGGTGG  
 TTGGAGATGAATTAAGATTCCATTCTGTTCTTGGCATGATATTGAGAAAAGATAATGCACGATTCATGAA  
 AAAGGCAACCAAGCACTATGAGATGCTTGCTAATCGAACAGCTGCAAATGGTCACTGCATTGATATTTAT  
 GCTTGTGCCCTTGATCAAACCTGGACTTTTGGAGATGAAGTGTGTGCAAATCTTACTGGAGGCTACATGG  
 TAATGGGAGATTCTTCAACACTTCTCTTCAAGCAGACATTCCAAAGAATCTTACTAAAGATTTTAA  
 TGGAGATTTCCGAATGGCATTGGTGTCTCTTGGACGTAAGACCTCTCGGGAACCTGAAGATTGCAGGA  
 GCCATTGGTCCATGCGTATCTCTGAATGTGAAAGGACCGTGTGTGTCAGAAAATGAGCTTGGTGTGGTG  
 GCACGAGTCAGTGGAAAATCTGTGGCCTAGATCCTACATCTACACTTGGCATCTATTTTGAAGTTGTCAA  
 TCAGCACAAACACCCCGATCCCCAAGGAGGAGGAGGCCATCCAGTTTGTACGCATTATCAGCACTCC  
 AGCACCCAGAGACGCATCCGCGTGACCACCATCGCCCGAAATGGGCAGATGTACAGAGTCAGCTCAGGC  
 ACATAGAAGCAGCATTGACCAGGAGGCTGCGGCAGTGTGATGGCACGGCTTGGGGTGTCCGAGCGGA  
 GTCAGAGGAGGGGCCGATGTGCTCCGGTGGCTGGACCGACAACCTCATCCGACTGTGTCAAAGTTTGA  
 CAGTATAACAAAGAAGACCCCACTCTTTTAGGTTATCAGATTCCTTTTCTATATCCTCAGTTTATGT  
 TCCATCTGAGAAGATCTCCATTTCTTCAAGTGTAAACAACAGTCTGATGAGTCGTATATTACAGACA  
 TCATTTTGCCCGGAGGACCTGACCCAGTCCCTCATCATGATCCAGCCATTCTCTACTCTTACTCCTTT  
 CATGGGCCACAGAGCCAGTACTCTTGGATAGCAGCAGCATTCTAGCTGACAGAATTTTGTGATGGATA  
 CTTTCTTTCAAATGTCAATTTATCTTGGTGGAGCCATAGCCAGTGGCGTAAAGCTGGCTACCAGGCAT  
 GCCCGAGTATGAAAACCTCAAGCACCTTCTGCAGGCACCACTGGATGATGCTCAAGAAATCTGCAAGCA  
 CGCTTCCCGATGCCACGTTACATCAACACGGAGCATGGAGGCAGTCAGGCTCGATTCTTTTGTCCAAAG  
 TGAACCCATCTCAGACACACAATAACCTGTATGCTTGGGACAGGAACTGGAGCACCCATCCTAACTGA  
 TGATGTTAGCCTGCAGGTGTTTATGGACCATTTGAAGAAGCTGGCTGTCTCCAGTGCCTGT

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG217911 representing NM\_006363  
 Red=Cloning site Green=Tags(s)


MATYLEFIQQNEERDGVRFSSWNVWVSSRLEATRMVVPLACLITPLKERPDLPVQYEPVLCRSRPTCKAVL  
 NPLCQVDYRAKLWACNFCFQRNQFPAYGGISEVNQPAELMPQFSTIEYVIQRGAQSPLIFLYVVDTCLE  
 EDDLQALKESLQMSLSLLPPDALVGLITFGRMVQVHELSCGEGISKSYVFRGKDLTAKQIQDMLGLTKPA  
 MPMQQAQPAQPQEHFFASSRFLQPVHKIDMNLTDLLGELQRDPWPVTQGKRPLRSTGVALSIAVGLLEGT  
 FPNTGARIMLFTGGPPTQGPVMVGDDELKIPIRSWHDIKDNARFMKKATKHYEMLANRTAANGHCIDIY  
 ACALDQTGLLEMCCANLTGGYVMGDSFNTSLFKQTFQRIFTKDFNGDFRMAFGATLDVKTSTRELKIAG  
 AIGPCVSLNVKGPCVSENELGVGGTSQWKICGLDPTSTLGIYFEVVNQHNTPIPGGGRGAIQFVTHYQHS  
 STQRRIRVTTIARNWADVQSQLRHIEAAFDQEAALVLMARLGVFRAESEEGPDVLRWLDRLIRLCQKFG  
 QYNKEDPTSFRLSDSFLYPQFMFHLRRSPFLQVFNNSPDESSYRHHFARQDLTQSLIMIPIILYSYSF  
 HGPPPEVLLDSSSILADRILLMDTFFQIVIYLGETIAQWRKAGYQDMPEYENFKHLLQAPLDDAQEILQA  
 RFPMPRYINTEHGGSQARFLLSKVNPSQTHNNLYAWGQETGAPILTDDVSLQVFMHLLKLAIVSSAC

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



## Plasmid Map:



ACCN: NM\_006363

ORF Size: 2301 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_006363.6](#)

RefSeq Size: 3429 bp

RefSeq ORF: 2304 bp

Locus ID: 10483

UniProt ID: [Q15437](#)

Cytogenetics: 20p11.23

Domains: zf-Sec23\_Sec24, Sec23\_trunk, Sec23\_helical, Gelsolin

**Gene Summary:** The protein encoded by this gene is a member of the SEC23 subfamily of the SEC23/SEC24 family, which is involved in vesicle trafficking. The encoded protein has similarity to yeast Sec23p component of COPII. COPII is the coat protein complex responsible for vesicle budding from the ER. The function of this gene product has been implicated in cargo selection and concentration. Multiple alternatively spliced transcript variants have been identified in this gene. [provided by RefSeq, Feb 2010]