

Product datasheet for **RG217913**

5HT4 Receptor (HTR4) (NM_001040171) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	5HT4 Receptor (HTR4) (NM_001040171) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	5HT4 Receptor
Synonyms:	5-HT4; 5-HT4R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG217913 representing NM_001040171 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACAACTTGATGCTAATGTGAGTTCTGAGGAGGTTTCGGGTCAGTGGAGAAGGTGGTGCTGCTCA
CGTTTCTCTCGACGGTTATCCTGATGGCCATCTGGGGAACCTGCTGGTATGGTGGCTGTGTGCTGGGA
CAGGCAGCTCAGGAAAAATAAAACAAATTATTTTCATTGTATCTCTTGCTTTTGGGATCTGCTGGTTTCG
GTGCTGGTATGCCCTTTGGTGCCATTGAGCTGGTTCAAGACATCTGGATTTATGGGAGGTGTTTTGTC
TTGTTCCGACATCTCTGGACGCTCTGCTACAACGGCATCGATTTTTCACCTGTGCTGCATTCTCTGGA
TAGGTATTACGCCATCTGCTGCCAGCCTTTGGTCTATAGGAACAAGATGACCCCTCTGCGCATCGCATT
ATGCTGGGAGGCTGCTGGGTCATCCCCACGTTTATTTCTTTTCTCCCTATAATGCAAGGCTGGAATAACA
TTGGCATAAATGATTTGATAGAAAAGAGGAAGTTCAACCAGAAGCTAACTCTACGTACTGTGCTTTCAT
GGTCAACAAGCCCTACGCCATCACCTGCTCTGTGGTGGCCTTCTACATCCCATTCTCTCATGGTGCTG
GCCTATTACCGCATCTATGTCACAGCTAAGGAGCATGCCATCAGATCCAGATGTTACAACGGGCAGGAG
CCTCCTCCGAGAGCAGGCCTCAGTCGGCAGACCAGCATAGCACTCATCGCATGAGGACAGAGACCAAAGC
AGCCAAGACCCTGTGCATCATCATGGGTTGCTTCTGCCTCTGCTGGGCACCATTCTTTGTACCAATATT
GTGGATCCTTTTATAGACTACACTGTCCCTGGCAGGTGTGGACTGCTTTCTCTGGCTCGGCTATATCA
ATTCCGGTTGAACCTTTTCTCTACGCCTTCTTGAATAAGTCTTTTAGACGTGCCTTCTCATCATCCT
CTGCTGTGATGATGAGCGCTACCGAAGACCTTCCATTCTGGGCCAGACTGTCCCTTGTTCACCAACAACC
ATTAATGGATCCACATGTACTAAGATTT

ACCGTACCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG217913 representing NM_001040171
 Red=Cloning site Green=Tags(s)

MDKLDANVSSEEGFGSVEKVVLLTFLSTVILMAILGNLLVMVAVCWDRQLRKIKTNYFIVSLAFADLLVS
 VLVMPFGAIELVQDIWIYGEVFLVRTSLDLVLLTTASIFHLCCISLDRYYAICQPLVYRNKMTPLRIAL
 MLGGCWVPTFFISFLPIMQGWNNIGIIDLIEKRKFNQNSNSTYCVFMVKNPYAITCSVAFYIPFLMLVL
 AYYRIYVTAKEHAHQIQLMRAGASSESRPQSADQHSRMRTEKAAKTLCIIMGCFCLCWAPFFVTNI
 VDPFIDYTVPGQVWTAFLWLGYINSGLNPFYAFLNKSFRRRAFLIILCCDDERYRRPSILGQTVPCSTTT
 INGSTHVLRF

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001040171

ORF Size: 1080 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001040171.1](#), [NP_001035261.1](#)

RefSeq Size: 1355 bp

RefSeq ORF: 1082 bp

Locus ID: 3360


Cytogenetics: 5q32

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Calcium signaling pathway, Neuroactive ligand-receptor interaction

Gene Summary: This gene is a member of the family of serotonin receptors, which are G protein coupled receptors that stimulate cAMP production in response to serotonin (5-hydroxytryptamine). The gene product is a glycosylated transmembrane protein that functions in both the peripheral and central nervous system to modulate the release of various neurotransmitters. Multiple transcript variants encoding proteins with distinct C-terminal sequences have been described. [provided by RefSeq, May 2010]

Product images:



Circular map for RG217913