

Product datasheet for **RG218145**

CAMK2B (NM_172078) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CAMK2B (NM_172078) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CAMK2B
Synonyms:	CAM2; CAMK2; CAMKB; CaMKIIbeta; MRD54
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG218145 representing NM_172078
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCACCACGGTGACCTGCACCCGCTTCACCGACGAGTACCAGCTCTACGAGGATATTGGCAAGGGG
 CTTTCTCTGTGGTCCGACGCTGTGTCAAGCTCTGCACCGCCATGAGTATGCAGCCAAGATCATCAACAC
 CAAGAAGCTGTCAGCCAGAGATCACCAGAAGCTGGAGAGAGAGGCTCGGATCTGCCGCCCTCTGAAGCAT
 TCCAACATCGTGCCTCCACGACAGCATCTCCGAGGAGGGCTTCCACTACCTGGTCTTCGATCTGGTCA
 CTGGTGGGAGCTCTTTGAAGACATTGTGGCGAGAGAGTACTACAGCGAGGCTGATGCCAGTCACTGTAT
 CCAGCAGATCCTGGAGGCCGTTCTCCATTGTACCAAATGGGGTCTCCACAGAGACCTCAAGCCGGAG
 AACCTGCTTCTGGCCAGCAAGTCAAAGGGGCTGCAGTGAAGCTGGCAGACTTCGGCCTAGCTATCGAGG
 TGCAGGGGACCAGCAGGCATGGTTTGGTTTCGCTGGCACACCAGGCTACCTGTCCCTGAGGTCCTTCG
 CAAAGAGGCGTACGGCAAGCCCGTGACATCTGGGCATGTGGGTGATCCTGTACATCCTGCTCGTGGGC
 TACCCACCTTCTGGGACGAGGACCAGCACAAGCTGTACCAGCAGATCAAGGCTGGTGCCTATGACTTCC
 CGTCCCTGAGTGGGACACCGTCACTCCTGAAGCCAAAAACCTCATCAACCAGATGCTGACCATCAACCC
 TGCCAAGCGCATCACAGCCATGAGGCCCTGAAGCACCCGTGGTCTGCCAACGCTCCACGGTAGCATCC
 ATGATGCACAGACAGGAGACTGTGGAGTGTCTGAAAAAGTTCAATGCCAGGAGAAAGCTCAAGGGAGCCA
 TCCTCACCACCATGCTGGCCACACGGAATTTCTCAGCCAAGAGTTTACTCAACAAGAAAGCAGATGGAGT
 CAAGCCCCAGACGAATAGCACAAAAACAGTGCAGCCGCCACCAGCCCCAAAGGGACGCTTCTCCTGCC
 GCCCTGGAGCCTCAAACCACCGTATCCATAACCCAGTGGACGGGATTAAGGAGTCTTCTGACAGTGCCA
 ATACCACCATAGAGGATGAAGACGCTAAAGCCCCGGAAGCAGGAGATCATTAAAGACCACGGACAGCTCAT
 CGAGGCCGTCAACAACGGTGACTTTGAGGCTACGCGAAAAATCTGTGACCCAGGGCTGACCTCGTTTGAG
 CCTGAAGCACTGGCAACCTGGTTGAAGGGATGGACTTCCACAGATTCTACTTCGAGAACCTGCTGGCCA
 AGAACAGCAAGCCGATCCACACGACCATCCTGAACCCACACGTCACGTCATTGGAGAGGATGCCGCTG
 CATCGTTACATCCGGCTCACGCAGTACATTGACGGGCAGGGCCGCCCCGACCAGCCAGTCTGAGGAG
 ACCCGCGTGTGGCACCGCCGCGACGGCAAGTGGCAGAACGTCACCTTCCACTGCTCGGGCGCCTGTGG
 CCCCCTGCAG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG218145 representing NM_172078
 Red=Cloning site Green=Tags(s)

MATTVTCTRFDTDEYQLYEDIGKGAFSVVRRCVKLCTGHEYAAKIINTKKLSARDHQKLEREARICRLLKH
 SNIVRLHDSISEEGFHYLVFDLVTGGELFEDIVAREYYSEADASHCIQQILEAVLHCHQMGGVVRDLKPE
 NLLLASKCKGAAVKLADFLAIEVQGDQAWFGFAGTPGYLSPEVLRKEAYGKPVDIWACGVILYILLVG
 YPPFWDEDQHKLYQIKAGAYDFPSPEWDTVTPEAKNLINQMLTINPAKRITAHEALKHPWVCQRSTVAS
 MMHRQETVECLKKFNARRKLGAILTTMLATRNFSAKSLLNKKADGVKPKQTNSTKNSAAATSPKGLTTPA
 ALEPQTTVIHNPVDGIKESSDSANTTIEDEDAKARKQEIKTTEQLIEAVNNGDFEAYAKICDPGLTSFE
 PEALGNLVEGMDFHRFYFENLLAKNSKPIHTTILNPHVHVIGEDAACIAYIRLTQYIDGQGRPRTSQSEE
 TRVWHRDGGKQNVVHFHCSGAPVAPLQ

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_172078

ORF Size: 1551 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_172078.1](#), [NP_742075.1](#)

RefSeq Size: 2066 bp

RefSeq ORF: 1629 bp


Locus ID: 816

UniProt ID: [Q13554](#)

Cytogenetics: 7p13

Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Calcium signaling pathway, ErbB signaling pathway, Glioma, GnRH signaling pathway, Long-term potentiation, Melanogenesis, Neurotrophin signaling pathway, Olfactory transduction, Oocyte meiosis, Wnt signaling pathway
Gene Summary:	The product of this gene belongs to the serine/threonine protein kinase family and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. Calcium signaling is crucial for several aspects of plasticity at glutamatergic synapses. In mammalian cells, the enzyme is composed of four different chains: alpha, beta, gamma, and delta. The product of this gene is a beta chain. It is possible that distinct isoforms of this chain have different cellular localizations and interact differently with calmodulin. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014]

Product images:



Circular map for RG218145