

Product datasheet for **RG218232**

SYPL1 (NM_006754) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SYPL1 (NM_006754) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SYPL1
Synonyms:	H-SP1; SYPL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG218232 representing NM_006754 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCGCCCAACATCTACTTGGTTCGCCAGCGGATCAGTCGACTCGGCCAGAGGATGTCCGGCTTCCAGA
 TCAACCTCAACCCGCTCAAGGAGCCACTCGGCTTCATCAAGGTCCTCGAGTGGATTGCTTCTATCTTTGC
 TTTTGCCACCTGTGGAGGTTTTAAGGGCCAAACAGAAATTCAGTGAATTGCTCCTGCAGTTACTGAG
 AATAAACTGTTACAGCTACTTTTGGTTATCCATTAGGTTGAATGAGGCATCATTCAGCCACCTCCAG
 GTGTAACATATGTGATGTAATTGGAAGATTACGTCCTCATAGGCGATTACTTCTTCTGCACAATT
 CTATGTTACCTTTGCAGTCTTTGTGTTCTGTACTGCATTGCTGCCCTTCTGCTTTATGTTGGCTACACG
 AGTCTGTATCTGGATAGTCGTAACCTTCATGATAGACTTTGTTGTACACTTGTGCCACTTTTTTGT
 GGTTGGTGAGCACTTCAGCCTGGGCTAAAGCTCTGACAGATATTTAAATAGCTACTGGTCACAATATTAT
 TGATGAACTTCGCCTTGTAAGAAGAAAGCAGTACTGTGTTACTTTGGCTCTGTGACCAGTATGGGATCC
 CTAAATGTATCTGTGATATTTGGCTTTCTAAATATGATACTCTGGGGAGGAAATGCTTGGTTTGTGTACA
 AGGAGACCAGCTACACAGTCCATCAAATACATCTGCCCTCATAGCCAAGGAGGTATCCACCTCCTAC
 CGGAATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >RG218232 representing NM_006754
 Red=Cloning site Green=Tags(s)

MAPNIYLVQRISRLGQRMSGFQINLNPLKEPLGFIKVLWIASIFAFATCGGFKGQTEIQVNCPPAVTE
 NKTVTATFGYPFRLNEASFQPPPGVNICDVNWKDYVLIGDYSSSAQFYVTFVAVFVFLYCIAALLLYVGYT
 SLYLDSRKLPIDFVVTLVATFLWL VSTSAWAKALTDIKIATGHNIIDELPPCKKKAVLCYFGSVTSMGS
 LNVSVIFGFLNMILWGGNAWFVYKETSLHSPSNTSAPHSQGGIPPPTGI

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_006754

ORF Size: 777 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_006754.5](#)

RefSeq Size: 2335 bp

RefSeq ORF: 780 bp

Locus ID: 6856


UniProt ID: [Q16563](#)

Cytogenetics: 7q22.3

Domains: Synaptophysin

Protein Families: Transmembrane

Product images:



Circular map for RG218232