

Product datasheet for **RG218417**

ZNF626 (NM_001076675) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF626
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG218417 representing NM_001076675
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGACCATTGCAATTTAGAGATGTGCCATAGAATTCTCTCTGGAGGAGTGGCATTGCCTGGACACTG
 CACAGCGGAATTTATATAGGAATGTGATGTTAGAGAACTACAGTAACCTGGTCTTCCTTGGTATTACTGT
 TTCTAAGCCAGACCTGATCACCTGTCTGGAGCAAGGAAGAAAACCTTTGACCATGAAGAGAAATGAGATG
 ATAGCCAAACCTCAGTAATGTGTTCTATTTTGCCTAAGACCTTTGGCCAGAGCAGAGCATGAAGATT
 CTTTCCAAAAAGTGGTACTGAGAAGATGAAAAATGTGAACATGACAATTTACAGTTAAAAAAGGATG
 TATAAGTGTGGATGAGTGAAGGTGCACAAAGAAGTTATAATGAACCTTAACCAATGTTTGACAACCTACC
 CCAAGAAAAATATGTCAATGTGATAAATATGTGAAAGTCCTTCATCAATTTCCAAATCAAACGGACAAA
 AGAGAGGACATACTGGGAAAAACCTTTCAAATATATAGAATGTGGCAAAGCTTTTAAGCAGTTCTCAAC
 TCTTACTACACATAAGAAAATTCATACTGGAGGGAACCTTACAAATGTGAAGAATGTGGCAAAGCCTTT
 AACCCTCTTGTAGCCTTACTAGACATAAGAAAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAAT
 GTGGCAAAGCCTTTAAGCACTCCTCCACTTACTACACATAAGAGAAATCATACTGGAGAGAAACCTTA
 CAAGTGTGATAAATGTGGCAAAGCCTTTATGTCATCCTCAACCTTAGTAAACATGAGATAATTCATACG
 GAAAAGAAACCTTACAAATGTGAAGAATGTGGCAAAGCCTTCAACCGGTCTCAACCTTACTACACATA
 AGATAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAATGTGACAAAGCCTTTAAGTACTCCTATAC
 CCTTACTACACATAAAAGAATTCATACTGAAGACAAACCTTACAAATGTGAAGAATGTGGCAAAGCCTTT
 AAGTACTCCTTACCCTTACTACACATAAGAGAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAAT
 GTGGCAAAGCCTTCAAGCGGTCTTCAAGCCTTACTACGCATAAGATAATTCATACTGGAGAGAAACCTTA
 CAAATGTGAAGAATGTGGCAAAGCCTTTAAGTACTCCTTAACTTACTACACATAAGAAAATTCATACT
 GGAGAGAGACCTTACAAATGTGAAGAATGTGGCAAAGCCTTCAACCAAGTCTCAATCCTTACTACACATA
 GGAGAATTCATACTGGAGAGAAATTCATAATGTGAAGAATGTGGCAAAGCCTTTAAGTGTCTCTCTAA
 CCTTACTACACATAAGAAAATTCATACTGGAGAGAGACCTTACAAATGTGAAGAATGTGGCAAAGCCTTC
 AACGAGTCTCAATCCTTACTACACATGAGAGAATCATACTGGAGAGAAATTCATAATGTGAAGAATG
 TGGCAAAGCCTTCAAGCGGTCTTACACCTTACTACACATAAGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG218417 representing NM_001076675
 Red=Cloning site Green=Tags(s)

MGPLQFRDVAIEFSLEEWHLDTAQRNL YRNVMLNYSNLVFLGITVSKPDLITCLEQGRKPLTMKRNEM
 IAKPSVMCSHFAQDLWPEQSMKDSFQKVVLRRYEKCEHDNLQLKKGCISVDECKVHKEGYNELNQCLTTT
 PRKICQCDKYVKVLHQFPNSNGQKRGHTGKKPFKYIECGKAFKQFSTLTTHKKIHTGGKPYKCEECGKAF
 NHSCSLTRHKIHTGEKPYKCEECGKAFKHSSTLTTHKRNTGEKPYKCDKCGKAFMSSSTLSKHEIHT
 EKKPYKCEECGKAFNRSSLTTHKKIHTGEKPYKCEECDFKYSYTLTTHKRIHTEDKPYKCEECGKAF
 KYSSLTTHKRIHTGEKPYKCEECGKAFKRSSDLTTHKKIHTGEKPYKCEECGKAFKYSSNLTHKKIHT
 GERPYKCEECGKAFNQSSILTTHRRHTGEKPYKCEECGKAFKCSSNLTHKKIHTGERPYKCEECGKAF
 NQSSILTTHERRIILERNSTNVKNVAKPSSGPHLLHIR

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001076675

ORF Size: 1584 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001076675.3](#)

RefSeq Size: 6094 bp

RefSeq ORF:	1587 bp
Locus ID:	199777
UniProt ID:	Q68DY1
Cytogenetics:	19p12
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for RG218417