

Product datasheet for **RG218455**

KRTAP12-1 (NM_181686) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: TurboGFP

Symbol: KRTAP12-1

Synonyms: KAP12.1; KRTAP12.1

Mammalian Cell Selection: Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG218455 representing NM_181686
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGTGCCACACCAGCTGCTCCTCGGGCTGCCAGCCAGCCTGCTGCGGCCAGCCCCTGCCAGGCATCCT
GTTACATCCCCGTGGGCTGCCAGTCCTCCGTGTGCGTGCCCGTGAGCTTCAAGCCAGCCGTGTGTGTGCC
CGTGAGATGCCAGTCCTCTGTGTGCGTGCCCGTGAGCTGCAGGCCCGTCGTGTATGCGGCTCCCTCCTGC
CAGTCCTCTGGGTGCTGCCAGCCTTCTGCACCAGCGTCCTCTGCAGACCATCTCCTGCAGACCCCTT
CCTGCTGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG218455 representing NM_181686
Red=Cloning site Green=Tags(s)


MCHTSCSSGCQPACCAPSPCQASCYIPVGCQSSVCVPVSFKPAVCVPVRCQSSVCVPVSCRPPVYAAPSC
QSSGCCQPSCTSVLCRPISCSTPSCC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_181686

ORF Size: 288 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_181686.2
RefSeq Size:	604 bp
RefSeq ORF:	291 bp
Locus ID:	353332
UniProt ID:	P59990
Cytogenetics:	21q22.3
Gene Summary:	In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG218455