

## Product datasheet for **RG218488**

### Proline rich membrane anchor 1 (PRIMA1) (NM\_178013) Human Tagged ORF Clone

#### Product data:

|                    |                                |
|--------------------|--------------------------------|
| Product Type:      | Expression Plasmids            |
| Tag:               | TurboGFP                       |
| Symbol:            | Proline rich membrane anchor 1 |
| Synonyms:          | PRIMA                          |
| Mammalian Cell     | Neomycin                       |
| Selection:         |                                |
| Vector:            | pCMV6-AC-GFP (PS100010)        |
| E. coli Selection: | Ampicillin (100 ug/mL)         |

**ORF Nucleotide Sequence:** >RG218488 representing NM\_178013  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGCTCCTCCGGGACTTGGTGCTGCGCCGTGGCTGCTGCTGGTCTCGCTGCTGCTGCACTGCGCGCTCC  
ACCCGCTCTGGGGCTTCGTGCAGGTGACGCATGGTGAGCCCCAGAAGTCTGCTCCAAAGTACTGACAG  
CTGCCGACAGCTGCCAGTGCCGGCCCCCTCCCCGCTGCCCCGCGCCCCACCCCGCCACCTCCC  
AGACTCCTCTCCGCCCCAGTCCCAACTCTACCTCTTGCCCCACTGAGGAAAGCTGGTGGTGGGGCTGG  
TGATCATCATTGCCGATGCTGTGCCTCCCTGGTGTCTGACTGTGCTTGTGCATCATTGTCTACAAAGC  
CATAAAAAGGAAACCACTGAGAAAAGACGAAAATGGCACCAGCGTTGCTGAGTATCCCATGAGTGCTTCG  
CAGAGCAACAAAGGAGTAGACGTGAACAACGAGTGGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG218488 representing NM\_178013  
Red=Cloning site Green=Tags(s)


MLLRDLVLRGCCWSSLLHLCALHPLWGFVQVTHGEPQKSKSVTDSCRHVCQCRPPPLPPPPPPPP  
RLLSAPAPNSTSCPTEESWWSGLVIIIVACASLVFLTVLVIICYKAIKRKPLRKDENGTSVAEYPMAS  
QSNKGVDVNNNAV

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI



## Cloning Scheme:



**ACCN:** NM\_178013

**ORF Size:** 459 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_178013.2](#)

**RefSeq Size:** 910 bp

**RefSeq ORF:** 462 bp

**Locus ID:** 145270

**UniProt ID:** [Q86XR5](#)

**Cytogenetics:** 14q32.12

**Protein Families:** Transmembrane

**Gene Summary:** The product of this gene functions to organize acetylcholinesterase (AChE) into tetramers, and to anchor AChE at neural cell membranes. [provided by RefSeq, Nov 2008]