

Product datasheet for **RG218998**

TTC14 (NM_001042601) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TTC14
Synonyms:	DRDL5813; PRO19630
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG218998 representing NM_001042601
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACCGGGACCTTTTGCAGCTCGCTAAATTGCCACGGGTCGTCTTTGCTCTCTACTTCGGAGCG
AACAGCAGGACAATCCACACTTCCTGAGCCTCCTGGGGTCGGCCGCCGAGCCAGCCCGGGGCCCGCCGCC
CCAGCACCCGTTGCAGGCGAGAAAAGAGAGAGTTGACAACATCGAGATACAGAAATTCATCTCCAAA
AAAGCGGATCTGCTTTTGCACCTTCTGGAAATCAGATGCACCTGCAACTTCTGAAATTAATGAAGACA
GTGAAGATCATTATGCAATCATGCCACCTTTAGAGCAATTCATGGAGATACCTAGTATGGATCGGAGAGA
GCTGTTTTCCGAGATATTGAGCGTGGTGATATAGTGATTGGAAGAATTAGTCTATTTCGGGAATTCGGT
TTTTTCATGGTGTTGATCTGTTAGGAAGTGGTATCATGAGAGATATAGCCCACTTAGAAATCACAGCTC
TTTGTCCTTAAGAGATGTGCCTTCTCACAGTAACCATGGGGATCCTTTATCATATTACCAAACTGGTGA
CATATTTCGAGCTGGAATCAAGGATATTGACAGATACCATGAAAAGCTAGCAGTATCTCTGTATAGCTCT
TCTCTTCCACCACCTATCTGGTATTAATTAAGGTGAATTAGCTCTGAAGAGCTTCCTTTATACTACA
GGAGAAGTGTGAGCTAAATAGCAATTCCTTGGAGTCCTATGAAAATGTCATGCAGAGTTCTTGGGATT
TGTTAATCCAGGAGTAGTTGAATTCCTCTAGAAAACTAGGAATAGATGAATCTAATCCACCATCTTTA
ATGAGAGGCCCTACAAAGCAAAATTTCTCTGAAGATGATTTTGCTTTCGATTGAGAAAAAACAATCCG
CATCTTGGGCTTTAAATGTGTGAAGATCGGAGTTGACTATTTTAAAGTTGGACGCCATGTGGATGCTAT
GAATGAATAACAATAAGCTTTGGAATAGACAAACAAACGTGGAAGCTTTGGTAGCTCGTGAGCATT
TATGCGACAAAAGGAAGTTTGAACAAAGCAATAGAAGATTTTGGCTTGCATTAGAAAAGTCCAACTC
ACAGAAATGCAAGAAAATACCTCTGCCAGACACTTGTAGAGAGAGGAGGACAGTTAGAAGAAGAAGAAA
GTTTTTAAATGCTGAAAGTTACTATAAGAAAGCTTTGGCTTTGGATGAGACTTTTAAAGATGCAGAGGAT
GCTTTGCAGAACTTCATAAATATATGCAGGTGATTCCTTATTTCTCTTAGAAATT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG218998 representing NM_001042601

Red=Cloning site Green=Tags(s)


MDRDLRLQSLNCHGSSLLSLLRSEQQDNPHFRSLLGSAAPARGPPPQHPLQGRKEKRVDNIEIQKFISK
KADLLFALSWKSDAPATSEINEDSEDHYAIMPPLQFMEIPSMRRELFFRDIERGDIVIGRISSIREFG
FFMVLICLGSIMRDIAHLEITALCLRDVPSHSHNGDPLSYYQTGDIIRAGIKDIDRYHEKLAVSLYSS
SLPPHLSGIKLGVISSEELPLYRRSVELNSLSLESYENVMQSSLGFVNPVGVEFLLEKLGIDESNPPSL
MRGLQSKNFSEDDFASALRKKQSASWALKCVKIGVDYFKVGRHVDAMNEYNKALEIDKQNVREALVARGAL
YATKGSLSNKAIEDFELALENCPTHNRKYLQCTLVERGGQLEEEKFLNAESYKKALALDETFKDAED
ALQKLHKYMQVIPYFLEI

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001042601

ORF Size:

1317 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001042601.3

RefSeq Size: 4758 bp


RefSeq ORF: 1320 bp

Locus ID: 151613

UniProt ID: Q96N46

Cytogenetics: 3q26.33

Product images:



Circular map for RG218998