

Product datasheet for **RG219043**

KRTAP4-7 (NM_033061) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: KRTAP4-7
Synonyms: KAP4.7; KRTAP4.7
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG219043 representing NM_033061
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTCAGCTCCTGTTGTGGCTCCGTGTGCTCTGACCAGGGCTGCGGCCAAGTCCTCTGTCAGGAGACCT
GCTGCCGCCCCAGCTGCTGTCAGACCACCTGTTGCAGGACCACCTGCTACCGCCCCAGCTGTTGTGTGTC
CAGCTGCTGCAGGCCCCAGTGTGCCAGTCTGTGTGCTGCCAACCCACCTGCTGTCGCCCCAGCTGCTGT
GAGACGACCTGCTGCCACCCTAGGTGCTGCATCTCCAGCTGCTGCCGCCCCAGCTGCTGTATGTCCAGCT
GCTGCAAGCCCCAGTGTGCCAGTCTGTGTGCTGCCAGCCACCTGCTGCCGCCCCAGCTGCTGCATCTC
CAGCTGCTGTGCCCCAGCTGCTGTGTGCTCCAGGTGCTGCAGGCCCCAGTGTGCCAGTCTGTGTGCTGC
CAGCCAACCTGCTGCCGTCCAGCTGCTGCATCTCCAGCTGCTGCCGCCCTCTTGTGTGAATCCAGCT
GCTGCCGCCCATGCTGCTGCCGCCCTGCTGCTGCCTGCGTCCAGTCTGTGGCCGAGTCTCCTGCCACAC
CACTTGCTATCGCCCAACCTGTGTCATCTCCACCTGTCCCCGCCCTTGTGCTGTGCCTCCTCTTGCTGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG219043 representing NM_033061
Red=Cloning site Green=Tags(s)

MVSSCCGSVCSDQCGQVLCQETCCRPSCCQTTCRRTTCYRPSCCVSSCCRPQCCQSVCCQPTCCRPSCC
ETTCCHPRCCISSCCRPSCCMSSCKPQCCQSVCCQPTCCRPSCCISCCRPSCCVSRCCRPQCCQSVCC
QPTCCRPSCCISCCRPSCCESSCCRPCCCRPCCCLRPVCGRVSCHTTCYRPTCVISTCPRLCCASSCC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN:	NM_033061
ORF Size:	630 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_033061.2 , NP_149050.2
RefSeq Size:	1202 bp


RefSeq ORF: 468 bp

Locus ID: 100132476

Cytogenetics: 17q21.2

Gene Summary: This protein is a member of the keratin-associated protein (KAP) family. The KAP proteins form a matrix of keratin intermediate filaments which contribute to the structure of hair fibers. KAP family members appear to have unique, family-specific amino- and carboxyl-terminal regions and are subdivided into three multi-gene families according to amino acid composition: the high sulfur, the ultrahigh sulfur, and the high tyrosine/glycine KAPs. This protein is a member of the ultrahigh sulfur KAP family and the gene is localized to a cluster of KAPs at 17q12-q21. [provided by RefSeq, Mar 2009]

Product images:



Circular map for RG219043