

## Product datasheet for **RG219078**

### **ALS2CR13 (FAM117B) (NM\_173511) Human Tagged ORF Clone**

#### **Product data:**

|                           |                                                       |
|---------------------------|-------------------------------------------------------|
| Product Type:             | Expression Plasmids                                   |
| Product Name:             | ALS2CR13 (FAM117B) (NM_173511) Human Tagged ORF Clone |
| Tag:                      | TurboGFP                                              |
| Symbol:                   | ALS2CR13                                              |
| Synonyms:                 | ALS2CR13                                              |
| Mammalian Cell Selection: | Neomycin                                              |
| Vector:                   | pCMV6-AC-GFP (PS100010)                               |
| E. coli Selection:        | Ampicillin (100 ug/mL)                                |



[View online »](#)

**ORF Nucleotide Sequence:**

>RG219078 representing NM\_173511  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTCCCAGCGGGTGAGGCGCAATGGGTCCCCACGCGCGCGGCTCCCTTGGGGGTGGTGCGGTGGCCA  
 CGGCCGGGGGACCGGGAGCCGCTTGACGCCCATGAGGGCGACGGTTCGGTTCAGCTGAAGCAGCAGCA  
 GCAGCAGCAACATGGCAGCCCCACGCGGAGCGGCGGCGGCGGCGGCAACAACAACGGTGCTGCTGT  
 GGTGGCGCCTCAGGCCCGCAGGCGGCGGCGGCGGCGGTGGCCGCGCACCGCCTCGCGCAGCACGAGCC  
 CCACGCGCGGCGGGAACGCGGCCGCGCACCAGCCCCACGGTGGCCACGCAGAGCGGCGCGTCCGC  
 GACGTCCACGCGAGGCACAGCCCCACGCGCAGCGCGCGCTGGAGCTCGCGGAGCCCCCACGGCCG  
 CCGCGCGCGCGCGCTGCTGGGACCGTGTCTCGCCAGCTCGTCGCCACCCACCTGTGGACCGGCG  
 AGGTGAGCGCGGCCCCACCCACGCCCGCTCCGGCATCGGAGGAGGTCTCCGGAGCAGAGCCGAAGCTC  
 GCCGGAAGAGGAGCCCCAGCGCCCCGTTTCAAAGCAGGTGACAAAACACGACAGCCTTCTTCAAGC  
 CCTCCAGTATTATCCGACGCACTTCTCCCTGGATACTCTTGCTGCACCGTATCTTGCTGGACACTGGC  
 CTCGGGATAGCCATGGGCAAGCTGCACCTTGCATGAGGGACAAAGCTACACAGACAGAGAGTGCATGGGC  
 TGAAGAATACTCTGAAAAGAAGAAAGGGTCTCACAAGCGCTCAGCATCTTGGGGCAGTACAGATCAACTT  
 AAGGAGATTGCAAAATTACGCCAGCAGTTGCAGAGAAGTAAACACAGCAGTCGGCATCATCGAGATAAAG  
 AAAGACAGTCTCCATTTTCATGGCAACCATGCAGCTATTAACAGTGTGAGGCTCTGTTCCTCAAGAGTGC  
 ACTTATTCCTGTAATCCCATCACCAATCAACAGGCTCCCGTTCCGGAATAGCGTGGAAGGATTGAAT  
 CAGGAGATTGAAATAATAATTAAGAGACTGGGAAAAGGAAGAGCAACTTATACCGCAAGATATCCAG  
 ATGGCCATCGTGCTCCACCCCCCTTGACAGAGAAGTAGCAGCACGCGCAGCATTGACACACAGACGCC  
 TGGTGGGCGAGACAGGGGAAGCAACAACAGCAGCGGTTCCAGTCCGTGTCCCAACATCGTTCCTCACC  
 ATTTCCAATGAAGGTAGCGAGGAGAGTCTTGCTCAGCGGATGACCTGCTTGTGATCCCAGAGATAAAG  
 AGAATGGGAACAATTCTCTTTGCCCAATATGCAACCTACCAAAACCTAAACAAGTTATATGTTCAA  
 AAGGGAACCTCTGAGGGCTGTGAAAGGTCAAAGTCTTTGAGGAATGCTCGCCAAAACAGCTTCATGAA  
 ATCCCAGCCTTTTATTGCTGACAAAACAAGGTGAATTTATTCTTAAAGTGGATCTGCTTTCTGCC  
 TCGTCAGCATCTCAAGCCACTCTTCTACCCGGATCTTACTCAAGGGCTCTGGCCACAGCCTGAC  
 AGTCACCACTGGCATGACAACCACTCTGCTGCAGCCTATTGCTGTGGCCTCCCTGTCTACAAACACAGAG  
 CAAGACCGAGTCTCTGAGGAACAAGTACAGTCATGCCATCAGCTTCTCTACTCCACCACCAGAACCA  
 TTGAGGAAGCTGAAGGA

**ACGCGTACGCGGCGGCTCGAG** – GFP Tag – GTTTAA

**Protein Sequence:**

>RG219078 representing NM\_173511  
 Red=Cloning site Green=Tags(s)


MSQVRVRNGSPTPAGSLGGGAVATAGGPGSRLQPMRATVPFQLKQQQQQHQSPTRSGGGGGNNNGGCC  
 GGASGPAGGGGGGPRTASRSTSPTRGGGNAARTSPTVATQTGASATSTRGTSPTRSAAPGARGSPPRP  
 PPPPPLLGTVSSPSSSPHLWTGEVSAAPPPARVRHRRRSPEQSRSSPEKRSPSAPVCKAGDKTRQPS  
 PSSIIIRRTSSDLTAAPYLGHWPDRSHQAAPCMRDKATQTESAWAEEYSEKKKGSHKRSASWGSTDQL  
 KEIAKLQQQLQSKHSSRHHRDKERQSPFHGNHAAINQCQAPVPKSALIPVIPITKSTGSRFRNSVEGLN  
 QEIEIIKETGEKEEQLIPQDIPDGHRAAPPLVQRSSSTRSIDTQTPGGADRGSNSSRSQSVSPTSFLT  
 ISNEGSEESPCSADDLLVDRDKENGNNSPKPYATSPKPNNSYMFKREPPEGCERVKVFECCSPKQLHE  
 IPAFYCPDKNKVNFIKSGSAFCLVSIKPLLPTPDLTLKSGHSLTVTTGMTTLLQPIAVASLSTNTE  
 QDRVSRGTSTVMPSASLLPPPEPIEEAEG

**TRTRPLE** – GFP Tag – V

**Restriction Sites:**

SgfI-MluI

## Cloning Scheme:



ACCN: NM\_173511

ORF Size: 1767 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM\_173511.4  
 RefSeq Size: 5795 bp  
 RefSeq ORF: 1770 bp  
 Locus ID: 150864  
 UniProt ID: Q6P1L5  
 Cytogenetics: 2q33.2

## Product images:



Circular map for RG219078