

Product datasheet for **RG219522**

PEX5 related protein (PEX5L) (NM_016559) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX5 related protein (PEX5L) (NM_016559) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PEX5 related protein
Synonyms:	PEX5R; PEX5RP; PXR2; PXR2B; TRIP8b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG219522 representing NM_016559
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTACCAGGGACACATGCAGAAAAGTAAAGAACAAGGATATGGAAACTAAGCAGTGATGAAGACCTCG
 AAATAATTGTTGATCAAAAGCAGGGAAAAGGCTCTAGGGCGGCAGATAAGGCTGTTGCCATGGTGATGAA
 GGAGATACCGAGGGAGGAGTCTGCTGAAGAAAAGCCCTCCTTACTATGACATCAGAGCTGGTGAATGAG
 CAACAAGAAAGCAGACCCCTCTGAGTCCCTCCATCGATGACTTTCTCTGTGAAACCAAAATCGGAAGCAA
 TAGCAAGGCCAGTAACATCCAATACAGCTGTATTGACCACTGGCTTAGATCTCCTCGACCTGAGTGAACC
 AGTCTCTCAAACCCAAACCAAGCAAGAAGTCAGAGCCCTCATCAAAACCTCATCCCTCAAGAAAAAG
 GCCGATGGATCTGACCTCATCAGCACGGATGCTGAGCAGAGAGGCCAGCCTCTCAGAGTCCCGGAGACTT
 CATCCTTAGATCTAGACATTCAACACAACCTGGAAAAATGGGACGATGTTAAGTTTCATGGAGATCGAAA
 TACCAAGGGACATCCAATGGCAGAGAGAAAATCATCCTCATCTAGAACTGGATCAAAAGAGCTCTTATGG
 TCCTCAGAACACAGATCTCAACCAGAACTGAGTGGTGGAAAAAGCGCCCTCAACTCTGAGTCGGGCTTCAG
 AATTGGAATTAGTGGCTCCGACTCAGGCTCGACTGACCAAGAACATCGTGGGGAAGCGCATTACTTTTC
 TAGAAACCACTCCTTAGAAGAAGAGTTTGAAAGGGCAAAAGCAGCAGTGGAGTCAGATACAGAGTTTGG
 GATAAGATGCAAGCAGAAATGGGAAGAAATGGCTCGGAGGAAGTGGATATCTGAGAACCAAGAAGCCGAGA
 ACCAAGTAACCATCTCGGCTAGTGAGAAGGGATATTACTTTCACACTGAAACCCCTTCAAGGACTGGCC
 TGGAGCATTGGAAGAAGGCTTAAAAAGGCTGAAGGAAGGGGATCTGCCAGTCACCATCCTGTTTCATGGAA
 GCAGCAATTCTTCAGGACCTGGAGATGCAGAGGCATGGCAGTTCCTCGGATAACCCAGGCGGAGAATG
 AAAATGAACAAGCAGCTATTGTGCCCCTCAGAGGTGCTTAGAATTACAGCCCAACAACCTTAAAGCTTT
 GATGGCCTTGGCTGTGAGTTATACTAACACTGGCCATCAGCAGGATGCCTGTGACGCTCTGAAGAATTGG
 ATTAAGCAAAATCCAAGTACAATACCTTGTAAGGAAGCAAGAAGGGATCTCCAGGCCTCACCCGCGGA
 TGTCTAAGTCCCCAGTTGATAGCTCTGTTCTGGAAGGGGTGAAGGAATTATCTGGAAGCTGCCACCA
 AAATGGAGATATGATCGACCCAGACCTGCAGACAGGTCTAGGGGTTCTGTTCCACCTGAGTGGAGAATTT
 AATAGAGCAATAGATGCATTTAACGCTGCCTTAAGTCTCGGCCAGAGGACTATTCATATGGAACCGCC
 TCGGGGCGACCTTGGCGAACGGAGACCGCAGCGAGGAAGCCGTGGAGGCCTATACGCGAGCACTGGAGAT
 TCAGCCAGGATTCATCCGGTCCAGATACAACCTAGGAATAAGCTGCATCAACCTGGGCGCCTACAGAGAA
 GCGGTGAGCAATTTCTCACTGCCCTCAGTTTGCAAAGAAAGAGCAGGAATCAGCAGCAAGTTCTCATC
 CTGCAATCTCTGGAATATCTGGGCTGCCCTCAGAATTGCGCTCTCTCTGATGGACCAACCAAGAACTCTT
 CCAGGCGGCTAATCTTGGTGACCTGGATGTCCTCTTAAGAGCTTTCAACTTGGATCCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG219522 representing NM_016559
 Red=Cloning site Green=Tags(s)


MYQGHMQSKEQGYGKLSSDEDLEIIVDQKQGGKSRAADKAVAMVMKEIPREESAEKPLLMTMSQLVNE
 QQESRPLLSPSIDDFLCETKSEAIARVTSNTAVLTTGLDLLDLSEPVSTQTKAKKSEPSSTSSSLKKK
 ADGSDLISTDAEQRGQPLRVPETSSLDLDIQTQLEKWDDVKFHGDRNTKGHPMAERKSSSRTGSKELLW
 SSEHRSQPELSSGKSALNSESASELELVAPTQARLTKEHRWGSALLSRNHSLEEEFERAKAAVESDTEFW
 DKMQAEWEEMARRNWISENQEAQNQVTISASEKGYFHTENPFKDWPGAFEEGLKRLKEGDLPTVILFME
 AAILQDPGDAEAWQFLGITQAENENEQAAIVALQRCLELQPNNLKALMALAVSYTNTGHQQDACDALKNW
 IKQNPKYKYLKSKKSGPLTRRMSKSPVDSSVLEGVKELYLEAAHQNGMDIDPLQTGLGVLFHLSGEF
 NRAIDAFNAALTVRPEDYSLWNRLGATLANGDRSEEAVEAYTRALEIQPGFIRSRYNLGISCINLGAYRE
 AVSNFLTALSLQRKSRNQVPHPAISGNIWAALRIALSLMDQPELFQAANLGDLVDVLLRAFNLDP

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_016559

ORF Size: 1878 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_016559.3](#)


RefSeq Size: 3577 bp

RefSeq ORF: 1881 bp

Locus ID: 51555

UniProt ID:	Q8IYB4
Cytogenetics:	3q26.33
Domains:	TPR
Protein Families:	Druggable Genome
Gene Summary:	Accessory subunit of hyperpolarization-activated cyclic nucleotide-gated (HCN) channels, regulating their cell-surface expression and cyclic nucleotide dependence.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG219522