

Product datasheet for **RG219611**

THOC5 (NM_001002879) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	THOC5 (NM_001002879) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	THOC5
Synonyms:	C22orf19; Fmip; fSAP79; PK1.3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG219611 representing NM_001002879
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCATCAGAATCGAGCAAAAAACGGAAGCCCAAAGTGATCCGAAGCGATGGAGCCCCAGCTGAAGGAA
AGCGGAATCGATCTGACACCGAGCAGGAAGGTAATACTACAGTGAGGAGCCGAGGTGGATCTGCGGGA
CCCTGGCAGAGACTATGAGTTATACAAGTACACCTGCCAGGAGCTACAGAGGCTGATGGCTGAGATCCAA
GACCTGAAGAGCAGGGGTGGCAAGGATGTGGCAATAGAAATAGAAGAACGGAGGATCCAGAGCTGTGTGC
ATTTTCATGACTCTAAAGAAGCTTAACCGATTAGCCACATCAGGTTGAAGAAAGGAGAGATCAGACCCA
CGAGGCTAAGCAGAAAAGTAGATGCCTATCATCTGCAGCTCCAGAACCTGTTGTATGAGGTGATGCACCTA
CAGAAGGAGATCACCAAATGTTTGGAGTTAAGTCAAAGCATGAAGAAATTGATCTGGTCAGTTTAGAGG
AGTTTTATAAGGAGGCTCCACCAGATATCAGCAAGGCCGAAGTCACCATGGGAGACCCCTCACCAGCAAAC
ACTGGCACGTCTGGACTGGGAGCTGGAGCAGCGAAAAGGCTGGCAGAGAAGTACCGAGAGTGCCTATCT
AACAAAGGAGAAGATTCTCAAGGAGATTGAGGTGAAGAAGGAGTACCTGAGCAGCCTCCAGCCCCGCCTCA
ACAGCATCATGCAGGCTTCCCTTCCGGTGCAGGAGTACCTGTTTATGCCATTCGACCAAGGCTCACAAAGCA
GTATGAGACAGCCAGACACCTGCCGCCTCCCTCTATGTCCTCTTTGTTCAAGCCACTGCGTATGGGCGAG
GCCTGTGATAAGACGTTATCTGTGGCAATCGAAGGCAGTGGATGAAGCCAAGGCTCTGTTCAAACCTC
CAGAGGACTCCCAAGATGACGAGAGTGAATCAGATGCCGAGGAGGAGCAGACTACGAAGCGCCGGAGACC
CACACTGGGGTTTCAGTTGGACGACAAACGCAAGGAGATGCTGAAGAGGCACCCACTGTCTGTGATGCTC
GACCTGAAGTGCAAAGATGACAGTGTGCTTACCTGACTTTCTACTACCTCATGAACCTCAACATCATGA
CAGTAAAAGCCAAAGTGACAACCTGCCATGGAGCTGATCACCCCATCAGTGCAGGTGACTTGTGTCTCC
TGACTCAGTCCTGAGTTGCTTGTATCCTGGGGATCATGGAAAGAAAACCTCCGAATCCAGCCAATCAGTAT
CAGTTTGATAAAGTTGGCATCCTGACTTTGAGCGACTATGTACTTGAGCTAGGTCACCCCTATTTGTGGG
TGCAGAAGCTGGGTGGCCTCCACTTCCCAAAGAGCAGCCCAAGCAACAGTGATTGCTGACCACTCGCT
GAGCGCCAGCCACATGGAGACCACCATGAAACTTCTGAAGACCAGGGTGCAGTCCCGCCTGGCCCTCCAC
AAACAGTTTGCATCCCTAGAACATGGCATTGTGCCAGTTACCAGTGATTGCCAGTACCTCTTCCCTGCCA
AGGTTGTCTCTCGCCTGGTGAATGGGTGACAGTTGCCCATGAGGATTACATGGAGCTGCACTTCAACAA
AGACATTGTGGATGCGGGACTGGCTGGGACACCAATCTCTACTACATGGCGCTCATCGAAAGGGGCACA
GCCAAACTGCAGGCCGCTGTGGTGTGAACCTGGCTACTCCTCCATCCCACCTGTTTTCCAGCTCTGTT
TGAAGTGGAAAGGGGAGAAAACCAACAGCAACGATGACAACATTCGGGCCATGGAGGGCGAAGTCAATGT
GTGCTACAAGGAGCTGTGTGGCCCTTGGCCAGCCACCAGCTGTTGACCAACCAGCTGCAGCGGCTGTGT
GTGCTGTGGATGTTTACCTGGAGACCGAGAGCCATGACGACAGTGTGGAGGGGCCCAAGGAATTTCCCC
AGGAGAAGATGTGTCTCGGCTCTTCAGGGTCTTAGCAGGATGAAGCCATTTAAATACAACCATCTCA
GGGATTCTTCAGCCATCGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG219611 representing NM_001002879
 Red=Cloning site Green=Tags(s)


MSESSKKRKP VIRSDGAPAEGKRNRSDEQEGKYYSEEAEVDLRDPGRDYEL YKYTCQELQRLMAEIQ
 DLKSRGGKDVAIEIEERRIQSCVHFM LKLNRLAHIRLKKGRDQTHEAKQKVDAYHLQLQNLLYEV MHL
 QKEITKCLEFKSKHEEIDLVSLEEFYKEAPPDISKA EVTMGDPHQQT LARLDWELEQRKRLAEKYRECLS
 NKEKILKEIEVKKEYLSSLQPR LNSIMQASLPVQEYLFMPFDQAHKQYETARHLPPPLYVLFVQATAYGQ
 ACDKTL SVAIEGSVDEAKALFKPPEDSQDDESDSDAE EEQTTKRRRPTLGVQLDDKRKEM LKRHPLSVML
 DLKCKDSDV LHLTFYLLMNLNIMTVKAKVTTAMELITPISAGDLLSPDSVLSCLYPGDH GKKT PNPANQY
 QFDKVGIL T LSDYVLELGH PYLWVQKLGGLHFPKEQPQQT VIADHSL SASHMETT MLLKTRVQSRLALH
 KQFASLEHG IVPVTSDCQYLFPAKVVSRLVKWVTVAHEDY MELHFTKDIVDAGLAGDTNL YYMALIERGT
 AKLQAAVVL NPGYSSIPPV FQLCLNWKGEK TNSNDNIRAMEGEVNV CYKELCGPWP SHQLL TNQLQRLC
 VLLDVYLETESH DSV EGPKEFPQEKMCLRLFRG P SRMKPFKYNHPQGFFSHR


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:


ACCN: NM_001002879

ORF Size: 2049 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001002879.1](#), [NP_001002879.1](#)

RefSeq Size: 2578 bp

RefSeq ORF: 2052 bp

Locus ID: 8563

UniProt ID: [Q13769](#)

Cytogenetics: 22q12.2

Gene Summary: Acts as component of the THO subcomplex of the TREX complex which is thought to couple mRNA transcription, processing and nuclear export, and which specifically associates with spliced mRNA and not with unspliced pre-mRNA. TREX is recruited to spliced mRNAs by a transcription-independent mechanism, binds to mRNA upstream of the exon-junction complex (EJC) and is recruited in a splicing- and cap-dependent manner to a region near the 5' end of the mRNA where it functions in mRNA export to the cytoplasm via the TAP/NFX1 pathway. The TREX complex is essential for the export of Kaposi's sarcoma-associated herpesvirus (KSHV) intronless mRNAs and infectious virus production. THOC5 in conjunction with ALYREF/THOC4 functions in NXF1-NXT1 mediated nuclear export of HSP70 mRNA; both proteins enhance the RNA binding activity of NXF1 and are required for NXF1 localization to the nuclear rim. Involved in transcription elongation and genome stability. Involved in alternative polyadenylation site choice by recruiting CPSF6 to 5' region of target genes; probably mediates association of the TREX and CFIm complexes.[UniProtKB/Swiss-Prot Function]