

## Product datasheet for **RG219798**

### IGF2 (NM\_000612) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** IGF2 (NM\_000612) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** IGF2  
**Synonyms:** C11orf43; GRDF; IGF-II; PP9974; SRS3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG219798 representing NM\_000612  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGAATCCCAATGGGGAAGTCGATGCTGGTCTTCTCACCTTCTTGGCCTTCGCCTCGTCTGCATTG  
 CTGCTTACCGCCCCAGTGAGACCCTGTGCGCGGGGAGCTGGTGGACACCCTCCAGTTCGCTGTGGGA  
 CCGCGCTTCTACTTCAGCAGGCCCGAAGCCGTGTGAGCCGTCGAGCCGTGGCATCGTTGAGGAGTGC  
 TGTTCCGCAGCTGTGACCTGGCCCTCCTGGAGACGTACTGTGCTACCCCGCCAAGTCCGAGAGGGACG  
 TGTCGACCCTCCGACCGTCTCCGGACAACCTCCCAAGATACCCCGTGGGCAAGTTCTCCAATATGA  
 CACCTGGAAGCAGTCCACCCAGCGCCTGCGCAGGGGCTGCTGCCCTCCTGCGTCCCGCGGGGTAC  
 GTGCTCGCCAAGGAGCTCGAGGCGTTCAGGGAGGCCAAACGTACCCGTCCTGATTGCTCTACCCACC  
 AAGACCCCGCCACGGGGCGCCCCCAGAGATGGCCAGCAATCGGAAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG219798 representing NM\_000612  
 Red=Cloning site Green=Tags(s)

MGIPMGKSMVLVLLTFLAFASCCIAAYRPSETLCGGELVDTLQFVCGDRGFYFSRPASRVSRRSRGIVEEC  
 CFRSCDLALLETYCATPAKSERDVSTPPTVLPDNFPRYPVGKFFQYDTWKQSTQRLRRLPALLRARRGH  
 VLAKELEAFREAKRHRPLIALPTQDPAHGGAPPEMASNRK

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



**Cloning Scheme:**


**ACCN:** NM\_000612

**ORF Size:** 540 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_000612.6](#)

**RefSeq Size:** 1640 bp

**RefSeq ORF:** 543 bp

**Locus ID:** 3481

**UniProt ID:** [P01344](#)


**Cytogenetics:** 11p15.5

**Domains:** IIGF

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein

**Gene Summary:** This gene encodes a member of the insulin family of polypeptide growth factors, which are involved in development and growth. It is an imprinted gene, expressed only from the paternal allele, and epigenetic changes at this locus are associated with Wilms tumour, Beckwith-Wiedemann syndrome, rhabdomyosarcoma, and Silver-Russell syndrome. A read-through INS-IGF2 gene exists, whose 5' region overlaps the INS gene and the 3' region overlaps this gene. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2010]

**Product images:**



Circular map for RG219798