

Product datasheet for **RG219809**

CAMKK2 (NM_172216) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CAMKK2 (NM_172216) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CAMKK2
Synonyms:	CAMKK; CAMKKB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG219809 representing NM_172216
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGTCATCATGTGTCTCTAGCCAGCCAGCAGCAACCGGGCCGCCCCAGGATGAGCTGGGGGCGAGGG
 GCAGCAGCAGCAGCGAAAGCCAGAAGCCCTGTGAGGCCCTGCGGGGCCTCTCATCCTTGAGCATCCACCT
 GGGCATGGAGTCTTCATTGTGGTCACCGAGTGTGAGCCGGGCTGTGCTGTGGACCTCGGCTTGGCGCGG
 GACCGGCCCTGGAGGCCGATGGCCAAGAGGTCCCCCTTGACACCTCCGGGTCCAGGCCCGGCCACC
 TCTCCGGTCGCAAGCTGTCTCTGCAAGAGCGGTCCAGGGTGGGCTGGCAGCCGGTGGCAGCCTGGACAT
 GAACGGACGCTGCATCTGCCCGTCCCTGCCCTACTCACCCGTGAGTCCCGCAGTCTCGCCTCGGCTG
 CCCCAGCGGCCGACAGTGGAGTCTACCACGTCCATCACGGGTATGCAGGACTGTGTGCAGCTGAATC
 AGTATACCCTGAAGGATGAAATTGGAAAGGGCTCCTATGGTGTGCTCAAGTTGGCCTACAATGAAAATGA
 CAATACCTACTATGCAATGAAGGTGCTGTCCAAAAGAAGCTGATCCGGCAGGCCGGCTTCCACGTCGC
 CCTCCACCCGAGGCACCCGGCCAGCTCCTGGAGGCTGCATCCAGCCAGGGGCCCCATTGAGCAGGTGT
 ACCAGGAAATTGCCATCCTCAAGAAGCTGGACCACCCCAATGTGGTGAAGCTGGTGGAGGTCCTGGATGA
 CCCCAATGAGGACCATCTGTACATGGTGTTCGAACTGGTCAACCAAGGGCCCGTATGGAAGTGGCCACC
 CTCAAACCACTCTCTGAAGACCAGGCCGTTTCTACTTCCAGGATCTGATCAAAGGCATCGAGTACTTAC
 ACTACCAGAAGATCATCCACCGTACATCAAACCTTCAAACCTCCTGGTCGGAGAAGATGGGCACATCAA
 GATCGTGACTTTGGTGTGAGCAATGAATCAAGGGCAGTGACGCGCTCCTCTCAAACACCGTGGGCACG
 CCCGCCTCATGGCACCCGAGTCTCTCTGAGACCCGCAAGATCTTCTGGAAGGCCCTGGATGTTT
 GGGCCATGGGTGTGACACTATACTGCTTTGTCTTTGGCCAGTGCCATTATGGACGAGCGGATCATGTG
 TTTACACAGTAAGATCAAGAGTCAAGCCCTGGAATTTCCAGACCAGCCGACATAGCTGAGGACTTGAAG
 GACCTGATCACCCGATGCTGGACAAGAACCCGAGTCAAGGATCGTGGTCCCGGAAATCAAGATCCTGG
 TGAAGACCATGATACGTAACGCTCCTTTGGGAACCCATTGAGGGCAGCCGGCGGGAGGAACGCTCACT
 GTCAGCGCCTGGAAACTTGCTCACCAAAAAACCAACCAGGGAATGTGAGTCCCTGTCTGAGCTCAAGGAA
 GCAAGGCAGCGAAGACAACCTCCAGGGCACCGACCCGCCCCCGTGGGGGAGGAGGAAAGTCTTGTGA
 GAGGCAGTCCCTGCGTGGAAAGTTGCTGGGCCCCCGCCCCGGCTCCCCGCACGCATGCATCCACTGCG
 GCCGGAGGAGCCATGGAGCCCGAG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG219809 representing NM_172216
 Red=Cloning site Green=Tags(s)

MSSCVSSQPSSNRAAPQDELGGRGSSSESQKPCALRGLSSL SIHLGMESFIVVTECEPGCAVDLGLAR
 DRPLEADGQEVPLDTSQSARPHLSGRKLSLQERSQGLAAGGSLDMNGRIPCPSLPYSPVSSPQSSPRL
 PRRPTVESHVSIITGMQDCVQLNQYTLKDEIGKGSYGVVVKLAYNENDNTYYAMKVL SKKKLIRQAGFPRR
 PPPRGTRPAPGGCIQPRGPIEQVYQEIAILKKLDHPNVVKLVEVLDDPNEDHLYMVFELVNQGPVMEVPT
 LKPLSEDQARFYFDL IKGIEYLHYQKIIHRDIKPSNLL VGEDGHIKIADFGVSNFKGSDALLSNTVGT
 PAFMAPESLSETRKIFSGKALDVWAMGVTLYCFVFGQCPFMDERIMCLHSKIKSQALEFPDQPDIAEDLK
 DLITRMLDKNPESRIVVPEIKILVKTMIKRSFGNPFEGSRREERSLSAPGNLLTKKPTRECESLSELKE
 ARQRRQPPGHRPAPRGGGSALVRGSPCVESCWAPAPGSPARMHPLRPEEAMEPE


TRTRPLE – GFP Tag – V

Restriction Sites:


Sgfl-Mlul

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN:	NM_172216
ORF Size:	1635 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172216.1 , NP_757365.1
RefSeq Size:	5491 bp
RefSeq ORF:	1638 bp
Locus ID:	10645
UniProt ID:	Q96RR4
Cytogenetics:	12q24.31
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways:	Adipocytokine signaling pathway
Gene Summary:	The product of this gene belongs to the Serine/Threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. The major isoform of this gene plays a role in the calcium/calmodulin-dependent (CaM) kinase cascade by phosphorylating the downstream kinases CaMK1 and CaMK4. Protein products of this gene also phosphorylate AMP-activated protein kinase (AMPK). This gene has its strongest expression in the brain and influences signalling cascades involved with learning and memory, neuronal differentiation and migration, neurite outgrowth, and synapse formation. Alternative splicing results in multiple transcript variants encoding distinct isoforms. The identified isoforms differ in their ability to undergo autophosphorylation and to phosphorylate downstream kinases. [provided by RefSeq, Jul 2012]