

Product datasheet for **RG219911**

PRB1 (NM_199354) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PRB1 (NM_199354) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PRB1
Synonyms: PM; PMF; PMS; PRB1L; PRB1M
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG219911 representing NM_199354
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTGTTGATTCTGCTGTCTGCTGAGTGGCCTTGTGGCCCTGAGCTCAGCTCAGAACTTAAATGAAGATGTCA
 GCCAGGAAGAATCTCCCTCCCTAATAGCAGGAAATCCACAAGGACCATCCCCACAAGGAGGCAACAAGCC
 CCAAGGCCCCACCTCCTCCAGGAAAGCCACAAGGACCACCCCAAGGAGGCAACAACCTCAAGGT
 CCCCCACCTCCAGGAAAGCCACAAGGACCACCCCAAGGGACAAGTCCCGAAGTCCCGATCTCCTC
 CAGGAAAACCACAAGGACCACCCCAAGGAGGAAAGCCACAAGGACCACCCGACAAGGAGGCAGCAA
 GTCCCAAAGTGCCCGATCTCCTCAGGAAAGCCACAAGGACCACCCCAACAAGAAGGCAACAATCTCAA
 GGTCCCCACCTCCAGCAGGAGGCAATCCCCAGCAGCCTCAGGCACCTCTGCTGGACAGCCCCAGGGAC
 CACCACGCCCTCTCAAGGGGGCAGACCTTCCAGACCTCCCCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG219911 representing NM_199354
 Red=Cloning site Green=Tags(s)

MLLILLSVALLALSSAQNLDVDSQEEPSLIAGNPQGPSQGGNKPQGPPPPGKPPGPPQGGNKPQG
 PPPGKPPGPPQGDKSRSPRSPGKPPGPPQGGKPPGPPAQGGSKSQSARSPPGKPPGPPQEGNPNPQ
 GPPPPAGGNPQQPAPPAGQPQGPPRPPQGGRPSRPPQ

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_199354

ORF Size: 534 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199354.1](#), [NP_955386.1](#)

RefSeq Size: 714 bp

RefSeq ORF: 537 bp

Locus ID: 5542


UniProt ID: [P04280](#)

Cytogenetics: 12p13.2

Protein Families: Druggable Genome

Gene Summary: This gene encodes a member of the heterogeneous family of basic, proline-rich, human salivary glycoproteins. The encoded preproprotein undergoes proteolytic processing to generate one or more mature peptides before secretion from the parotid glands. Multiple alleles of this gene exhibiting variations in the length of the tandem repeats have been identified. The reference genome encodes the "Medium" allele. This gene is located in a cluster of closely related salivary proline-rich proteins on chromosome 12. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar proteolytic processing. [provided by RefSeq, Nov 2015]

Product images:



Circular map for RG219911