

Product datasheet for **RG219936**

ZNF586 (NM_001077426) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: ZNF586 (NM_001077426) Human Tagged ORF Clone
 Tag: TurboGFP
 Symbol: ZNF586
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-AC-GFP (PS100010)
 E. coli Selection: Ampicillin (100 ug/mL)
 ORF Nucleotide Sequence: >RG219936 representing NM_001077426
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGCAGCAGCCGCTCTGAGGGCGCCTGCTCAGGTTGTTGGCATGGAGGGGAAGATGAGGCAGCACC
 TTCTAAGCAGAGCACGTGTATACATATATACAAAGACCAGGGAGGTCATAGTGGAGAAAGGCCTTATGAG
 TGTGGGGAATATAGGAAATTATTTAAGAACAAGTCCTGCCTCACTGAACCCAGAAGAGATCACAAACACA
 GGAATGTTGCGACTGGAGAAAGGCCTTATGAGTGCAGCAAATATGGGAAATTATTCACCAAAAGCCTAC
 ACTCCATATTCATGAGAGATTTTCACTGGGCAAAAGACCTATGAGTGCAGTGAGTGTGGAAAATCATT
 CACCAAAAGCTCTTCACTCTTGCAGCGTCAGACACTTCACACTAGAGAAAGGCCTTATGAGTGTATTGAAT
 GTGGGAAAGCCTTTGCTGAAAAGTCCAGTCTCATTAACCACAGGAAAGTCACTCTGGAGCAAAGCGTTA
 TGAATGAATGAATGTGGGAAGTCCTTTGCTTATACATCTAGTCTCATTAAACACAGGAGGATTCACACT
 GGAGAGAGGCCTTATGAGTGCAGTGAATGTGGGAGATCCTTTGCTGAAAAGTCCAGTCTTATTAAACACT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG219936 representing NM_001077426
 Red=Cloning site Green=Tags(s)


MAAAAALRAPAQVVGMEGKMRQHLLSRARVYIYTKTREVIVEKGLMSVGNIGNYLRTSPASLNPEEITNT
 GMFALEKGLMSAANMGNYFTKSLHSIFMRDFILGKRPMASVSVENHFTKALHSCSVRHFTLEKGLMSVLN
 VGKPLLKSPVSLTTGKFTLEQSVNMAMNVGSPLLIHLVSLNTGGFTLERGLMSAVNVGDPLLKTPVLLNT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001077426

ORF Size: 630 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001077426.2](#)

RefSeq Size: 2088 bp

RefSeq ORF: 633 bp

Locus ID: 54807
UniProt ID: [Q9NXT0](#)
Cytogenetics: 19q13.43
Protein Families: Transcription Factors
Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG219936