

## Product datasheet for **RG220149**

### **P2Y10 (P2RY10) (NM\_198333) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	P2Y10 (P2RY10) (NM_198333) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	P2Y10
Synonyms:	LYPSR2; P2Y10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG220149 representing NM_198333 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTAACCTTGACAAATACACTGAAACATTCAAGATGGGTAGCAACAGTACCAGCACTGCTGAGATTT  
ACTGTAATGTCACTAATGTGAAATTTCAACTCCCTCTATGCAACCACCTATATCCTCATATTCATTCC  
TGGTCTTCTGGCTAACAGTGCAGCCTTGTGGTCTGTGCCGCTTCATCAGCAAGAAAAATAAGCCATC  
ATTTTCATGATCAACCTCTCTGTGGCTGACCTTGCTCATGTATTATCTTTACCCCTCCGGATTTACTATT  
ACATCAGCCACCACTGGCCTTTCCAGAGAGCCCTTGCCTGCTCTGCTTCTACCTGAAGTATCTCAACAT  
GTATGCCAGCATTGTTTCCTGACGTGCATCAGTCTTCAAAGGTGCTTTTTCTCCTCAAGCCCTTCAGG  
GCCAGAGACTGGAAGCGTAGGTACGATGTGGGCATCAGTCTGCCATCTGGATCGTTGTGGGGACTGCAT  
GTTTGCCATTTCCCATCCTGAGAAGCACAGACTTAAACAACAAGTCCTGCTTTGCTGATCTTGGATA  
CAAGCAAATGAATGCAGTTGCGTTGGTCGGGATGATTACAGTTGCTGAGCTTGCAGGATTTGTGATCCCA  
GTGATCATCATCGCATGGTGTACCTGAAAACTACTATATCCTTGAGACAGCCACCAATGGCTTTCCAAG  
GGATCAGTGAGAGGCAGAAAGCACTGCGGATGGTGTTCATGTGTGCTGCAGTCTTCTTCATCTGCTTAC  
TCCCTATCATATTAACCTTATTTTTTACACCATGGTAAAGGAAACCATCATTAGCAGTTGCCCCGTTGTC  
CGAATCGCACTGTATTTCCACCCTTTTTGCCTGTGCCTTGAAGTCTCTGCTGCCTTTTGGATCCAATTC  
TTTATTACTTTATGGCTTCAGAGTTTCGTGACCAACTATCCC GCCATGGCAGTTCTGTGACCCGCTCCCG  
CCTCATGAGCAAGGAGAGTGGTTCATCAATGATTGGC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG220149 representing NM\_198333  
Red=Cloning site Green=Tags(s)

MANLDKYTETFKMGSNSTSTAEIYCNVTNVKFQYSLYATTYILIFIPGLLANSAALWVLCRFISKKNKAI  
 IFMINLSVADLAHVLSLPLRIYYIYSHHWPFQRALCLLCFYLYLNMYSICFLTCISLQRCFFLLKPF  
 ARDWKRRYDVGISAAIWI VVGTA CLPFPILRSTDLN NNKSCFADLGYKQMNAVALVGMITVAELAGFVIP  
 VIIIAWCTWKTITISLRQPPMAFQGISERQKALRMVFMCAAVFFICFTPYHINIFIFYTMVKETI ISSCPV  
 RIALYFHPFC LCLASLCLLDPILYYFMASEFRDQLSRHGSSVTRSR LMSKESGSSMIG

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_198333

**ORF Size:** 1017 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_198333.1</a> , <a href="#">NP_938147.1</a>
<b>RefSeq Size:</b>	1679 bp
<b>RefSeq ORF:</b>	1020 bp
<b>Locus ID:</b>	27334
<b>UniProt ID:</b>	<a href="#">O00398</a>
<b>Cytogenetics:</b>	Xq21.1
<b>Protein Families:</b>	Druggable Genome, GPCR, Transmembrane
<b>Protein Pathways:</b>	Neuroactive ligand-receptor interaction
<b>Gene Summary:</b>	The protein encoded by this gene belongs to the family of G-protein coupled receptors that are preferentially activated by adenosine and uridine nucleotides. There is a pseudogene for this gene nearby on chromosome X. Multiple alternatively spliced transcripts have been observed. [provided by RefSeq, Apr 2016]

## Product images:



Circular map for RG220149