

Product datasheet for **RG220192**

ZNF433 (NM_001080411) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF433 (NM_001080411) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF433
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG220192 representing NM_001080411
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGATGTTTCAGGATTCAGTGGCCTTTGAGGATGTGGCTGTGACCTTCACCCAAGAGGAGTGGGCTTTGC
 TGGATCCTTCCCAGAAAAATCTCTGTAGAGATGTGATGCAAGAAACCTTCAGGAACCTGGCCTCTATAGG
 GAAAAATGGAACCCCAAGACATATATGTAGAGTACGAAAACTAAGGAGAAACCTAAGAATTGTGGGA
 GAGAGACTCTTTGAAAGTAAAGAAGTCACTCAGCATGGAGAAATTTGACCCAGGTTCCAGATGACATGC
 TGAAGAAAAACAACACTACTGGAGTAAATCATGCGAAAGCAGTGTGTATGGAGAAGTAGGCAGTGCTATTC
 ATCTCTTAATAGGCACATCAGAGATGACACTGGACACAAGGCATATGAGTATCAAGAATATGGACAGAAA
 CCATATAAATGTAATACTGTAAAAACCTTTCAACTGTCTCTCCTCTGTTCCAGACATGAAAGGGCTC
 ATAGTGAAGGAACTCTATGTTTGTGAGGAATCGGAAAAACATTTATTTCCATTCAAACCTTCAAAG
 ACACAGGATAATGCACCGTGGAGATGGACCTTATAAGTGTAAATTTGTGGGAAAGCCTTGATGTTTCTC
 AGTTTGTATCTTATCCACAAACGAACTCACACTGGAGAGAAACCATATCAATGTAAACAGTGTGGTAAAG
 CCTTTAGTCATTCTAGTAGCCTTCGAATACATGAAAGAACTCACACTGGGGAGAAGCCTTATAAATGTAA
 TGAATGTGGGAAAGCATTCCATAGTTCCACATGCCTTCATGCTCATAAAGAACTCACACTGGGGAGAAG
 CCATATGAATGTAACAGTGTGGGAAAGCCTTCAGCTCTTCCATTCTTTCAAATACATGAAAGAACTC
 ACACGGGGGAGAAGCCATATGAATGTAAGGAATGTGGAAGCATTCAAGTGTCCAGTTCTGTTTCGCAG
 ACATGAAAGAACCCTCTAGGAAAAACCTATGAATGTAAACATTGTGGGAAAGTATTATCTTATCTT
 ACCAGCTTTCAAACCACTTGGGAATGCACACTGGAGAGATATCTCATAAATGTAAGATATGTGGGAAAG
 CCTTTTATCTCCAGTTCACCTCAAACACATGAAAAAACTCACACTGGAGAGAAACCTATAAATGCAA
 CCAATGTGGTAAAGCCTTAATTCTTCCAGTTCCTTCCGATATCATGAAAGAACTCACACTGGAGAGAAA
 CCTTACGAGTGTAAAGCAATGTGGGAAAGCCTTCAGATCTGCCTCACTCCTTCAAACACATGGTAGGACTC
 ACACGGGAGAGAAACCTATGCATGTAAGGAATGTGGAAGCATTAGTAATTTCTCTTTCTTTCAAAT
 ACATGAAAGGATGCACAGAGAAGAGAAGCCGTATGAATGTAAGGGTTATGGGAAACATTCAAGTTGCCC
 AGTTTATTTATAGACATGAAAGGACTCACACTGGAGGAAAAACCTATGAATGCAAGCAGTGTGGCAGAT
 CCTTCACTGTTTCGAGCTCCTTTCGATATCATGGAAGGACTCACACTGGAGAGAAACCTATGAATGCAA
 GCAATGTGGGAAAGCCTTCAGATCTGCCTCACAGCTTCAAATTCATGGAAGGACTCACACTGGAGAGAAA
 CCTTATGAATGTAAGCAGTGTGGGAAAGCCTTTGGATCTGCCTCACACCTTCAAATGCATGGAAGGACTC
 ACACTGGAGAGAAACCTATGAATGTAAGCAGTGTGGGAAAGTCTTTGGATGTGCCTCGCGACTTCAAAT
 GCATGGAAGGACTCACACTGGAGAGAAACCGTATAAATGTAAGCAATGTGGGAAAGCTTTTGGATGTCCC
 TCAAACCTTGAAGGCATGGAAGGACTCACACTGGAGAGAAACCTATAAATGTAACCAATGTGGTAAAG
 TCTTTAGATGTTCTTCAAACTTCAAGTGCATGGAAGGGCTCACTGCATAGACACCCCA

ACGCGTACGCGGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG220192 representing NM_001080411
 Red=Cloning site Green=Tags(s)


MMFQDSVAFEDVAVTFTQEEWALLDPSQKNLCRDVMQETFRNLASIGKKWKPNIIYVEYENLRNLRIVG
 ERLFESKEGHQHGEILTQVPDDMLKTTTGVKSCCESSVYGEVSAHSSLNRHIRDDTGHKAYEYQYEQK
 PYKCKYCKKPFNCLSSVQTHERAHSGRKLYVCEECGKTFISHSNLQRHRIMHRGDGPYCKFCGKALMFL
 SLYLIIHKRTHTEKPYQCKQCGKAFSHSSSLRIHERHTHTEKPYKCNECGKAFHSSTCLHAHKRTHTEK
 PYECKQCGKAFSSSHSFQIHERHTHTEKPYECKECGKAFKCPSSVRRHERTHSRKKPYECKHCGKVL SYL
 TSFQNHLMGMTGEISHKCKICGKAFYSPSSLQTHEKTHTEKPYKCNQCGKAFNSSSSFRYHERHTHTEK
 PYECKQCGKAFRSASLLQTHGRTHTEKPYACEKCGKPF SNF SFFQIHERMHREEKPYECKGYGKTFSLP
 SLFHRHERHTHTEKPYECKQCGRSFNCSFFRYHGRTHTEKPYECKQCGKAFRSASQLQIHGRTHTEK
 PYECKQCGKAFGSASHLMHGRTHTEKPYECKQCGKSFGCASRLQMHGRTHTEKPYCKQCGKAFGCP
 SNLRRHGRTHTEKPYKCNQCGKVFRCSQLQVHGRAHCIDTP

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001080411

ORF Size: 2019 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080411.2](#), [NP_001073880.1](#)

RefSeq Size: 2321 bp

RefSeq ORF: 2022 bp

Locus ID: 163059


UniProt ID: [Q8N7K0](#)

Cytogenetics: 19p13.2

Protein Families: Transcription Factors

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG220192