

## Product datasheet for **RG220204**

### FGL1 (NM\_147203) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	FGL1 (NM_147203) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FGL1
Synonyms:	HFREP1; HP-041; HPS; LFIRE-1; LFIRE1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG220204 representing NM_147203 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAAAGGTGTTTCAGTTTCATCCTTGTACCACCGCTCTGACAATGGGCAGGGAAATTCGGCGCTCG  
AGGACTGTGCCAGGAGCAGATGCGGCTCAGAGCCAGGTGCGCCTGCTTGAGACCCGGTCAAACAGCA  
ACAGGTCAAGATCAAGCAGCTTTTGCAGGAGAATGAAGTCCAGTTCCTTGATAAAGGAGATGAGAATACT  
GTCATTGATCTTGAAGCAAGAGGCAGTATGCAGATTGTTTCAGAGATTTTCAATGATGGGTATAAGCTCA  
GTGGATTTTACAAAATCAAACCTCTCCAGAGCCCAGCAGAATTTTCTGTTTATTGTGACATGTCCGATGG  
AGGAGGATGGACTGTAATTCAGAGACGATCTGATGGCAGTGAAAACCTTTAACAGAGGATGGAAAGACTAT  
GAAAATGGCTTTGGAAATTTGTCCAAAACATGGTGAATATTGGCTGGCAATAAAAATCTTCACTTCT  
TGACCACTCAAGAAGACTACACTTTAAAATCGACCTTGAGATTTTGAAGAAATAGCCGTTATGCACA  
ATATAAGAATTTCAAAGTTGGAGATGAAAAGATTTCTACGAGTTGAATATTGGGGAATATTCTGGAACA  
GCTGGAGATTCCTTGCAGGGAATTTTCATCCTGAGGTGCAGTGGTGGGCTAGTCACCAAAGAATGAAAT  
TCAGCACGTGGGACAGAGATCATGACAACTATGAAGGAACTGCGCAGAAGAAGATCAGTCTGGCTGGT  
GTTTAAACAGGTGTCACCTGCAAACCTGAATGGTGTATACTACAGCGGCCCTACACGGCTAAAACAGAC  
AATGGGATTGTCTGGTACACCTGGCATGGTGGTGGTATTCTCTGAAATCTGTGGTTATGAAAATTAGGC  
CAAATGATTTTATTCAAATGTAATT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG220204 representing NM\_147203  
 Red=Cloning site Green=Tags(s)

MAKVFSFILVTTALTMGREISALEDCAQEQMRLRAQVRLLETRVKQQQVVIKQLLQENEVQFLDKGDENT  
 VIDLGSKRQYADCSEIFNDGYKLSGFYKIKPLQSPAEFVSYCDMSDGGGWTVIQRRSDGSENFNRGWKDY  
 ENGFGNFVQKHGEYWLGNKNLHFLTQEDYTLKIDLADFEKNSRYAQYKNFKVKGDEKNFYELNIGEYSGT  
 AGDSLAGNFHPEVQWASHQRMKFSTWDRDHDNYEGNCAEEDQSGWWFNRCHSANLNGVYYSGPYTAKT  
 NGIVWYTWHGWWYSLKSVVMKIRPNDFIPNVI

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_147203

**ORF Size:** 936 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_147203.2</a> , <a href="#">NP_671736.2</a>
<b>RefSeq Size:</b>	1373 bp
<b>RefSeq ORF:</b>	939 bp
<b>Locus ID:</b>	2267
<b>UniProt ID:</b>	<a href="#">Q08830</a>
<b>Cytogenetics:</b>	8p22
<b>Protein Families:</b>	Druggable Genome, Secreted Protein
<b>Gene Summary:</b>	Fibrinogen-like 1 is a member of the fibrinogen family. This protein is homologous to the carboxy terminus of the fibrinogen beta- and gamma- subunits which contains the four conserved cysteines of fibrinogens and fibrinogen related proteins. However, this protein lacks the platelet-binding site, cross-linking region and a thrombin-sensitive site which are necessary for fibrin clot formation. This protein may play a role in the development of hepatocellular carcinomas. Four alternatively spliced transcript variants encoding the same protein exist for this gene. [provided by RefSeq, Jul 2008]