

Product datasheet for **RG220350**

Cytokeratin 4 (KRT4) (NM_002272) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cytokeratin 4 (KRT4) (NM_002272) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cytokeratin 4
Synonyms:	CK-4; CK4; CYK4; K4; WSN1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG220350 representing NM_002272
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**C

ATGACTTCTGTAGGTGTGTTTAGTGACATGCTCAACGGGTGCGGAAGGATGGGCTTGTGCCAAGGGCCA
 AGCCAGAGATGTTTCAGATTTTCCCTTTATGCCCTGCAACCAAGCCCTGCTGCTCCAGGACATATAA
 GAGACGAAGGCTGAGGGCTCCAGCACTACCGGCTGGGCCCTGTCACTTCTCTGATAGCTCCAGCTCG
 CTCTCTGCAGCCATGATTGCCAGACAGCAGTGTGTCCGAGGCGGGCCCGGGGCTTCAGCTGTGGCTCGG
 CCATTGTAGGCGGTGGCAAGAGAGGTGCCTTCAGCTCAGTCTCCATGTCTGGAGGTGCTGGCCGATGCTC
 TTCTGGGGGATTGGCAGCAGAAGCCTCTACAACCTCAGGGGAACAAAAGCATCTCCATGAGTGTGGCT
 GGGTCACGACAAGGTGCCTGCTTTGGGGGTGCTGGAGGCTTTGGCACTGGTGGCTTTGGTGGTGGATTTG
 GGGGCTCCTTCAGTGGTAAGGGTGGCCCTGGCTTCCCCGTCTGCCCGCTGGGGGAATTCAGGAGGTAC
 CATCAACCAGAGCTTGCTACCCCCCTCCAGTGGAGATTGACCTGAGATCCAGAAAGTCCGGACGGAA
 GAGCGCAACAGATCAAGCTCCTCAACAACAAGTTTGCCTCCTTCATCGACAAGGTGCAGTTCTTAGAGC
 AACAGAATAAGGTCTTGAGACCAAAATGGAACCTGCTCCAGCAGCAGACGACCACCTCCAGCAAAAA
 CCTTGAGCCCCTCTTTGAGACCTACCTCAGTGTCTGAGGAAGCAGCTAGATACCTTGGGAATGACAAA
 GGGCGCTGCAGTCTGAGCTGAAGACCATGCAGGACAGCGTGGAGGACTTCAAGACTAAGTATGAAGAGG
 AGATCAACAAACGCACAGCAGCCGAGAATGACTTTGTGGTCTAAAGAAGGACGTGGATGCTGCCTACCT
 GAACAAGGTGGAGTTGGAGGCCAAGGTGGACAGTCTTAATGACGAGATCAACTTCTGAAGGTCTCTAT
 GATGCGGAGCTGTCCAGATGCAGACCCATGTCAGCGACACGTCCGTGGTCTTTCCATGGACAACAACC
 GCAACCTGGACCTGGACAGCATTATTGCCGAGGTCCGTGCCAGTACGAGGAGATTGCCAGAGGAGCAA
 GGCTGAGGCTGAAGCCCTGTACCAGACCAAGGTCCAGCAGCTCCAGATCTCGGTTGACCAACATGGTGAC
 AACCTGAAGAACACCAAGAGTGAATTGCAGAGCTCAACAGGATGATCCAGAGGCTGCGGCAGAGATCG
 AGAACATCAAGAAGCAGTGCCAGACTCTTCAGGTATCCGTGGCTGATGCAGAGCAGCGAGGTGAGAATGC
 CCTTAAAGATGCCACAGCAAGCGCTAGAGCTGGAGGCTGCCCTGCAGCAGGCCAAGGAGGAGCTGGCA
 CGAATGCTGCGTGAGTACCAGGAGCTCATGAGTGTGAAGCTGGCCTTGGACATCGAGATCGCCACCTACC
 GCAAACCTGCTGGAGGGCAGGAGTACAGAATGTCTGGAGAATGCCAGAGTGCCGTGAGCATCTCTGTGGT
 CAGCGGTAGCACCACTGGAGGCATCAGCGGAGGATTAGGAAGTGGCTCCGGGTTTGGCTGAGTAGT
 GGCTTTGGCTCCGGCTCTGGAAGTGGCTTTGGGTTTGGTGGCAGTGTCTCTGGCAGTTCCAGCAGCAAGA
 TCATCTCTACCACCACCTGAACAAGAGACGA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG220350 representing NM_002272
 Red=Cloning site Green=Tags(s)


MTSVGVFSDMLNGCGKDG LVPRAKPRDVSD FSLYAPATKPCCSR TYKRRRLRAPALTGLGPVTS LIAPSS
 LSAAMIARQQCVRG GPRGFSCGSAIVGGGKRGAFSSVSM SGGAGRCSSGGFGSRSLYNLRGNKSISMSVA
 GSRQGACFGGAGGFGTGGFGGGFSGSGGPGFPVCPAGGIQEVTINQSLLTPLHVEIDPEIQKVRTE
 EREQIKLLNNKFASFIDKVQFLEQQNKVLETKWNLLQQTTTTSSKNLEPLFETYLSVLRKQLDTLGNDK
 GRLQSELKTMQDSVEDFKTKYEEEINKRTAAENDFVVLKKDVDAAYLNKVELEAKVDSLNDINFLKVL
 DAELSQM QTHVSDTSVVL SMDNNRNL DLSIIAEVRAQYEEIAQRSKAEAEALYQTKVQQLQISVDQHGD
 NLKNTKSEIAELNRM IQR LRAE IENIKKQCQTLQVSVADAEQRGENALKDAHSKRVELEAALQQAKEELA
 RMLREYQELMSVKLALDIEIATYRKLLGEEYRMSGECQSAVSISVVS GSTSTGGISGGLSGSGFGLSS
 FGSGSGSGFGFGGSVSGSSSKIISTTLNKR

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_002272

ORF Size: 1782 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_002272.2](#), [NP_002263.2](#)

RefSeq Size: 2576 bp

RefSeq ORF: 1563 bp

Locus ID: 3851


UniProt ID: [P19013](#)

Cytogenetics: 12q13.13

Domains: filament

Gene Summary: The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in differentiated layers of the mucosal and esophageal epithelia with family member KRT13. Mutations in these genes have been associated with White Sponge Nevus, characterized by oral, esophageal, and anal leukoplakia. The type II cytokeratins are clustered in a region of chromosome 12q12-q13. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG220350