

Product datasheet for **RG220479**

RPLP1 (NM_213725) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPLP1 (NM_213725) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: RPLP1
Synonyms: LP1; P1; RPP1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG220479 representing NM_213725
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCTCTGTCTCCGAGCTCGCCTGCATCTACTCGGCCCTCATTCTGCACGACGATGAGGTGACAGTCA
CGGCCCTGGCCAACGTCAACATTGGGAGCCTCATCTGCAATGTAGGGCCGGTGGACCTGCTCCAGCAGC
TGGTGTGCACCAGCAGGAGTCTGCCCTCCACTGCTGCTGCTCCAGCTGAGGAGAAGAAAGTGGAA
GCAAAGAAAGAAGAAATCCGAGGAGTCTGATGATGACATGGCTTTGGTCTTTTGGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG220479 representing NM_213725
Red=Cloning site **Green**=Tags(s)

MASVSELACIYSALILHDEVTVTALANVNIGSLICNVGAGGPAPAAGAAPAGGPAPSTAAAPAEKKVE
AKKEESEESDDDMGFLFD


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_213725

ORF Size: 267 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_213725.2
RefSeq Size:	437 bp
RefSeq ORF:	270 bp
Locus ID:	6176
UniProt ID:	P05386
Cytogenetics:	15q23
Protein Pathways:	Ribosome
Gene Summary:	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal phosphoprotein that is a component of the 60S subunit. The protein, which is a functional equivalent of the E. coli L7/L12 ribosomal protein, belongs to the L12P family of ribosomal proteins. It plays an important role in the elongation step of protein synthesis. Unlike most ribosomal proteins, which are basic, the encoded protein is acidic. Its C-terminal end is nearly identical to the C-terminal ends of the ribosomal phosphoproteins P0 and P2. The P1 protein can interact with P0 and P2 to form a pentameric complex consisting of P1 and P2 dimers, and a P0 monomer. The protein is located in the cytoplasm. Two alternatively spliced transcript variants that encode different proteins have been observed. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]