

Product datasheet for **RG220549**

BRCC45 (BRE) (NM_004899) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRCC45 (BRE) (NM_004899) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	BABAM2
Synonyms:	BRCC4; BRCC45; BRE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG220549 representing NM_004899 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCGCCAGAAAGTGGCCTTGAACCGAATATCTCCAATGCTCTCCCTTTCATATCTAGCGTGGTCCGGA
ATGGAAAAGTGGGACTGGATGCTACAACTGTTTGAGGATAACTGACTTAAAATCTGGCTGCACATCATT
GACTCTGGGCCCAACTGTGACCGATTTAACTGCACATACCATATGCTGGAGAGACATTAAGTGGGAT
ATCATTTCATGCCCAATACCCAGAAGTGCCTCCGATTTTATCTTTGGAGAAGATGCTGAATTCCTGC
CAGACCCCTCAGCTTTGCAGAATCTTGCTCCTGGAATCCTTCAAATCCTGAATGTCTTACTTGTGGT
GAAGGAACTTGTGCAACAATATACCAATTCAAATGTAGCCGCTCCGGGAGAGCTCCCGCTCATGTTT
GAATACCAGACATTACTGGAGGACCCAGTATGGAGAGAATGGAAATTTATGCTGGGAAAAAACA
ACTGGACTGGTGAATTTTCAGCTCGTTTCTTTGAAAGCTGCCCGTAGATTTTCAGCAATATCCCCACATA
CCTTCTCAAGGATGTAATGAAGACCCTGGAGAAGATGTGGCCCTCCTCTCTGTTAGTTTGGAGCACT
GAAGCCACCCAGGTGTACCCCAAGCTGTACTTGTACCTCGAATTGAGCATGCACCTGGAGGCTCCTCAG
CTCTTCATATCCAGCTTTTCCAGGAGGAGGATGTCTCATTGATTACGTTCCCTCAAGTATGCCACCTGCT
CACCAACAAGGTGCAGTACGTGATTCAAGGGTATCACAAAAGAAGAGAGTATATTGCTGCTTTTCTCAGT
CACTTTGGCACAGGTGCTGGAATATGATGCAGAAGGCTTTACAAAACACTCACTGCTGATGTGGA
AAGATTTTGTGTTTCTGTACACATTGACCTGCCTCTGTTTTTCCCTCGAGACCAGCCAACCTCACATT
TCAGTCCGTTTATCACTTTACCAACAGTGGACAGCTTTACTCCAGGCCCAAAAAATTATCCGTACAGC
CCCAGATGGGATGAAATGAAATGGCCAAAAGAGCAAAGGGATGCCAAGGGAGCAGAGATGCCTGCAGCC
CGTGGGAGCAAGTCTGGCCTTTGCAGTTGCAAAAACCTGGCTGCAAGCTGCTCCAGCCCCAGAGAACTG
GCCAAGCTCCAGAGGGCCTCCTTGGAGGGCCTCAGAGGGAGAGAACTGCTCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG220549 representing NM_004899
 Red=Cloning site Green=Tags(s)

MSPEVALNRISPMLSPFISSVVRNGKVGLDATNCLRITDLKSGCTSLTPGPNCDFKLHIPPYAGETLKWD
 IIFNAQYPELPPDFIFGEDAEFLPDPALQNLASWNPSNPECLLLVVKELVQYHQFQCSRLRESSRLMF
 EYQTLLEEPQYGENMEIYAGKKNWTGEF SARFLKLPVDFSNIPTYLLKDVNEDPGEDVALLSVSFEDT
 EATQVYPKLYLSPRIEHALGGSSALHIPAFPPGGCLIDYVPQVCHLLTNKVQYVIQGYHKRREYIAAFLS
 HFGTGVVEYDAEGFTKLTLLLMWKDFCFLVHIDLPLFFPRDQPTLTFQSVYHFTNSGQLYSQAQKNYPYS
 PRWDGNEMAKRAKGCQGSRDACSPWEQVLAFAVAKTGCKLLQPQRNWPSSRGPWRASEGERTAQ


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_004899

ORF Size:	1245 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_004899.4 , NP_004890.2
RefSeq Size:	1978 bp
RefSeq ORF:	1248 bp
Locus ID:	9577
UniProt ID:	Q9NXR7
Cytogenetics:	2p23.2
Protein Families:	Druggable Genome
Gene Summary:	This gene encodes an anti-apoptotic, death receptor-associated protein that interacts with tumor necrosis factor-receptor-1. The encoded protein acts as an adapter in several protein complexes, including the BRCA1-A complex and the BRISC complex. The BRCA1-A complex possesses ubiquitinase activity and targets sites of double strand DNA breaks, while the BRISC complex exhibits deubiquitinase activity and is involved in mitotic spindle assembly. This gene is upregulated in several types of cancer. [provided by RefSeq, Jun 2016]