

Product datasheet for **RG220607**

GM CSF Receptor alpha (CSF2RA) (NM_172247) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GM CSF Receptor alpha (CSF2RA) (NM_172247) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CSF2RA
Synonyms:	alphaGMR; CD116; CDw116; CSF2R; CSF2RAX; CSF2RAY; CSF2RX; CSF2RY; GM-CSF-R-alpha; GMCSFR; GMCSFR-alpha; GMR; GMR-alpha; SMDP4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG220607 representing NM_172247 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCTTCTCCTGGTGACAAGCCTTCTGCTCTGTGAGTTACCACACCCAGCATTCTCCTGATCCCAGAGA
AATCGGATCTGCGAACAGTGGCACCAGCCTCTAGTCTCAATGTGAGGTTTGACTCCAGGACGATGAATTT
AAGCTGGGACTGCCAAGAAAACACAACCTTCAGCAAGTGTTCCTAACTGACAAGAAGAAGAGTTCGTG
GAACCCAGGCTCAGTAACAACGAATGTTTCGTGCACATTTTCGTGAAATTTGTCTGCATGAAGGAGTCACAT
TTGAGGTTACGTGAATACTAGTCAAAGAGGATTTCAACAGAACTGCTTTATCCAAATTCAGGAAGGGA
GGTACCGCTGCTCAGAATTTCTCCTGTTTCATCTACAATGCGGATTTAATGAACTGTACCTGGGCGAGG
GGTCCGACGGCCCCCGTGACGTCCAGTATTTTTGTACATACGAACTCAAAGAGAAGGAGGGAGATCC
GGTGTCTTATTACATACAAGACTCAGGAACCCATGTGGGATGTCACCTGGATAACCTGTCAGGATTAAC
GTCTCGCAATTAATTTCTGGTTAACGGAACCCAGCCGAGAAATGGCATCCAATTTCTTGATTCACTTTTG
GACACAAAGAAAATAGAACGATTCAACCCCTCCCAGCAATGTCACCGTACGTTGCAACACGACGCACTGCC
TCGTACGGTGGAAACAGCCAGGACCTATCAGAAGCTGTCGTACCTGGACTTTCAGTACCAGCTGGACGT
CCACAGAAAGAATACCCAGCCTGGCAGGAAAACCTACTGATTAATGTTTCTGGTGATTTGGAAAATAGA
TACAATTTCCAAGCTCTGAGCCAGAGCAAAACACAGTGTGAAGATCAGAGCTGCAGACGTCCGCATCT
TGAATTTGAGCTCCTGGAGTGAAGCCATTGAATTTGGTTCCTTAGGATACAGCGGCTGTTCCCGCCAGTT
CCACAGATCAAAGACAAAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG220607 representing NM_172247
 Red=Cloning site Green=Tags(s)

MLLVTSLLLCELPHPAFLLIPEKSDLRTVAPASSLNVRFDSTRMNL SWDCQENTTF SKCFLTDKKNRVV
 EPRLSNNECSCTFREICLHEGVTFEVHVNTSQRGFQKLLYPNSGREGTAAQNFSCFIYNADLMNCTWAR
 GPTAPRDVQYFLYIRNSKRRREIRCPYYIQDSGTHVGHLDNL SGLTSRNYFLVNGTSREIGIQFFDSSL
 DTKKIERFNPPSNVTVRCNTTHCLVRWKQPRTYQKL SYLDFQYQLDVHRKNTQPGTENLLINVSGLLENR
 YNFPSSPRAKHSVKIRAADVRILNWSWSEAIEFGSLGYSGCSRQFHRSKTN


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_172247

ORF Size:	999 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172247.3
RefSeq Size:	1587 bp
RefSeq ORF:	1002 bp
Locus ID:	1438
UniProt ID:	P15509
Cytogenetics:	X;Y
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway, Pathways in cancer
Gene Summary:	The protein encoded by this gene is the alpha subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Multiple transcript variants encoding different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble. [provided by RefSeq, Jul 2008]