

Product datasheet for **RG220906**

DNMT3A (NM_153759) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DNMT3A (NM_153759) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DNMT3A
Synonyms:	DNMT3A2; HESJAS; M.HsaIIIA; TBRS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG220906 representing NM_153759
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGGATCCTGGAGCGGTTGTGAGAAGGAATGGGCGCTGGATCGTAGCCTGAAAGACGAGTGTGATA
 CGGCTGAGAAGAAAGCCAAGGTCATTGCAGGAATGAATGCTGTGGAAGAAAACCAGGGGCCCGGGAGTC
 TCAGAAGGTGGAGGAGGCCAGCCCTCCTGCTGTGCAGCAGCCACTGACCCCGCATCCCCACTGTGGCT
 ACCACGCCTGAGCCCGTGGGGTCCGATGCTGGGACAAGAATGCCACAAAGCAGGCGATGACGAGCCAG
 AGTACGAGGACGGCCGGGCTTTGGCATTGGGGAGCTGGTGTGGGGAAACTGCGGGGCTTCTCCTGGTG
 GCCAGGCCGATTGTGTCTTGGTGGATGACGGGCCGAGCCGAGCAGCTGAAGGCACCCGCTGGGTCATG
 TGGTTCGAGACGGCAAATTCAGTGGTGTGTGTTGAGAAGCTGATGCCGCTGAGCTCGTTTTGCAGTG
 CGTTCCACCAGGCCACGTACAACAAGCAGCCCATGTACCGCAAAGCCATCTACGAGTCTCGAGGTGGC
 CAGCAGCCGCGCGGGGAAGCTGTCCCGGTGTGCCACGACAGCGATGAGAGTGACTGCCAAGGCCGTG
 GAGGTGCAAGAACAAGCCATGATTGAATGGGCCCTGGGGGGCTTCCAGCCTTCTGGCCCTAAGGGCCTGG
 AGCCACCAAGAAGAAGAGAAGAATCCCTACAAGAAGTGTACACGGACATGTGGGTGGAACCTGAGGCAGC
 TGCCACGACACCCTCCACCAGCCAAAAGCCCGGAAGAGCACAGCGGAGAAGCCCAAGGTCAAGGAG
 ATTATTGATGAGCGCACAAAGAGACGGCTGGTGTACGAGGTGCGGAGAGTGCAGGAACTTGGAGACA
 TCTGCATCTCCTGTGGGAGCCTCAATGTTACCCGTGAACACCCCTCTTTCGTTGGAGGAATGTGCCAAAA
 CTGCAAGAACTGCTTTCTGGAGTGTGCTACCAGTACGACGACGACGGCTACCAGTCTACTGCACCATC
 TGCTGTGGGGCCGTGAGGTGCTCATGTGCGGAAACAACAACCTGCTGCAGGTGCTTTTTCGTTGGAGTGTG
 TGGACCTTTGGTGGGGCCGGGGCTGCCAGGCAGCCATTAAGGAAGACCCCTGGAACCTGCTACATGTG
 CGGGCACAAAGGTACCTACGGGCTGCTGCGGCGGAGAGGACTGGCCCTCCCGGCTCCAGATGTTCTTC
 GCTAATAACCACGACCAGGAATTTGACCCTCAAAGGTTTACCCACCTGTCCCAGCTGAGAAGAGGAAGC
 CCATCCGGGTGCTGTCTCTTTGATGGAATCGCTACAGGGCTCCTGGTGTGAAGGACTTGGGCATTCA
 GGTGGACCGCTACATTGCCTCGGAGGTGTGTGAGGACTCCATCACGGTGGGCATGGTGGCCACCAGGGG
 AAGATCATGTACGTCGGGGACGTCGCAGCGTCACACAGAAGCATATCCAGGAGTGGGGCCATTTCGATC
 TGGTGATTGGGGCAGTCCCTGCAATGACCTCTCCATCGTCAACCCTGCTCGAAGGGCCTTACGAGGG
 CACTGGCCGGCTCTTCTTTGAGTTCTACCGCTCTGATGATGCGCGGCCCAAGGAGGGAGATGATCGC
 CCCTTCTTCTGGCTCTTTGAGAATGTGGTGGCCATGGGCGTTAGTGACAAGAGGGACATCTCGGATTTT
 TCGAGTCCAACCTGTGATGATTGATGCCAAAGAAGTGTGAGCTGCACACAGGGCCCGCTACTTCTGGGG
 TAACCTTCCCGGTATGAACAGGCCGTTGGCATCCACTGTGAATGATAAGCTGGAGCTGCAGGAGTGTCTG
 GAGCATGGCAGGATAGCCAAGTTCAGCAAAGTGGAGACCATTACTACGAGGTCAAACCTCCATAAAGCAGG
 GCAAAGACCAGCATTTTCTGTCTTCAATGAATGAGAAAGAGGACATCTTATGGTGCAGTGAATGGAAAG
 GGTATTTGGTTTCCAGTCCACTATACTGACGTCTCAACATGAGCCGCTTGGCGAGGCAGAGACTGCTG
 GGCCGGTCATGGAGCGTGCCAGTCATCCGCCACCTCTTCGCTCCGCTGAAGGAGTATTTTTCGCTGTGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG220906 representing NM_153759
 Red=Cloning site Green=Tags(s)

```


MGILERVRRNGRVDRLKDECDTAEKKAKVIAGMNAVEENQGPGESQKVEEASPPAVQQPTDPASPTVA
TTPEPVGSDAGDKNATKAGDDEPEYEDGRGFGIGELVWGKLRGFSWPGRIVSWWMTGRSRAAEGTRWVM
WFGDGKFSVVCVEKLMPLSSFFCSAFHQATYNKQPMYRKAIYEVLQVASSRAGKLPVCHDSDESDTAKAV
EVQNKPMIEWALGGFQPSGPKGLEPPEEENPYKEVYTDMMWVEPEAAAYAPPPPAKKPRKSTAEKPKVKE
IIDERTRERLVYEVQRKCRNIEDICISGSLNVTLEHPLFVGGMCQNCKNCFLECAQYDDDDGYQSYCTI
CCGGREVL MCGNNCCRCFCVECDVLLVGPAAQAAIKEDPWNCYMGCHKGTYGLLRREDWPSRLQMFF
ANNHDQEFDPKVPVPAEKRPPIRVLSLFDGIATGLLVKDLGIQVDRYIASEVCEDSITVGMVRHQG
KIMVYGDVRSVTQKHIQEWGPFDLVIGGSPCNDLSIVNPARKGLYEGTGRLFFEFYRLLHDARPEGDGR
PFFWLFENVVAMGVSDKRDISRFLSNPVMIDAKEVSAHRARYFWGNLPGMNRPLASTVNDKLELQECL
EHGRIAKFSKVRTITTRSNSIKQKQDQHFVFMNEKEDILWCTEMERVFVGFVHYTDVSNMSRLARQRL
GRSWSVPVIRHLFAPLKEYFACV
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:




ACCN: NM_153759

ORF Size: 2169 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_153759.3 , NP_715640.2
RefSeq Size:	3608 bp
RefSeq ORF:	2172 bp
Locus ID:	1788
UniProt ID:	Q9Y6K1
Cytogenetics:	2p23.3
Domains:	PWWP, DNA_methylase
Protein Families:	Druggable Genome
Protein Pathways:	Cysteine and methionine metabolism, Metabolic pathways
Gene Summary:	CpG methylation is an epigenetic modification that is important for embryonic development, imprinting, and X-chromosome inactivation. Studies in mice have demonstrated that DNA methylation is required for mammalian development. This gene encodes a DNA methyltransferase that is thought to function in de novo methylation, rather than maintenance methylation. The protein localizes to the cytoplasm and nucleus and its expression is developmentally regulated. [provided by RefSeq, Mar 2016]

Product images:



Circular map for RG220906