

Product datasheet for **RG220932**

LILRB4 (NM_001081438) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LILRB4 (NM_001081438) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LILRB4
Synonyms:	CD85K; HM18; ILT3; LILRB5; LIR-5; LIR5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG220932 representing NM_001081438
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATCCCCACCTTCACGGCTCTGCTCTGCCTCGGGCTGAGTCTGGGCCCCAGGACCCACATGCAGGCAG
 GGCCCTCCCAAACCCACCCTCTGGGCTGAGCCAGGCTCTGTGATCAGCTGGGGAACTCTGTGACCAT
 CTGGTGTAGGGGACCCTGGAGGCTCGGGAGTACCGTCTGGATAAAGAGGAAAGCCAGCACCCCTGGGAC
 AGACAGAACCCTGGAGCCCAAGAACAAGGCCAGATTCTCCATCCCATCCATGACAGAGGACTATGCAG
 GGAGATACCGCTGTTACTATCGCAGCCCTGTAGGCTGGTACAGCCAGTGACCCCTGGAGCTGGTGTG
 GACAGGAGCTACAGTAAACCCACCCTTTCAGCCCTGCCGAGTCTCTTGTGACCTCAGGAAAGAGCGTG
 ACCCTGCTGTGTCAGTACGGAGCCCAATGGACACTTTCCTTCTGATCAAGGAGCGGGCAGCCCATCCCC
 TACTGCATCTGAGATCAGAGCACGGAGCTCAGCAGCACCAGGCTGAATCCCCATGAGTCTGTGACCTC
 AGTGCACGGGGGACCTACAGGTGCTTCAGCTCACACGGCTTCTCCACTACCTGCTGTACACCCCACT
 GACCCCTGGAGCTCATAGTCTCAGGATCCTTGGAGGATCCCAGGCCCTACCCACAAGGTCCTGCTCAA
 CAGCTGCAGGCCCTGAGGACCAGCCCTCATGCCTACAGGGTCACTCCCCACAGTGGTCTGAGAAGGCA
 CTGGGAGGTACTGATCGGGTCTTGGTGGTCTCCATCCTGCTTCTCTCCCTCCTCCTCTCTCCTCCTC
 CAACACTGGCGTCAGGAAAACACAGGACATTGGCCCAGAGACAGGCTGATTTCCAACGTCTCCAGGGG
 CTGCCGAGCCAGAGCCCAAGGACGGGGCTACAGAGGAGGTCCAGCCAGCTGCTGACGTCCAGGGAGA
 AAATTCTGTGCTGCCGTGAAGAACACACAGCCTGAGGACGGGTGGAATGGACTCGGAGCCACAC
 GATGAAGACCCCAAGCAGTACGTATGCCAAGGTGAAACTCCAGACCTAGGAGAGAAATGGCCTCTC
 CTCCTCCCACTGTCTGGGGAATTCCTGGACACAAAGGACAGACAGGCAGAGGAGGACAGACAGATGGA
 CACTGAGGCTGCTGCATCTGAAGCCCCCAGGATGTGACCTACGCCAGCTGCACAGCTTTACCCTCAGA
 CAGAAGGCAACTGAGCCTCCTCCATCCCAGGAAGGGCCTCTCCAGCTGAGCCAGTGTCTATGCCACT
 TGGCCATCCAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG220932 representing NM_001081438
 Red=Cloning site Green=Tags(s)

MIPTFTALLCLGLSLGPRTHMQAGPLPKPTLWAEPGSVISWNSVTIWCQGTLEAREYRLDKEESPAPWD
 RQNPLEPKNKARFSIPSMTEYAGRYRCYRSPVQSPDPLELVMTGAYSKPTLSALPSPLVTSKSV
 TLLCQSRSPMDTFLLIKERAHPLLHLRSEHGAQQHQAEPMPVTSVHGGTYRCFSSHGFSHYLLSHPS
 DPLELIVSGSLEDPRPSPTRSVSTAAGPEDQPLMPTGSPHSGLRHWEVLIGLVVVSILLLSLLLFLL
 QHWRQGHRTLAQRQADFQRPPGAAEPEPKDGLQRRSSPAADVQGENFCAAVKNTQPEDGVEMDTRSPH
 DEDPQAVTYAKVHSRPRREMASPPSPLSGEFLDTKDRQAEEDRQMDTEAAASEAPQDVTYAQLHSFTLR
 QKATEPPPSQEGASPAEPSVYATLAIH

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_001081438.1, NP_001074907.1</p>
RefSeq Size:	<p>2094 bp</p>
RefSeq ORF:	<p>1343 bp</p>
Locus ID:	<p>11006</p>
Cytogenetics:	<p>19q13.42</p>
Protein Families:	<p>Transmembrane</p>
Gene Summary:	<p>This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). The receptor is expressed on immune cells where it binds to MHC class I molecules on antigen-presenting cells and transduces a negative signal that inhibits stimulation of an immune response. The receptor can also function in antigen capture and presentation. It is thought to control inflammatory responses and cytotoxicity to help focus the immune response and limit autoreactivity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p>