

Product datasheet for **RG221136**

Nogo A (RTN4) (NM_207521) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nogo A (RTN4) (NM_207521) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RTN4
Synonyms:	ASY; Nbla00271; Nbla10545; NI220/250; NOGO; NSP; NSP-CL; RTN-X; RTN4-A; RTN4-B1; RTN4-B2; RTN4-C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Protein Sequence: >RG221136 representing NM_207521
 Red=Cloning site Green=Tags(s)


MDLKEQPGNTISAGQEDFPSVLLETAASLPSLSPLSAASFKEHEYLGNLSTVLPTEGTLQENVSEASKEV
 SEKAKTLLIDRDLTEFSELEYSEMGSSFVSPKAESAVIVANPREEIIIVKNKDEEEKLVSNILHNQOEL
 PTALTKLVKEDVVSSEKAKDSFNEKRVAVEAPMREEYADFKPFERVWEVKDSKEDSDMLAAGGKIESNL
 ESKVDKCFADSLEQTNHEKDSSESSNDDTSFPSTPEGIKDRSGAYITCAPFNPAATESIATNIFPLLGD
 TSENKTDEKKIEEKKAQIVTEKNTSTKTSNPFLVAAQDSETDYVTTDNLTKVTEEVVANMPEGLTPDLVQ
 EACESELNEVTGKIAIYETKMDLVQTSVMEQESLYPAAQLCPSFEESEATPSPVLPDIVMEAPLNSAVPS
 AGASVIQPSSSPLEASSVNYESIKHEPENPPPYEEAMSVLKKVSGIKEEIKEPENINAALQETEAPYIS
 IACDLIKETKLSAEPAPDFSDYSEMAKVEQVPDPHSELVEDSSPDSEPVDLFSDSIPDVPQKQDETVML
 VKESLTETSFESMIEYENKEKLSALPPEGGKPYLESFKLSLNTKDTLLPDEVSTLSKKEKIPLQMEELS
 TAVYSNDDLFIKSAQIRETETFSOSSPIEIIDEFPTLISSKTDSFSLKAREYTDLEVSHKSEIANAPDG
 AGSLPCTELPHDLSLKNIQPKVEEKISFSDDFSKNGSATSKVLLLPDVSALATQAEIESIVKPKVLVKE
 AEKLPDTEKEDRSPSAIFSAELSKTSVVDLLYWRDIKKTGVVFGASLFLLLSLTVFSIVSVTAYIALA
 LLSVTISFRIYKGIQAIQKSDGHPFRAYLESEVAISEELVQKYSNSALGHVNCIKELRRLFLVDDL
 DSLKFAVLMWVFTYVGFALFNGLTLLILALISLFSVPVIYERHQAQIDHYLGLANKNVKDAMAKIQAKIPG
 LKRKAE

TRTRPLE – GFP Tag – V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_207521

ORF Size: 2958 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_207521.2](#)

RefSeq Size: 4129 bp

RefSeq ORF: 2961 bp

Locus ID: 57142

UniProt ID: [Q9NQC3](#)

Cytogenetics: 2p16.1

Protein Families: Transmembrane

Gene Summary: This gene belongs to the family of reticulon encoding genes. Reticulons are associated with the endoplasmic reticulum, and are involved in neuroendocrine secretion or in membrane trafficking in neuroendocrine cells. The product of this gene is a potent neurite outgrowth inhibitor which may also help block the regeneration of the central nervous system in higher vertebrates. Alternatively spliced transcript variants derived both from differential splicing and differential promoter usage and encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]