

Product datasheet for **RG221330**

GJA5 (NM_181703) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GJA5 (NM_181703) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	GJA5
Synonyms:	ATFB11; CX40
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG221330 representing NM_181703 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCGATTGGAGCTTCTGGGAAATTCCTGGAGGAAGTACACAAGCACTCGACCGTGGTAGGCAAGG
TCTGGCTCACTGTCTTTCATATTCGGTATGCTCGTGGGCACAGCTGCTGAGTCTTCTGGGGGA
TGAGCAGGCTGATTTCCGGTGTGATACGATTACGCTGGCTGCCAGAATGTCTGCTACGACCAGGCTTTC
CCCATCTCCACATTCGCTACTGGGTGCTGCAGATCATCTTCGCTCCACGCCCTCTCTGGTGTACATGG
GCCACGCCATGCACACTGTGCGCATGCAGGAGAAGCGCAAGCTACGGGAGGCCGAGAGGGCCAAAGAGGT
CCGGGGCTCTGGCTCTTACGAGTACCCGGTGGCAGAGAAGGCAGAACTGTCTGCTGGGAGGAAGGGAAT
GGAAGGATTGCCCTCCAGGGCACTCTGCTCAACACCTATGTGTGCAGCATCCTGATCCGACCACCATGG
AGGTGGGCTTCATTGTGGGCCAGTACTTCATCTACGGAATCTTCTGACCACCCTGCATGTCTGCCGAG
GAGTCCCTGTCCCAACCCGGTCAACTGTTACGTATCCCGGCCACAGAGAAGAATGTCTTTCATTGTCTTT
ATGCTGGCTGTGGCTGCACTGTCCCTCCTCCTTAGCCTGGCTGAACTCTACCACCTGGGCTGGAAGAAGA
TCAGACAGCGATTTGTCAAACCGCGGCAGCACATGGCTAAGTGCCAGCTTCTGGCCCTCTGTGGGCAT
AGTCCAGAGCTGCACACCACCCCGACTTTAATCAGTGCCTGGAGAATGGCCCTGGGGGAAAATTCCTC
AATCCCTCAGCAATAATATGGCCTCCCAACAAAACACAGACAACCTGGTCACCGAGCAAGTACGAGGTC
AGGAGCAGACTCCTGGGAAGGTTTCATCCAGGTTGTTATGGCCAGAAGCCTGAGGTGCCCAATGGAGT
CTCACCAGGTCACCGCCTTCCCATGGCTATCATAGTGACAAGCGACGTCTTAGTAAAGGCCAGCAGCAAG
GCAAGGTCAGATGACCTATCAGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG221330 representing NM_181703
 Red=Cloning site Green=Tags(s)

MGDWSFLGNFLEEVEHKHSTVVGKVLTVLFI FRMLVLGTA AESSWGDEQADFRCDTIQPGCQNV CYDQAF
 PISHIRYWVLQIIFVSTPSLVYMGHAMHTVRMQEKRLREAERAKEVRGSGSYEYPVAEKAE LSCWEEGN
 GRIALQGTLLNTYVCSILIRTTMEVGFIVGQYFIYGIFL TTLHVCRRSPCPHPVNCYVSRPTEKNV FIVF
 MLAVAALSLLLSLAELYHLGWKKIRQRFVKPRQHMAKQ LSGPSVGI VQSC TPPPDFNQCLENGPGGKFF
 NPFSNNMASQONTDLNLTVEQVRGQEQTPEGFIQVRYGQKPEV PNVSPGHRLPHGYHSDKRRLSKASSK
 ARSDDL SV

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_181703

ORF Size: 1074 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181703.4](#)

RefSeq Size: 2566 bp

RefSeq ORF: 1077 bp

Locus ID: 2702


UniProt ID: [P36382](#)

Cytogenetics: 1q21.2

Protein Families: Transmembrane

Gene Summary: This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. Mutations in this gene may be associated with atrial fibrillation. Alternatively spliced transcript variants encoding the same isoform have been described. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG221330