

Product datasheet for **RG221348**

CTAGE5 (MIA2) (NM_203355) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTAGE5 (MIA2) (NM_203355) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CTAGE5
Synonyms:	CTAGE5; MEA6; MGEA; MGEA6; MGEA11; TALI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG221348 representing NM_203355
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGAGCCCGGGTTACCCCTCAACCGTATTTGGGGCTGCTCCTGGAGGAGCTACGCAGGGTTGTGG
 CAGCACTGCCTGAAGGTATGAGACCAGATTCTAATCTTTATGGTTTTCCATGGGAATTGGTGATATGTGC
 AGCTGTTGTTGGATTTTTGCTGTTCTCTTTTTTTGTGGAGAAGTTTTAGATCGGTTAGGAGTCGGCTT
 TATGTGGGACGAGAGAAAAAGCTTGCTCTAATGCTTTCTGGACTAATTGAAGAAAAAGTAACTACTTG
 AAAAAATTAGCCTTGTTCAAAAAGAGTATGAAGGCTATGAAGTAGAGTCATCTTTAAAGGATGCCAGCTT
 TGAGAAGGAGGCAACAGAACACAAAGTTGGAGGCAACCTGTGAAAAGCTGAACAGGTCCAATTCTGAA
 CTTGAGGATGAAATACTGTCTAGAAAAAGATTAAAAGAAGAGAAATCCAACATTCTGAACAAGATG
 AATTGATGGCGGATATTTCAAAAAGGATACAGTCTCTAGAAGATGAGTCAAAATCCCTCAAATCACAAAGT
 AGCTGAAGCCAAAATGACCTTCAAGATATTTCAAATGAATGAAGAACGACTGAAGATAGCAATAAAAGAT
 GCTTTGAATGAAAATTCTCAACTTCAGGAAAGCCAGAAACAGCTTTTGCAAGAAGCTGAAGTATGAAAAG
 AACAAAGTGAGTGAACCTAATAAACAGAAAGTAACATTTGAAGACTCCAAAGTACATGCAGAACAAAGTCT
 AAATGATAAAGAAAGTCACATCAAGACTCTGACTGAACGCTTGTTAAAGATGAAAGATTGGGCTGCATG
 CTTGGAGAAGACATAACGGATGATGATAACTTGGAAATAGAAAAGAACAGTGAATCGGAAAATGGTGCTT
 ACTTAGATAATCCTCCAAAAGGAGCTTTGAAGAAACTGATTATGCTGCTAAGTTAAATGCTTCTTTAAA
 AACCTTAGAAGGAGAAAGAAACCAATTTATATTCAGTTGTCTGAAGTTGATAAAAACAAAGGAAGAGCTT
 ACAGAGCATATTAATACTTCAGACTGAACAAGCATCTTGCAGTCAGAAAACACACATTTTAAAAATG
 AGAATCAGAAGCTTCAACAGAACTTAAAGTAATGACTGAATTATCAAGAAAAATGAAATGAAACTCCA
 CAGGAAAATTAACAGTAGAGGAAAATTATCGGTTAGAGAAAAGAAGAGAAACTTTCTAAAGTAGATGAAAAG
 ATCAGCCATGCCACTGAAGAGCTGGAGACCTATAGAAAAGCGAGCCAAAAGATCTTGAAGAAGAAATGGAGA
 GAACTATTCATTCTATCAAGGGCAGATTATTTCCCATGAGAAAAAGCACATGATAATTGGTTGGCAGC
 TCGGAATGCTGAAAGAACTCAATGATTTAAGGAAAGAAAATGCTCACAAACAGACAAAAATTAAGTAA
 ACAGAGCTTAAATTTGAACTTTTAGAAAAAGATCCTTATGCACTCGATGTTCCAAATACAGCATTGGCA
 GAGGCTCAGGAGCCAGGGAATCCTCTGGACCATCAGATTACCAATGAAAGAGGAGAATCAAGCTGTGA
 TAGGTTAACCGATCCTCATAGGGCTCCCTGACTGGGTCTGTGTACCTCCATGGGACCAGGACCGT
 AGGATGATGTTTCTCCGCCAGGACAATCATATCCTGATTACGCCCTTCTCCACAAAGGCAAGACAGAT
 TTTGTTCTAATCTGGTAGACTGTCTGGACCAGCAGAAGTCAAGATTTTAAATATGCCTTCTTTGGATAA
 AATGGATGGGTCAATGCCTTCAGAAATGGAATCCAGTAGAAATGATACCAAAGATGATCTTGGTAATTTA
 AATGTGCCTGATTCATCTCTCCCTGCTGAAAATGAAGCCACTGGCCCTGGCTTTGTTCTCCACCTCTTG
 CTCCAATCAGAGGTCCATTGTTTCCAGTGGATGCAAGAGGCCATTTCTTGAGAAGAGGACCTCCTTTCCC
 CCCACCTCCTCCAGGAGCATGTTTGGAGCTTCTCGAGATTATTTTCCACCAGGGGATTTCCAGGTCCA
 CCACCTGCTCCATTTGCAATGAGAAATGTCTATCCACCAGGGGTTTTCTCCTTACCTTCCCCCAAGAC
 CTGGATTTTTCCCCCACCCACATTTGAAGGTAGAAGTGAAGTCCCCTCAGGTTTGATTCCACCTTC
 AAATGAGCCTGCTACTGAACATCCAGAACCACAGCAAGAAACC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG221348 representing NM_203355
Red=Cloning site Green=Tags(s)

```
MEEPGVTPQPYLGLLLEELRRVVAALPEGMRPDSNLYGFPWELVICA AVVGVFFAVLFFLWRSFRSVRSRL
YVGREKKLALML SGLIEEKSLLLEKFSLVQKEYEGYEVESSLK DASFKEATEAQSLEATCEKLNRSNSE
LEDEILCLEKELKEEKSKHSEQDEL MADISKRIQSLEDESKSLKSQVAEAKMTFKIFQMNEERLKI AIKD
ALNENSQ LQESQKQLLQEA EVWKEQVSELNKQKVTFFEDSKVHAEQV LNDKESH IKT LTERLLKMKDWAAM
LGEDITDDNLELEMNSENGAYLDNPPKGALKKLIHAAKLNASLKTLEGERNQIYIQLSEVDKTEEL
TEHIKNLQTEQASLQSENTHFENENQKLQKLVMTL YQENEMKLHRKLTVEENYRLEKEEKLSKVDEK
ISHATEELETYRKRKADLEEEELERTIHSYQGQII SHEKKAHDNWL AARNAERNLNDLRKENAHRQKLTE
TELKFELLEKDPYALDVPNTAFGRGSRGPNLDHQITNERGESSCDRL TDPHRAPSDTGSLSPPWDQDR
RMMFPPPGQSYPSALPPQRQDRFCNSNGRLSGPAELRSFNMPSLDKMDGSMPSSEMSSRNDTKDDLGNL
NVPDSSLPAENEATGPGFVPPPLAPIRGPLFPVDARGPFLRRGPPPPPPPPGAMFGASRDYFPPGDFPGP
PPAPFAMRNVYPPRGFPYLP RPFGFPFPPPHSEGRSEFP SGLIPPSNEPATEHPEPQQT
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_203355

ORF Size: 2283 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_203355.3](#)

RefSeq Size: 3551 bp

RefSeq ORF: 2286 bp


Locus ID: 4253

Cytogenetics: 14q21.1

Protein Families: Transmembrane

Gene Summary: This gene encodes s receptor in the endoplasmic reticulum, which plays a role in the export of large pre-chylomicrons and pre-very low density lipoproteins (pre-VLDLs). Three major classes of transcripts are generated from this gene- melanoma inhibitory activity 2-specific transcripts, cTAGE family member 5-specific transcripts and transcripts that include exons from both these transcript species (TANGO1-like or TALI). Additionally, alternative splicing in these transcripts results in multiple transcript variants encoding multiple isoforms. [provided by RefSeq, Sep 2016]

Product images:



Circular map for RG221348