

## Product datasheet for **RG221452**

### HSD11B1L (NM\_198707) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** HSD11B1L (NM\_198707) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** HSD11B1L  
**Synonyms:** 11-beta-HSD3; 11-DH3; HSD1L; HSD3; SCDR10; SCDR10B; SDR26C2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG221452 representing NM\_198707  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCTCCCCTGAGGCGCCGAGAGCGTGGTGCAGTTTTCGCTGGACAAGCTGGCGGGCTGGACTACC  
TCGTGCTGAACCACATCGGCGGCCCGCCGCGCACGCGAGCCCGCAGCCCCAGGCAACTCGCTGGCT  
CATGCAGTAACTTTGTGAGCTACGTGCAACTGACGTCGCGGGCGCTGCCAGCCTGACGGACAGCAAG  
GGCTCCCTGGTGGTGGTGTCTCGCTGCTCGGCCGCTGCCACGTCGTTCTCCACTCCCTACTCGGCGG  
CCAAGTTTTCGCTGGACGGCTTCTCGGCTCCCTGCGGCGGGAGCTGGACGTGCAGGACGTGAACGTGGC  
CATCACCATGTGCGTCTGGGCCTCCGAGATCGCGCCTCCGCGCCGAGGAGTCAAGGGAGTACAGGAGG  
GTCAAGGCGGCCCGGGGCCAAGGCAGCCCTGGCCGTGATCCGCGGCGGCCACGCGCGCGGCCGCGG  
TCTTCTACCCGTGGCGTTTCCGCTGCTGTGCTTCTCGGCGCTGGCTACCGCGCCCGCGGGCCTGGTT  
TATCCGCCAGGAGCTCAACGTACCGCCGCGGACGCC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >RG221452 representing NM\_198707  
Red=Cloning site Green=Tags(s)

MASPEAPESVVQFALDKLGGLDYLVLNHIGGAPAGTRARSPQATRWLMQVNFVSYVQLTSRALPSLTDSK  
GSLVVVSSLLGRVPTSFSTPYSAAKFALDGFFGSLRRELDVQDVNVAITMCVGLRDRASAAEAVRGVTR  
VKAAPGPKAALAVIRGGATRAAGVFYPWRFRLCLLRRWLPRPRAWFIRQELNVTAAAA

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_198707

**ORF Size:** 597 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_198707.3](#)

**RefSeq Size:** 1897 bp

**RefSeq ORF:** 600 bp

**Locus ID:** 374875


**UniProt ID:** [Q7Z5J1](#)

**Cytogenetics:** 19p13.3

**Protein Families:** Druggable Genome

**Gene Summary:** This gene is a member of the hydroxysteroid dehydrogenase family. The encoded protein is similar to an enzyme that catalyzes the interconversion of inactive to active glucocorticoids (e.g. cortisone). Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2012]

**Product images:**



Circular map for RG221452