

## Product datasheet for **RG221589**

### Caspase-7 (CASP7) (NM\_033338) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Caspase-7 (CASP7) (NM_033338) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CASP7
Synonyms:	CASP-7; CMH-1; ICE-LAP3; LICE2; MCH3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG221589 representing NM_033338 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGACTGTGTTGGTTGGCCTCCAGGCAGGAAGTGGCACTTGAAAAGAACACCAGCTGCGGTGGTAGCA  
GTGGGATTTGTGCTTCTTATGTTACCCAGATGGCAGATGATCAGGGCTGTATTGAAGAGCAGGGGTTGA  
GGATTCAGCAAATGAAGATTCAGTGGATGCTAAGCCAGACCGGTCCTCGTTGTACCGTCCCTCTTCAGT  
AAGAAGAAGAAAATGTCACCATGCGATCCATCAAGACCACCCGGGACCGAGTGCCTACATATCAGTACA  
ACATGAATTTTAAAAGCTGGGCAAATGCATCATAATAAACAACAAGAAGCTTTGATAAAGTGACAGGTAT  
GGGCGTTGAAAACGGAACAGACAAAGATGCCGAGGCGCTCTTCAAGTGCCTCCGAAGCCTGGGTTTTGAC  
GTGATTGTCTATAATGACTGCTCTTGTGCCAAGATGCAAGATCTGCTTAAAAAGCTTCTGAAGAGGACC  
ATACAAATGCCGCTGCTTCGCCTGCATCCTCTTAAGCCATGGAGAAGAAAATGTAATTTATGGGAAAGA  
TGGTGTACACCAATAAAGGATTTGACAGCCACTTTAGGGGGGATAGATGCAAAACCCCTTTAGAGAAA  
CCCAAATCTTCTTCAATCAGGCTTGCCGAGGGACCGAGCTTGATGATGGCATCCAGGCCGACTCGGGGC  
CCATCAATGACACAGATGCTAATCCTCGATAACAAGATCCCAGTGGAAAGCTGACTTCTCTTCGCCTATTC  
CACGGTTCCAGGCTATTACTCGTGGAGGAGCCAGGAAGAGGCTCCTGGTTTGTGCAAGCCCTCTGCTCC  
ATCCTGGAGGAGCACGGAAAAGACCTGGAAATCATGCAGATCCTCACCAGGTGAATGACAGAGTTGCCA  
GGCACTTTGAGTCTCAGTCTGATGACCCCACTTCCATGAGAAGAAGCAGATCCCCTGTGTGCTCCAT  
GCTCACCAAGGAAGTCTACTTCAGTCAA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG221589 representing NM\_033338  
Red=Cloning site Green=Tags(s)

MDCVGGWPPGRKWHLEKNTSCGGSSGICASYVTQMADDQGCIEEQVEDSANEDSVDAKPDRSSFVPSLFS  
 KKKKNVTMRSIKTTRDRVPTYQYNMNFELGKCIINNNKFDKVTGMGVRNGTDKDAEALFKCFRSLGFD  
 VIYVNDSCAKMQDLLKKASEEDHTNAACFACILLSHGEENVYIGKDGVTPIKDLTAHFRGDRCKTLLEK  
 PKLFFIQACRGTLEDGDIQADSGPINDTDANPRYKIPVEADFLFAYSTVPGYYSWRSPGRGSWFVQALCS  
 ILEEHGKDLIMQILTRVNDRVARHFESQSDDPHFHEKKQIPCVVSMMLTKELYFSQ


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_033338

**ORF Size:** 1008 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_033338.6</a>
<b>RefSeq Size:</b>	2712 bp
<b>RefSeq ORF:</b>	1011 bp
<b>Locus ID:</b>	840
<b>UniProt ID:</b>	<a href="#">P55210</a>
<b>Cytogenetics:</b>	10q25.3
<b>Protein Families:</b>	Druggable Genome, Protease
<b>Protein Pathways:</b>	Alzheimer's disease, Apoptosis
<b>Gene Summary:</b>	This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. The precursor of the encoded protein is cleaved by caspase 3 and 10, is activated upon cell death stimuli and induces apoptosis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, May 2012]