

Product datasheet for **RG221653**

CAMK2G (NM_172171) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CAMK2G (NM_172171) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CAMK2G
Synonyms:	CAMK; CAMK-II; CAMKG; MRD59
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG221653 representing NM_172171
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCACCACCGCCACCTGCACCCGTTTCACCGACGACTACCAGCTCTTCGAGGAGCTTGCAAGGGTG
 CTTTCTCTGTGGTCCGCAGGTGTGTGAAGAAAACCTCCACGCAGGAGTACGCAGCAAAAATCATCAATAC
 CAAGAAGTTGTCTGCCCGGGATCACCAGAACTAGAACGTGAGGCTCGGATATGTCGACTTCTGAAACAT
 CCAAACATCGTGCCTCCATGACAGTATTTCTGAAGAAGGGTTTCACTACCTCGTGTGGACCTTGTTA
 CCGGCGGGGAGCTGTTTGAAGACATTGTGGCCAGAGAGTACTACAGTGAAGCAGATGCCAGCCACTGTAT
 ACATCAGATTCTGGAGAGTGTAAACCACATCCACCAGCATGACATCGTCCACAGGGACCTGAAGCCTGAG
 AACCTGCTGTGGCGAGTAAATGCAAGGGTGCCCGCTCAAGCTGGCTGATTTTGGCCTAGCCATCGAAG
 TACAGGGAGAGCAGCAGGCTTGGTTTGGTTTTGCTGGCACCCAGGTTACTTGTCCCTGAGGTCTTGAG
 GAAAGATCCCTATGGAACCTGTGGATATCTGGCCTGCGGGGTATCTGTATATCCTCTGGTGGGC
 TATCCTCCCTTCTGGGATGAGGATCAGCACAAGCTGTATCAGCAGATCAAGGCTGGAGCCTATGATTTC
 CATCACCAGAATGGGACACGGTAACTCCTGAAGCCAAGAATTGATCAACCAGATGCTGACCATAAACCC
 AGCAAAGCGCATCACGGCTGACCAGGCTCTCAAGCACCCGTGGTCTGTCAACGATCCACGGTGGCATCC
 ATGATGCATCGTCAGGAGACTGTGGAGTGTTCGCAAGTTCAATGCCCGGAGAAAACCTGAAGGGTGCCA
 TCCTCACGACCATGCTTGTCTCCAGGAACCTCTCAGCTGCCAAAAGCCTATTGAACAAGAAGTCGGATGG
 CGGTGTCAAGCCACAGAGCAACAACAAAACAGTCTCGTAAGCCAGCCCAAGAGCCCGCCCTTGCGAG
 ACGGCCATGGAGCCAAACACTGTGGTACACAACGCTACAGATGGGATCAAGGGCTCCACAGAGAGCT
 GCAACACCACCACAGAAGATGAGGACCTCAAAGCTGCCCGCTCCGCACCTGGGAATGGCAGCTCGGTGCC
 TGAAGGACGGAGCTCCCGGGACAGAACAGCCCTCTGCAGGCATGCAGCCCGAGCCTTCTCTGCTCC
 TCAGCCATGCGAAAACAGGAGATCATTAAAGATTACAGAACAGCTGATTGAAGCCATCAACAATGGGGACT
 TTGAGGCCTACACGAAGATTTGTGATCCAGGCCTCACTTCCTTTGAGCCTGAGGCCCTTGTTAACCTCGT
 GGAGGGGATGGATTTCCATAAGTTTACTTTGAGAATCTCCTGTCCAAGAACAGCAAGCCTATCCATACC
 ACCATCTAAACCCACACGTCCACGTGATTGGGGAGGACGCAGCGTGCATCGCTACATCCGCTCACCC
 AGTACATCGACGGCAGGGTCGGCCTCGCACCAGCCAGTCAAGAGACCCGGTCTGGCACCGTCGGGA
 TGGCAAGTGGCTCAATGTCCACTATCACTGCTCAGGGGCCCTGCCGCACCGCTGCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG221653 representing NM_172171
 Red=Cloning site Green=Tags(s)

MATTATCTRFTDDYQLFEELGKGAFSVVRRRCVKKTSTQEYAAKIINTKKLSARDHQKLEREARICRLLKH
 PNIVRLHDSISEEGFHYLVFDLVTGGELFEDIVAREYYSEADASHCIHQILESVNHIHQHDIVHRDLKPE
 NLLLASKCKGAAVKLADFLAIEVQGEQQAWFGFAGTPGYLSPEVLRKDPYGKPVDIWACGVILYILLVG
 YPPFWDEDQHKLYQQIKAGAYDFPSPEWDTVTPEAKNLINQMLTINPAKRITADQALKHPWVCQRSTVAS
 MMHRQETVECLRKFNARRKLKGAAILTTMLVSRNFSAAKSLLNKSDGGVQPQSNNKNSLVSPAQEPAPLQ
 TAMEPQTTVVHNATDGIKSTESCENTTTEDELKAAPLRTGNGSSVPEGRSSRDRTAPSAGMQPQPSLCS
 SAMRKQEI IKITEQLIEAINNGDFEAYTKICDPGLTSFEPEALGNLVEGMDFHKFYFENLLSKNSKPIHT
 TILNPHVHVIGEDAACIAYIRLTQYIDGQGRPRTSQSEETRVVHRRDGGKWLNVVHYHCSGAPAAPLQ

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_172171

ORF Size: 1668 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_172171.2](#)
RefSeq Size: 3818 bp
RefSeq ORF: 1671 bp
Locus ID: 818
Cytogenetics: 10q22.2
Protein Families: Druggable Genome, Protein Kinase
Protein Pathways: Calcium signaling pathway, ErbB signaling pathway, Glioma, GnRH signaling pathway, Long-term potentiation, Melanogenesis, Neurotrophin signaling pathway, Olfactory transduction, Oocyte meiosis, Wnt signaling pathway
Gene Summary: The product of this gene is one of the four subunits of an enzyme which belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. Calcium signaling is crucial for several aspects of plasticity at glutamatergic synapses. In mammalian cells the enzyme is composed of four different chains: alpha, beta, gamma, and delta. The product of this gene is a gamma chain. Many alternatively spliced transcripts encoding different isoforms have been described but the full-length nature of all the variants has not been determined.[provided by RefSeq, Mar 2011]

Product images:



Circular map for RG221653