

Product datasheet for **RG221681**

MAP2 (NM_001039538) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAP2 (NM_001039538) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MAP2
Synonyms:	MAP-2; MAP2A; MAP2B; MAP2C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG221681 representing NM_001039538
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCAGATGAACGGAAAGATGAAGCAAAGGCACCTCACTGGACCTCAGCACCGCTAACAGAGGCATCTG
 CACACTCACATCCACCTGAGATTAAGGATCAAGGCGGAGCAGGGGAAGGACTTGTCCGAAGCGCCAATGG
 ATTCCCATACAGGGAGGATGAAGAGGGTGCCTTTGGAGAGCATGGGTACAGGGCACCTATTCAAATACC
 AAAGAGAATGGGATCAACGGAGAGCTGACCTCAGCTGACAGAGAAAACAGCAGAGGAGGTGTCTGCAAGGA
 TAGTTCAAGTAGTCACTGCTGAGGCTGTAGCAGTCTGAAAGGTGAACAAGAGAAAAGACTCAACATAA
 AGACCAGACTGCAGCTCTGCCTTAGCAGCTGAAGAAACAGCTAATCTGCCTCCTTCTCCACCCCATCA
 CCTGCCTCAGAACAGACTGTCACAGTGGAGGAAGCAGCAGGTGGGAATCAGCTCTGGCTCCCAGTGTAT
 TAAACAGGCAAAGGACAAAGTCTTAATTCTACCTTGCAAAGATTCTGCTTTACAGGGTAGCACAAA
 GTCCCCAAGATACAGCTCAGCCTGCCCTAGCACGACTAAAAGGGTACATTTTCTGACAGTTTATTAATA
 CAGCCCACCTCAGCAGGCTCCACAGACCGTTTGCATACCTCAAAAACAGGGAACAAGGACGGAGTAACCA
 AGAGCCCAGAAAAGCGCTCTTCTCTCCCAAGACCTTCTCCATTCTCCCTCCTCGGCGAGGTGTGTCAAG
 AGACAGAGATGAGAATTCCTTCTCTCAACAGTTCTATCTTCTTTCAGCACGGCGGACCACAGGTCA
 GAGCCAATTCGCAGAGCAGGGAAGAGTGGTACCTCAACACCCACTACCCCTGGGTCTACTGCCATCACTC
 CTGGCACCCACCAAGTTATTCTTACGCACACCAGGCACTCTGGAACCCCTAGCTATCCCAGGACCCC
 TCACACACCAGGAACCCCAAGTCTGCCATCTTGGTCCGAGTGAGAAGAAGGTCCGCATCATACGTA
 CCTCCAAAATCTCCTCGACTCCCAAGCAGCTTCGGCTTATTAACCAACCACTGCCAGACCTGAAGAATG
 TCAAAATCCAAAATCGGATCAACAGACAACATCAAATACCAGCCTAAAGGGGGCAGGTTAGGATTTAAA
 CAAGAAGATCGATTTTAGCAAAGTTTCAAGTCCAGATGTGGTTCCAAGGATAACATCAAACATTCGGCTGGG
 GCGGAAATGTACAAATTGTTACCAAGAAAATAGACCTAAGCCATGTGACATCCAATGTGGCTCTCTGA
 AGAACATCCGCCACAGGCCAGGTGGCGGACGTGTGAAAATTGAGAGTGTAAAAGTAAAGATTTCAAAGAAA
 GGCCCAAGCTAAAGTTGGTTCTCTTGATAATGCTCATCATGTACCTGGAGGTGGTAAATGTCAAGATTGAC
 AGCCAAAAGTTGAACTTCAGAGAGCATGCTAAAGCCCGTGTGGACCATGGGGCTGAGATCATTACACAGT
 CCCCAGGCAGATCCAGCTGGCATCACCCGACGACTCAGCAATGTCTCCTCGTCTGGAAGCATCAACCT
 GCTCGAATCTCCTCAGCTTGCCACTTTGGCTGAGGATGCTACTGCTGCACTCGCTAAGCAGGCTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG221681 representing NM_001039538
 Red=Cloning site Green=Tags(s)

MADERKDEAKAPHWTSAPL TEASAHSHPPEIKDQGGAGEGLVRSANGFPYREDEEGAFGEHGSQGTYSNT
 KENGINDELTSADRETAEEVSARIVQVVTAEAVAVLKGEQEKEAQHKDQT AALPLAAEETANLPPSPPPS
 PASEQTVTVVEEAAGGESALAPSVFKQAKDKVSNSTLSKIPALQGSTKSPRYSSACPSTTKRATFSDSLLI
 QPTSAGSTDRLPYKSGNKDGVTKSPEKRSSLPRPSSILPPRRGVSGDRDENSFSLNSSISSARRTTRS
 EPIRRAGKSGTSTPTTPGSTAITPGTPPSYSSRTPGTPGTPSYPRTPHTPGTPKSAILVPSEKKVAIIRT
 PPKSPATPKQLRLINQPLPDLKNVSKSIGSTDNIKYQPKGGQVIRLNKKIDF SKVQSRCSKDNIKHSAG
 GGNVQIVTKKIDLSHVTSKCGSLKNIRHRPGGGRVKIESVKLDFKEKAQAKVGLDNAHHVPGGNNVKID
 SQKLNFRHAKARVDHGAEIITQSPGRSSVASPRRLSNVSSSGSINLLESPQLATLAEDVTAALAKQGL

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001039538

ORF Size: 1677 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039538.2](#)

RefSeq Size: 5844 bp

RefSeq ORF: 1680 bp

Locus ID: 4133


UniProt ID: [P11137](#)

Cytogenetics: 2q34

Protein Families: Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS

Gene Summary: This gene encodes a protein that belongs to the microtubule-associated protein family. The proteins of this family are thought to be involved in microtubule assembly, which is an essential step in neurogenesis. The products of similar genes in rat and mouse are neuron-specific cytoskeletal proteins that are enriched in dendrites, implicating a role in determining and stabilizing dendritic shape during neuron development. A number of alternatively spliced variants encoding distinct isoforms have been described. [provided by RefSeq, Jan 2010]

Product images:



Circular map for RG221681