

Product datasheet for **RG221740**

CPT1A (NM_001876) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CPT1A (NM_001876) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CPT1A
Synonyms:	CPT1; CPT1-L; L-CPT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG221740 representing NM_001876
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGAAGCTACCAAGCTGTGGCCTTTCAGTTCACGGTCACTCCGGACGGGATTGACCTGCGGCTGA
 GCCATGAAGCTCTTAGACAAATCTATCTCTGGACTTCATTCTGGAAAAAGAAGTTCATCAGATTCAA
 GAACGGCATCATCACTGGCGTGTACCCGCAAGCCCTCCAGTTGGCTTATCGTGGTGGTGGCGTGATG
 ACAACGATGTACGCCAAGATCGACCCCTCGTTAGGAATAATTGCAAAAATCAATCGGACTCTGAAACGG
 CCAACTGCATGTCCAGCCAGACGAAGAAGCTGGTCAAGCGCGTGTGTTGGCACCGGCCTGTGGGTGGC
 CCTCATCGTCACCATGCGCTACTCCCTGAAAGTGTGCTCTCTACCACGGGTGGATGTTCACTGAGCAC
 GGCAAGATGAGTCGTGCCACCAAGATCTGGATGGGTATGGTCAAGATCTTTTCAGGCCGAAAACCCATGT
 TGTACAGCTTCCAGACATCGCTGCCTCGCTGCCGGTCCCAGGCTGTCAAAGACACTGTGAACAGGTATCT
 ACAGTCGGTGAGGCTCTTATGAAGGAAGAAGACTTCAAACGGATGACAGCACTTGTCAAGATTTTGTCT
 GTCGGTCTTGACCAAGATTACAGTGGTATTTGAAGTTAAAAATCTGGTGGGCTACAAAATACGTGAGCG
 ACTGGTGGGAGGAGTACATCTACCTCCGAGGACGAGGGCCGCTCATGGTGAACAGCAACTATTATGCCAT
 GGATCTGCTGTATATCTTCCAACCTCACATTCAGGCAGCAAGAGCCGGCAACGCCATCCATGCCATCCTG
 CTTTACAGGCGCAAACCTGGACCGGGAGGAAATCAAACCAATTCGTCTTTTGGGATCCACGATCCACTCT
 GCTCCGCTCAGTGGGAGCGGATGTTAATACTTCCCGGATCCCAGGAGAGGAGACAGACACCATCCAGCA
 CATGAGAGACAGCAAGCACATCGTCGTGTACCATCGAGGACGCTACTTCAAGGTCTGGCTCTACCATGAT
 GGGCGGCTGTGAAGCCCCGGGAGATGGAGCAGCAGATGCAGAGGATCCTGGACAATACCTCGGAGCCTC
 AGCCCCGGGAGGCCAGGCTGGCAGCCCTCACCGCAGGAGACAGAGTTCCCTGGCCAGGTCTGTCAGGC
 CTATTTTGGACGTGGGAAAAATAAGCAGTCTCTTGATGCTGTGGAGAAAGCAGCGTTCCTCGTGACGTTA
 GATGAAACTGAAGAAGGATACAGAAGTGAAGACCCGGATACGTCATGGACAGCTACGCCAAATCTCTAC
 TACACGGCCGATGTTACGACAGGTGGTTTGACAAGTCGTTACGTTTGTGTCTTCAAAAACGGGAAGAT
 GGGCCTCAACGCTGAACACTCCTGGGCAGATGCGCCGATCGTGGCCACCTTTGGGAGTACGTATGTCC
 ATTGACAGCCTCCAGCTGGGCTATGCGGAGGATGGGCACTGCAAAGGCGACATCAATCCGAACATCCGT
 ACCCCACCAGGCTGCAGTGGGACATCCCAGGGGAAATGTCAAGAGGTTATAGAGACCTCCCTGAACACCGC
 AAATCTTCTGGCAAACGACGTGGATTTCCATTCTTCCATTCTAGCCTTTGGTAAAGGAATCATCAAG
 AAATGTGCGACGAGCCCAGACGCTTTGTGACGCTGGCCCTCCAGCTGGCGCACTACAAGGACATGGGCA
 AGTTTTGCCTCACATACGAGGCCTCCATGACCCGGCTCTCCGAGAGGGGAGGACGGAGACCGTGCCTC
 CTGCAACACTGAGTCATGCGACTTCGTGCGGGCCATGGTGGACCCGGCCAGACGGTGAACAGAGGCTG
 AAGTTGTTCAAGTTGGCGTCTGAGAAGCATCAGCATATGTATCGCTCGCCATGACCCGGCTCTGGGATCG
 ATCGTCACCTCTTCTGCCTTACGTGGTGTCTAAATATCTCGCTGTGGAGTCCCCCTTCTTAAGGAAGT
 TTTATCTGAGCCTTGGAGATTATCAACAAGCCAGACCCCTCAGCAGCAAGTGGAGCTGTTGACTTGGAG
 AATAACCCAGAGTACGTGTCCAGCGGAGGGGCTTTGGACCGGTGTGTGACGGCTATGGTGTGTCGT
 ACATCCTTGTGGGAGAGAACCTCATCAATTTCCACATTTCTTCAAGTTCTTTGCCCTGAGACGGATTC
 TCATCGCTTTGAAGGCACCTGAAAGAAGCAATGACTGACATCATCACTTTGTTTGGTCTCAGTTCTAAT
 TCCAAAAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG221740 representing NM_001876
 Red=Cloning site Green=Tags(s)

MAEAHQAVAFQFTVTPDGIDLRLSHEALRQIYLSGLHSWKKKFIKFKNGIITGVYPASPSSWLIVVGVMT
 TTYAKIDPSLGIIAKINRTLETANCMSSQTKNVVSGVLFGTGLVVALIVTMRYSKLVLLSYHGWMFTEH
 GKMSRATKIWMGMVKIFSGRKPLYSFQTSPLPRLPVPVKDVTNRYLQSVRPLMKEEDFKRMTALAQDFA
 VGLGPRLQWYLKLSWWATNYVSDWWEYIYLRGRGPLMVNSNYAMDLLYILPHTHIQAARAGNAIHAIL
 LYRRKLDREEIKPIRLLGSTIPLCSAQWERMFNTRIPGEETDTIQHMRDSKHIVVYHRGRYFKVWLYHD
 GRLLKPREMEQQMQRILDNTSEPPGEARLAALTAGDRVPWARCRQAYFGRGKNKQSLDAVEKAFFVTL
 DETEEGYRSEDPDTSMDSYAKSLLHGRCYDRWFDKSFYVVFKNKMGMLNAEHSWADAPIVAHLWEYVMS
 IDSLQLGYAEDGHCKGDINPNIPYPTRLQWDIPGECQVEIETSLNTANLLANDVDFHSFPFVAFGKGIK
 KCRTSPDAFVQLALQLAHYKDMGKFCLEYEASMTRLFREGRTETVRSCTTESCDFVRAMVDAQTVEQRL
 KLFKLAASEKHQHMRYLAMTGSIDRHLFCLYVVS KYLAVESPFLEKVLSEPWRLSTSQTPOQQVELFDLE
 NNPEYVSSGGGFGPVADDGYGVS YILVGENLINFHIS SKFSCPETDSHRFGRHLKEAMTDIITLFLGSSN
 SKK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001876

ORF Size: 2319 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001876.3](#)

RefSeq Size: 2941 bp

RefSeq ORF: 2322 bp

Locus ID: 1374

UniProt ID: [P50416](#)

Cytogenetics: 11q13.3

Domains: Carn_acyltransf


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Adipocytokine signaling pathway, Fatty acid metabolism, PPAR signaling pathway

Gene Summary:

The mitochondrial oxidation of long-chain fatty acids is initiated by the sequential action of carnitine palmitoyltransferase I (which is located in the outer membrane and is detergent-labile) and carnitine palmitoyltransferase II (which is located in the inner membrane and is detergent-stable), together with a carnitine-acylcarnitine translocase. CPT I is the key enzyme in the carnitine-dependent transport across the mitochondrial inner membrane and its deficiency results in a decreased rate of fatty acid beta-oxidation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG221740