

## Product datasheet for **RG221777**

### **GPR64 (ADGRG2) (NM\_001079860) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	GPR64 (ADGRG2) (NM_001079860) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ADGRG2
Synonyms:	CBAVDX; EDDM6; GPR64; HE6; TM7LN2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)



Protein Sequence: >RG221777 representing NM\_001079860  
 Red=Cloning site Green=Tags(s)


MVFSVRQCGHVGRTVEVLLTFKIFLVIICLHVVLVTSLEEDTDNSSLSPPPAKLSVVSFAPSSNGTPDVT  
 LSLLPNETGVKPKQRNICNLSSICNDSAFFRGEIMFQYDKESTVPQNHITNGTLTGVLSLSELKRELN  
 KTLQTLSEYFIMCATAEAQSTLNCTFTIKLNNTMNACAVIAALERVKIRPMEHCCCSVRIPCSSPEEL  
 EKLQCDLQDPIVCLADHPRGPPFSSSQSIPVVPVPRATVLSQVPKATSFSAEPPDYSPVTHNVPSPIGEIQPL  
 SPQPSAPIASSPAIDMPPQSETISSPMPQTHVSGTPPPVKASFSSPTVSAPANVNTTSAPPVQTDIVNTS  
 SIDLENQVLQMEKALSLGSLEPNLAGEMINQVSRLHSPDMLAPLAQRLLKVVDDIGLQLNFSNTTIS  
 LTSPSLALAVIRVNASSFNNTTTFVAQDPANLQVSLETQAPENSIGTITLPSLMMNLPADMELASRVQF  
 NFFETPALFQDPSLENLSLISYVSSSVANLTVRNLTRNVTVTLKHINPSQDELTVRCVFWDLGRNGGRG  
 GWSDNCGSVKDRRLNETICTCSHLTSFGVLLDLRTSVLPAQMMALTFITYIGGLSSIFLSVTLVTYIA  
 FEKIRRDYPSKILIQCAALLLNLFLLDSWIALYKMQGLCISVAVFLHYFLLVSFTWMGLEAFHMYLA  
 LKVFNTYIRKYILKFCIVGWGVPVAVVTIILTISPNDYGLGSYGKFPNGSPDDFCWINNNAVFIYITVVG  
 YFCVIFLLNVSMFIVVLVQLCRIKKKKQLGAQRKTSIQDLRSIAGLTFLLGITWGF AFFAWGPVNVTFMY  
 LFAIFNTLQGGFFIFIFYCAKENVRKQWRRYLCCGKLRLENSDWSKTATNGLKKQTVNQGVSSSSNSLQ  
 SSSNSTNTLLVNDCSVHASGNNGASTERNVGSFVQNGDVCLHDFTKQKQHFNEKEDSCNGKGRMAL  
 RRTSKRGLHFIEQM


TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



**Plasmid Map:**


**ACCN:** NM\_001079860

**ORF Size:** 2985 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001079860.3](#)

**RefSeq Size:** 4829 bp

**RefSeq ORF:** 2988 bp

**Locus ID:** 10149

**UniProt ID:** [Q8IZP9](#)

**Cytogenetics:** Xp22.13

**Protein Families:** Druggable Genome, GPCR, Transmembrane

**Gene Summary:** This gene encodes a member of the G protein-coupled receptor family described as an epididymis-specific transmembrane protein. The encoded protein may be proteolytically processed as it contains a motif shown to be a protein scission motif in some members of this family (PMID: 11973329). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]