

## Product datasheet for **RG221790**

### DAZ2 (NM\_001005786) Human Tagged ORF Clone

#### Product data:

|                           |                                                                                |
|---------------------------|--------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                            |
| Product Name:             | DAZ2 (NM_001005786) Human Tagged ORF Clone                                     |
| Tag:                      | TurboGFP                                                                       |
| Symbol:                   | DAZ2                                                                           |
| Synonyms:                 | pDP1678                                                                        |
| Mammalian Cell Selection: | Neomycin                                                                       |
| Vector:                   | pCMV6-AC-GFP (PS100010)                                                        |
| E. coli Selection:        | Ampicillin (100 ug/mL)                                                         |
| ORF Nucleotide Sequence:  | >RG221790 representing NM_001005786<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAACTACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGC**

ATGCTGCTGCAAACTCTGAGACTCCAACTCAACCATCTCCAGAGAGGCCAGCACCCAGTCTTCATCAG  
 CTGCAGCTAGCCAAGGCTGGGTGTACCAGAAGGCAAAATCGTGCCAAACACTGTTTTGTTGGTGAAT  
 TGATGCTAGGATGGATGAACTGAGATTGGAAGCTGCTTTGGTAGATACGGTTCAGTGAAAGAAGTGAAG  
 ATAATCACGAATCGAACTGGTGTGTCAAAGGCTATGGATTTGTTTCGTTTGTTAATGACGTGGATGTCC  
 AGAAGATAGTAGGATCACAGATACATTTCCATGGTAAAAAGCTGAAGCTGGGCCCTGCAATCAGGAAACA  
 AAAGTTATGTGCTCGTCATGTGCAGCCACGTCTTTGGTAGTTAATCCTCCTCCTCCACCACAGTTTCAG  
 AACGTCTGGCGGAATCCAAACACTGAAACCTACCTGCAGCCCCAAATCACGCCGAATCCTGTAACCTCAGC  
 ACGTTCAGGCTTATTCTGCTTATCCACATTCACCAGGTCAGGTCATCACTGGATGTCAGTTGCTTGATATA  
 TAATTATCAGGAATATCCTACTTATCCCGATTCAGCATTTACAGTTCACCACTGGATATCAGTTGCCTGTA  
 TATAATTATCAGCCATTTCTGCTTATCCAAGTTCACCATTTACAGTTCAGTTCAGGATATCAGTTGCCTG  
 TATATAATTATCAGGCATTTCTGCTTATCCAAGTTCACCATTTACAGTTCACCACTGGATATCAGTTGCC  
 TGTATATAATTATCAGGCATTTCTGCTTATCCAAGTTCACCATTTACAGTTCACCACTGGATATCAGTTG  
 CCTGTATATAATTATCAGGCATTTCTGCTTATCCAAGTTCACCATTTACAGTTCACCACTGGATATCAGT  
 TGCTGTATATAATTATCAGGCATTTCTGCTTATCCAATTCAGCAGTTCAGGTCACCACTGGATATCA  
 GTTCCATGTATACAATTACCAGATGCCACCGCAGTGCCCTGTTGGGGAGCAAAGGAGAAATCTGTGGACC  
 GAAGCATACAATGGTGTATCTGTCTGTTTAAATCCAGAGAAGAGAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

**Protein Sequence:** >RG221790 representing NM\_001005786  
 Red=Cloning site Green=Tags(s)


MSAANPETPNSTISREASTQSSSAAASQGWVLPPEGKIVPNTVFVGGIDARMDTEIGSCFGRYGSVKEVK  
 IITNRTGVSKGYGFVSFVNDVDVQKIVGSQIHFGKLLKLGPAIRKQKLCAHVVQPRPLVNNPPPPQFQ  
 NVWRNPNTETYLQPIITPNPTQHVQAYSAYPHSPGQVITGCQLLVYNYQEYPTYPDSAFQVTTGYQLPV  
 YNYQFPAYPSSPFQVTAGYQLPVYNYQAFPAYPSSPFQVTTGYQLPVYNYQAFPAYPSSPFQVTTGYQL  
 PVYNYQAFPAYPSSPFQVTTGYQLPVYNYQAFPAYPNSAVQVTTGYQFHVYNYQMPPQCPVGEQRRNLWT  
 EAYKWWYLVCLIQRRD

TRTRPLE – GFP Tag – V

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001005786

**ORF Size:** 1098 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001005786.2](#), [NP\\_001005786.2](#)

**RefSeq Size:** 3277 bp

**RefSeq ORF:** 1101 bp


**Locus ID:** 57055

**UniProt ID:** [Q13117](#)

**Cytogenetics:** Yq11.223

**Gene Summary:** This gene is a member of the DAZ gene family and is a candidate for the human Y-chromosomal azoospermia factor (AZF). Its expression is restricted to premeiotic germ cells, particularly in spermatogonia. It encodes an RNA-binding protein that is important for spermatogenesis. Four copies of this gene are found on chromosome Y within palindromic duplications; one pair of genes is part of the P2 palindrome and the second pair is part of the P1 palindrome. Each gene contains a 2.4 kb repeat including a 72-bp exon, called the DAZ repeat; the number of DAZ repeats is variable and there are several variations in the sequence of the DAZ repeat. Each copy of the gene also contains a 10.8 kb region that may be amplified; this region includes five exons that encode an RNA recognition motif (RRM) domain. This gene contains one copy of the 10.8 kb repeat. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

## Product images:



Circular map for RG221790