

## Product datasheet for **RG221975**

### Myosin Phosphatase 2 (PPP1R12B) (NM\_002481) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Myosin Phosphatase 2 (PPP1R12B) (NM_002481) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PPP1R12B
Synonyms:	MYPT2; PP1bp55
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide  
Sequence:**

>RG221975 representing NM\_002481  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGCGGAACTGGAGCACCTAGGAGGGAAGCGGGCAGAGTCGGCGCAATCGGGCGGCAGAGCAGCTTC  
 GGCCTGGCGGGGCTCGCTGACAGAGCAGGAGCCTGCGGAGCGACGAGGCGGGGCGGCAGCCGCTGAC  
 CAGGCGCGGGAGCCCCAGGGTCCGCTTCGAGGACGGTGTGTCTTTCTGGCCGCCTGCTCTAGCGGGGAC  
 ACCGACGAGGTGAGAAAGCTTCTGGCAAGAGGTGCTGATATCAACACGGTCAACGTGGACGGCTTGACAG  
 CCCTGCACCAGGCATGTATTGATGAAAAATTTGGACATGGTGAAGTTTCTGGTGGAGAACAGAGCCAAATGT  
 AAACCAGCAAGACAACGAGGGCTGGACACCCCTTCATGCAGCAGCTTCTGTGGCTATCTCAACATAGCA  
 GAGTATTTCAATTAATCACGGAGCCAGTGTAGTATTGTCAATAGTGAAGGTGAAGTTCCTCTGACCTTG  
 CAGAAGAGCCAGCCATGAAGGATCTTCTTCTGGAGCAAGTAAAGAAGCAAGGAGTTGATCTAGAGCAGTC  
 AAGAAAAGAAGAAGAGCAGCAGATGTTGCAGGATGCCCGCAGTGGCTCAACAGTGGGAAAAATAGAGGAT  
 GTGAGGCAGGCTCGCTCAGGGGCTACAGCCCTTCATGTGGCTGTGCCAAGGGCTACTCTGAAGTCCCTCA  
 GACTTTTAATTCAGGCTGGCTATGAACTCAATGTTCAAGGATTATGATGGCTGGACTCCCCTCCATGCTGC  
 TGCACACTGGGGAGTGAAGGAGGCTTGCTCCATCCTGGCAGAAGCACTTTGTGACATGGATATTCGAAAT  
 AAAGTGGCCAGACACCATTTGATGTGGCTGATGAGGGTCTCGTGGAGCATTTGGAGTTGCTCCAGAAGA  
 AGCAGAATGTGCTTCAAGTGAAGGAGACACGGAATAAACTCATTGAGTCAGATCTGAACAGCAAGAT  
 TCAGAGTGGGTTCTTTAAGAACAAGAGAAGATGCTCTATGAGGAGGAGACACCTAAGTCCCAAGAAATG  
 GAGGAAGAAAATAAAGAATCTAGTAGCTCCAGCTCAGAGGAGGAGGAAGTGAAGATGAAGCTTCTGAGT  
 CAGAAAAGTGAAGGAGGAGATAAAAAGCCAGAAGCCTTTGTCAATCATTCCAACCTGAAAGCAAGAG  
 TAGTATCACAGAGCAGATACCAGCACAGCTCAAATACCTTCTCTGCCTTCTCTGCTAGGAGGTTCTCT  
 TCTGGCCTTTTTAAAGCCAGAAGAGCCAAAGATGAATCTCCTTCTCATGGAGATTGGGACTGAGAA  
 AAAGTGGCAGCCACAACATGCTGAGTGAAGTGGCCAATTCAGGGAACCTATAAGGGACCGAGGCTCTTC  
 CATCTATCGCTCCTTCAAGCCCTCGGATTTCTGCTCTACTGGACAACAAGATAAGGAGAGAGAAAAC  
 AAAAGCTATATTAGTTCCTAGCACCCCGAAGCTCAACAGCACAAGTATTTGAAGAAAAGGAGAAACA  
 GAGAATCAGCTGTTAATCTAGTGAAGTGGCTCCTATACCCGGCAGCTATGGAGGGATGAAGCAAAAGG  
 AAATGAAATCCCACAGACAATTGCTCCCTCCACCTATGTATCACTTACTGAAAAGGACTCCTCACAAA  
 TCCAGGCGGACACAACAGCAGAGAAAACAGCAGACAATGTCTTCTAGCACCCGCTCTGTGTGATCA  
 CCAATCGCCCTTCTAGCACTGCCAATGGGGTTACAGCTACTCCTGTGCTCTCCATTACTGGAACAGA  
 TTCTCTGTGGAAGCCAGGGAGAAGAGGAGTCTATCTGACTCCTGTACGGGATGAGGAAGCAGAGTCT  
 TTACGGAAGCAGCTCCAGACAAGCTCGGCAGACACGAAGTCTACTCAAGGTGTACCCTAACAGACC  
 TTCAAGAAGCAGAAAGGACATTCAGCCGGTCGAGGGCAGAGAGGCAAGCTCAGGAGCAGCCTCGTGAGAA  
 GCCACAGACACTGAAGGGCTTGAGGGGAGCCCTGAGAAGCATGAGCCCTCAGCAGTTCAGCAACAGAA  
 GCTGGGGAGGGCCAGCAGCCCTGGGGCAGGAGTCTGGATGAAGAGCCTATCTGTGATCGCTGAGGTGCC  
 CAGCTCAGCCAGACAAAACCCACCCAGCAGTCTCCTTCTACGTCAAGACCCCTACTCTACACCCAGTTC  
 CCACCTGCTATGGACAAAATAGATTTTCAGTCCCTGATTCTGAGAGTTCAGAGACTACCACAAAACACTACA  
 ACTGCAAAGGAAATGGACAAAAATGAGAATGAAGAAGCAGATTTGGATGAGCAGTCTCTAAGAGGCTGT  
 CCATCCGAGAGAGGAGGGGCCAAGGAACGACGAAGAGGCACAGGCATCAATTTCTGGACAAAAGGATGA  
 GGATGAAACTGATGGCTCTGAAGAGGTCAAAGAAACGTGGCATGAAAGACTTTCTAGGTTGGAATCGGGA  
 GGATAGTAATCCTACAACCAGTATTCTTACGGTGACCGGGCTTCAGCAAGAGCCCGTGGGAGGCCCGGG  
 AGGCCCGCTAGCCACCCTGACCAGCCGTGTAAGAAGACAGCAACAGAGATTAAAAAACTCTATGA  
 GAGTGTCTGACTGAAAACAAAAACTGAAAACAAAACTCAGGAAGCCAGCTAGAGCTAGCAGATATA  
 AAGTCAAAGCTTGAGAAGGTGGCCAGCAGAAAACAAGAAAAGACCTCTGACCGATCATCAGTGTGGAGA  
 TGGAGAAAACGGGAGAGGGCAGCCTTGGAGCGCAAATGTCAGAAATGGAGGAAGAAATGAAGGTGTTAAC  
 AGAACTGAAATCCGACAACCAGAGGCTGAAAGATGAAAATGGTGCCTCATCAGAGTCATCAGCAAACCTG  
 TCCAAG

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG221975 representing NM\_002481  
 Red=Cloning site Green=Tags(s)

MAELEHLGGKRAESARMRAEQLRRWRGSLTEQPAERRGAGRQPLTRRGSPRVRFEDGAVFLAACSSGD  
 TDEVKLLARGADINTVNV DGLTALHQACIDENLDMVKFLVENRANVNQDNEGWTPLHAAAACGYLNIA  
 EYFINHGASVGI VNSEGEVPSDLAE E PAMKDLLLEQVKKQGV DLEQSRKEEEQQMLQDARQWLN S G K I E D  
 VRQARSGATALHVAAAAGYSEVLRLLIQAGYELNVQDYDGTPLHAAAHWGVKEAC SILAEALCDMDIRN  
 KL G Q T P F D V A D E G L V E H L E L L Q K K Q N V L R S E K E T R N K L I E S D L N S K I Q S G F F K N K E K M L Y E E E T P K S Q E M  
 E E E N K E S S S S S E E E E G E D E A S E S E T E K E A D K K P E A F V N H S N S E S K S S I T E Q I P A P A Q N T F S A S S A R R F S  
 S G L F N K P E E P K D E S P S S W R L G L R K T G S H N M L S E V A N S R E P I R D R G S S I Y R S S S P R I S A L L D N K D K E R E N  
 K S Y I S S L A P R K L N S T S D I E E K E N R E S A V N L V R S G S Y T R Q L W R D E A K G N E I P Q T I A P S T Y V S T Y L K R T P H K  
 S Q A D T T A E K T A D N V S S T P L C V I T N R P L P S T A N G V T A T P V L S I T G T D S S V E A R E K R R S Y L T P V R D E E A E S  
 L R K A R S R Q A R Q T R R S T Q G V T L T D L Q E A E R T F S R S R A E R Q A Q E Q P R E K P T D T E G L E G S P E K H E P S A V P A T E  
 A G E G Q Q P W G R S L D E E P I C H R L R C P A Q P D K P T T P A S P S T S R P S L Y T S S H L L W T N R F S V P D S E S S E T T N T T  
 T A K E M D K N E N E E A D L D E Q S S K R L S I R E R R R P K E R R R G T G I N F W T K D E D E T D G S E E V K E T W H E R L S R L E S G  
 G S N P T T S D S Y G D R A S A R R E A R E A R L A T L T S R V E E D S N R D Y K K L Y E S A L T E N Q L K T K L Q E A Q L E L A D I  
 K S K L E K V A Q Q K Q E K T S D R S S V L E M E K R R R A L E R K M S E M E E M K V L T E L K S D N Q R L K D E N G A L I R V I S K L  
 S K


TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



**Plasmid Map:**


**ACCN:** NM\_002481

**ORF Size:** 2946 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<a href="#">NM_002481.4</a>
RefSeq Size:	11108 bp
RefSeq ORF:	2949 bp
Locus ID:	4660
UniProt ID:	<a href="#">O60237</a>
Cytogenetics:	1q32.1
Domains:	ANK
Protein Families:	Druggable Genome
Protein Pathways:	Vascular smooth muscle contraction
Gene Summary:	<p>Myosin phosphatase is a protein complex comprised of three subunits: a catalytic subunit (PP1c-delta, protein phosphatase 1, catalytic subunit delta), a large regulatory subunit (MYPT, myosin phosphatase target) and small regulatory subunit (sm-M20). Two isoforms of MYPT have been isolated--MYPT1 and MYPT2, the first of which is widely expressed, and the second of which may be specific to heart, skeletal muscle, and brain. Each of the MYPT isoforms functions to bind PP1c-delta and increase phosphatase activity. This locus encodes both MYTP2 and M20. Alternatively spliced transcript variants encoding different isoforms have been identified. Related pseudogenes have been defined on the Y chromosome. [provided by RefSeq, Oct 2011]</p>